900 series RELAX + upper arm blood pressure monitor

BPA-100AHBT-WT
Doc # L-03006, Rev.1

5 year limited warranty

IMPORTANT PRODUCT NOTICES AND SAFETY INSTRUCTIONS

When using your blood pressure monitor, basic precautions should always be followed. Please read and follow all instructions and warnings before using this product.

Save these instructions for future reference.

This device is intended for over-the-counter home use in adults 18 years and older with an arm circumference ranging from 9 inches to 17 inches (approx. 23cm to 43cm).

The device is equipped with an Advanced IHB (irregular heartbeat) detection feature to collect and analyze pulses. If the specific irregular heartbeats are detected and it may affect blood pressure reading with deviation, the device will give the user a warning signal. The feature can inform the user that the measured blood pressure reading may be inaccurate from the specific irregular heartbeats are detected.

This device features a built-in "Bluetooth® Data Transmission" function, which enables the device to automatically transmit measuring results to an auto-paired Bluetooth-enabled device. This function allows users to synchronize to current date and time, and check the battery status with the Homedics Health+ app.

When using your blood pressure monitor, basic precautions should always be followed. Please read and follow all instructions and warnings before using this product. Save these instructions for future reference.

• Please note that this is a home health care product only and this manual and device are not intended to serve as a substitute for the advice of a physician or medical professional.

• This device uses oscillometric method to measure systolic and diastolic blood pressure, as well as heart rate.

• DO NOT use this device for diagnosis or treatment of any health problem or disease. Measurement results are for reference only. Consult a health care professional for interpretation of pressure measurements. Contact your physician if you have or suspect any medical problem. Do not change your medications without the advice of your physician or health care professional.

• Proper cuff size is critical for accurate measurements. Follow the instructions in this manual and printed on the arm cuff to ensure the arm cuff is used properly.

• DO NOT take a measurement in a low (less than 41°F/5°C) and high (more than 104°F/40°C) temperature, nor in a place outside humidity ranges (15% ~ 93% R.H.), and altitude ranges (700 ~ 1060 hPa), or you may get inaccurate readings.

• This product is not suitable for:

- Pregnant women
- People undergoing intravenous injection on any limb
- People currently in a dialysis treatment
- Pregnant women in preeclampsia condition

• Those who have had a mastectomy surgery (especially when lymph nodes have been removed), it is recommended to take measurements on the unaffected side.

• This device may have difficulty determining proper blood pressure for users with irregular heartbeat, diabetes, liver disease, kidney disease, poor circulation of the blood or for users who have suffered a stroke. Please consult your healthcare professional before using this device. Excessive use may result in blood overflow interference, which is likely to cause uncomfortable sensations, such as partial subcutaneous hemorrhage, or temporary numbness to your arm. In general, these symptoms should not last long. However, if you continue to experience these sensations, please seek advice from a medical professional.

• When used along with other electronic medical equipment on the same arm, pressurization of the cuff may cause the other devices to temporarily malfunction.

• The pulse display is not suitable for checking the frequency of heart pacemakers.

• Electromagnetic interference: This device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. cell phones, microwave ovens). These may lead to temporary impairment of measurement accuracy.

• Consider the electromagnetic compatibility of the device (ex. power disturbance, radio frequency interference, etc.) Please use this device indoors only.

• Use blood pressure monitor only for its intended use.

• DO NOT wrap the cuff around body parts other than your arm.

• Not for use by or on persons under the age of 18.

• DO NOT use this device on infants, children, or those who cannot express their own intention.

• DO NOT plug or unplug the adapter power cord with wet hands.

• Please use only the specified AC adapter (not included, sold separately) with this monitor or 4.15V "AA" alkaline batteries for power supply.

• Please rest for at least 5-10 minutes before taking a measurement.

• To allow your blood vessels to return to the condition prior to taking the measurement, please wait at least 3-5 minutes between measurements. You may need to adjust the wait time according to your personal physiological situation.

• Wait 30-45 minutes before measurement if you've just consumed caffeinated beverages or smoked cigarettes.

• The applied part is the cuff.

• The patient is the intended operator.

Blood pressure measurements determined with this device are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultatory method, within the accuracy limits prescribed by the American National Standard, Manual, electronic, or automated sphygmomanometers.

In the event that the device needs to be checked for calibration, contact the distributor.

ABOUT BLOOD PRESSURE

What is blood pressure?
Blood pressure is the measurement of the force of blood pushing against the walls of the arteries. Arterial blood pressure is constantly fluctuating during the course of the cardiac cycle. The pressure measured when the heart contracts and sends blood out of the heart is systolic (highest) blood pressure. The pressure measured when the heart dilates with blood flowing back into the heart is called diastolic (lowest) blood pressure.

Why measure your blood pressure?

Among today's various health problems, those associated with high blood pressure are very common. High blood pressure dangerously correlates with cardiovascular diseases. Therefore, blood pressure monitoring is important for identifying those at risk.

Why do my readings vary?

Blood pressure is a body parameter that is subject to normal variations throughout the day. A single reading that is different from yours or your doctor's readings are not necessarily inaccurate. The average of several readings, taken under similar conditions, using the same arm is preferred for accurate blood pressure readings. Medical physicians generally recommend the "Rule of 3", where you are encouraged to take your blood pressure three times in a row (at 3 - 5 minute intervals), three times a day for three days. After three days you can average all the results and this will give you an accurate idea of what your blood pressure really is.

Why Are My Readings Different Than Those Taken at My Doctor's Office?

Many experience a phenomenon called "White Coat Hypertension" when measured by a doctor. White Coat Hypertension refers to blood pressure that rises above its usual level when measured in a clinical setting, such as a doctor's office.

BLOOD PRESSURE STANDARD

The table below contains defined levels for hypertension that are publicly available from the American Heart Association® (AHA 2017). Users can compare their own blood pressure readings against these defined levels to determine if they may be potentially at increased risk.

Blood Pressure Category	Systolic mm Hg (upper number)	Diastolic mm Hg (lower number)	Indicator Color
Normal	<120	and	<80
Elevated	120-129	and	<80
High Blood Pressure (Hypertension) Stage 1	130-139	or	80-89
High Blood Pressure (Hypertension) Stage 2	≥140	or	≥90
Hypertension Crisis (consult your doctor immediately)	>180	and/or	>120

*Source: American Heart Association® (AHA) 2017

Blood pressure tends to go up and down, even in people who normally don't have high readings. If your numbers stay above the "normal" range most of the time, you may be at increased risk and should consult your physician. Although one can easily find where their own blood pressure readings fall on this table, this monitor comes equipped with a Risk Category Index that automatically compares each reading to the defined levels and provides a helpful cue if your reading falls into one of the stages that could potentially indicate increased risk. See Risk Category Index section for more information on this feature.

Please note that cues provided by this monitor are only intended to assist you in using this table. The table and cues are only provided for convenience to help you understand your non-invasive blood pressure reading as it relates to the American Heart Association® (AHA 2017) information.

They are not a substitute for a medical examination or diagnosis by your physician. It is important for you to consult with your pressure range as well as the point at which you may actually be considered to be at risk.

HOW THIS BLOOD PRESSURE MONITOR WORKS

This monitor uses Smart Measure® inflation technology to detect your blood pressure. Once the cuff is wrapped around your upper arm, with a press of the START/STOP button, the monitor will turn on and the inflation will automatically start, creating pressure around the arteries inside the upper arm. Within the cuff is a gauge which senses the fluctuations (oscillations) in pressure. The fluctuation measured represents the degree of intensity that your arteries are contracting with each heartbeat, which is also a result of the pressure that the cuff has placed on the arm. The monitor measures these contractions and converts the information to a digital value. Once the measurement is complete, the cuff will automatically deflate, and the result is displayed on the screen. Please note that any muscle movement during inflation will cause measurement error. When measurement is complete, the monitor will display your systolic pressure, diastolic pressure, and pulse readings. Before measurement, it is suggested that you sit quietly for 15 minutes before measurement as measurements taken in a relaxed state have greater accuracy. The monitor automatically finds where your measurement result fall on the American Heart Association® (AHA 2017) table and provides a cue if your reading falls into one of the stages that could potentially indicate increased risk. Please refer to the Risk Category Index section for more information on this feature.

NAME/FUNCTION OF EACH PART

This monitor uses Smart Measure® inflation technology to detect your blood pressure. Once the cuff is wrapped around your upper arm, with a press of the START/STOP button, the monitor will turn on and the inflation will automatically start, creating pressure around the arteries inside the upper arm. Within the cuff is a gauge which senses the fluctuations (oscillations) in pressure. The fluctuation measured represents the degree of intensity that your arteries are contracting with each heartbeat, which is also a result of the pressure that the cuff has placed on the arm. The monitor measures these contractions and converts the information to a digital value. Once the measurement is complete, the cuff will automatically deflate, and the result is displayed on the screen. Please note that any muscle movement during inflation will cause measurement error. When measurement is complete, the monitor will display your systolic pressure, diastolic pressure, and pulse readings. Before measurement, it is suggested that you sit quietly for 15 minutes before measurement as measurements taken in a relaxed state have greater accuracy. The monitor automatically finds where your measurement result fall on the American Heart Association® (AHA 2017) table and provides a cue if your reading falls into one of the stages that could potentially indicate increased risk. Please refer to the Risk Category Index section for more information on this feature.

DISPLAY SYMBOL EXPLANATIONS

DISPLAY SYMBOL EXPLANATIONS

Date/Time Indicator

Memory Average

Excessive Body Motion Detector

Cuff Fitness-Detection

Pulse Rate Symbol

Heart Rate

Systolic Pressure

Risk Category Indicator

Diastolic Pressure

Low Battery Symbol: Appears when batteries should be replaced.

Pulse Symbol: Shows the heart rate per minute. Once pulse is detected, the symbol flashes with each pulse beat.

Risk Category Index: See Risk Category Index Section for more information. The arrow points out the specific Risk Category that your measurement reading falls in.

Avg. 3: Average of the Last 3 Measurements: This symbol appears when the LCD displays the average value of last 3 readings.

Excessive Body Motion Detector: Displayed if body movement is detected during measurement, especially, the movement on the arm the blood pressure monitor is worn on. If the cuff is worn improperly, or the shape of the upper arm is unusual for example, the circumference of the upper arm differs largely from the circumference of the forearm, excessive gap might exist between the arm cuff and the arm.

Note: The measured blood pressure reading may not be accurate if the icon is displayed.

Cuff Fitness Detection Symbol: Displayed if the cuff was wrapped incorrectly, which is too tight or too loose. This function aids in detecting if the cuff is wrapped properly.

Advanced IHB Detection Symbol: Displayed when there is Advanced IHB negative.

Irregular Heartbeat + Advanced Detection: See Irregular Heartbeat + Advanced Detection section for more information.

Bluetooth® Symbol: Appears when Bluetooth® is active, and when data is transmitting to your mobile device.

*These symbols are only generic symbols to indicate the detection (or non-detection) of IHB, not actual heartbeat waveforms detected as a function of time.

If any of the following letters and numbers appear in the area that systolic pressure should be displayed, an error has occurred with your reading. See Troubleshooting section of this manual for more information.

Display Explanation Symbols:

User 1: Appears when the monitor is operated by User 1.

User 2: Appears when the monitor is operated by User 2.

Low Battery Symbol: Appears when batteries should be replaced.

Pulse Symbol: Shows the heart rate per minute. Once pulse is detected, the symbol flashes with each pulse beat.

Risk Category Index: See Risk Category Index Section for more information. The arrow points out the specific Risk Category that your measurement reading falls in.

Avg. 3: Average of the Last 3 Measurements: This symbol appears when the LCD displays the average value of last 3 readings.

Excessive Body Motion Detector: Displayed if body movement is detected during measurement, especially, the movement on the arm the blood pressure monitor is worn on. If the cuff is worn improperly, or the shape of the upper arm is unusual for example, the circumference of the upper arm differs largely from the circumference of the forearm, excessive gap might exist between the arm cuff and the arm.

Note: The measured blood pressure reading may not be accurate if the icon is displayed.

Cuff Fitness Detection Symbol: Displayed if the cuff was wrapped incorrectly, which is too tight or too loose. This function aids in detecting if the cuff is wrapped properly.

Advanced IHB Detection Symbol: Displayed when there is Advanced IHB negative.

Irregular Heartbeat + Advanced Detection: See Irregular Heartbeat + Advanced Detection section for more information.

Bluetooth® Symbol: Appears when Bluetooth® is active, and when data is transmitting to your mobile device.

*These symbols are only generic symbols to indicate the detection (or non-detection) of IHB, not actual heartbeat waveforms detected as a function of time.

If any of the following letters and numbers appear in the area that systolic pressure should be displayed, an error has occurred with your reading. See Troubleshooting section of this manual for more information.

2. Connect the AC adapter with the AC adapter jack in the back of the unit as shown below.

3. Use only the AC adapter specified.

NOTE:

• When the AC adapter is your main power supply, make sure the adapter plug can be easily removed from the unit.

• Please unload the batteries when operating with an AC adapter for an extended period of time. Leaving the batteries in the compartment for a long time may cause leakage, which may lead to damage to the unit.

• No batteries are needed when operating with an AC adapter.

• Date and time will need to be reset if AC adapter is unplugged even if the unit has the batteries in the battery compartment.

If Advanced IHB is detected during measurement, the Advanced IHB positive symbol is displayed.

If Advanced IHB positive symbol appears frequently, we recommend the patient to seek professional medical advice. And please note that the device does not replace a cardiac examination, but serves to detect specific pulse irregularity at an early stage.

DATE & TIME SET PROCEDURE

Date & time can be set by two methods, either sync automatically using your Bluetooth® mobile device, or set manually using the USER/+ and buttons on the monitor.

If setting the date & time by syncing with your mobile device, it is important that this is done prior to taking any measurements to ensure the date & time are accurate.

To set using your device with Bluetooth® wireless technology:

1. Before taking the first measurement, open the Homedics Health+ app on your device. If you have not already installed the Homedics Health+ app onto your device, it is available to download on the App Store® and Google Play®. Make sure you have the app downloaded and open on your mobile device before trying to sync the date & time.

2. Press START/STOP button to switch on the monitor.

</

RISK CATEGORY INDEX

This monitor comes equipped with a Risk Category Index that automatically compares each reading to defined levels established by the American Heart Association® (AHA 2017) as described earlier in this manual, and provides a helpful cue if your reading falls into one of the stages that could potentially indicate increased risk. No cue is given if the reading falls in the normal range as defined by the AHA. Please note that cues provided by this monitor are only intended to assist you in using this table. The table and cues are only provided for convenience to help you understand your non-invasive blood pressure reading as it relates to the AHA information. They are not a substitute for a medical examination or diagnosis by your physician. It is important for you to consult with your physician regularly. Your physician will tell you your normal blood pressure range as well as the point at which you may actually be considered to be at risk.

Blood Pressure Category	Systolic mm Hg (upper number)		Dia stolic mm Hg (lower number)	Indicator Color
Normal	<120	and	<80	Green
Elevated	120-129	and	<80	Yellow
High Blood Pressure (hypertension) Stage 1	130-139	or	80-89	Red
High Blood Pressure (hypertension) Stage 2	≥140	or	≥90	
Hypertension Crisis (consult your doctor immediately)	>180	and/or	>120	

*Source: American Heart Association® (AHA) 2017

After measurement, the LCD displays the systolic and diastolic pressure, heart rate, date and time along with Risk Category Indicator symbols. Compare the symbol of LCD display, to know the classification of your blood pressure based on American Heart Association® standard (AHA 2017).

Note: When a person's systolic and diastolic pressures fall into different categories, the higher category should apply.

USING THE MEMORY FUNCTION

The monitor can store up to 240 memories for 2 users (120 each user), and automatically replaces the oldest data with new one.

After each measurement, the systolic and diastolic pressure, heart rate, pulse, Risk Category Indicator Bar, Advanced IHB Detection feature-related symbol with date and time will be automatically stored.

Note: If the Advanced IHB Detection measurement is taken without a blood pressure measurement, your data cannot be stored.

RECALLING VALUES FROM MEMORY

1. Press the **USER/+** button to select User 1 or User 2.
2. Press the **M** button to enter Memory Recall Mode. The LCD will display the calculated average of the last 3 measurement results first.

3. Press the **M** button again. The LCD will display the latest measurement. Press **M** or **8** buttons to scroll through all stored measurements.

Note: Each press of the **USER/+** button will show the prior reading, newest to oldest, in a cycling order. (example, avg 3, prior reading 1, prior reading 2, ...)

• Each press of the **8** button will show the more current reading, oldest to newest, in a cycling manner. (example, prior reading 2, prior reading 1, avg 3, ...)

4. To stop reading memories, press the START/STOP **STOP** button, and switch to standby mode.

CLEARING VALUES FROM MEMORY

1. Press the **USER/+** button to select User 1 or User 2.
2. Press the **M** button to enter Memory Recall Mode.
3. Press and hold the **USER/+** and the **8** buttons at the same time. The data will be erased automatically.

4. To confirm the data for the selected user has been erased, press the **M** button. No data should appear.

Note: Once deleted, your readings cannot be restored.

IMPORTANT NOTES REGARDING YOUR BLOOD PRESSURE MEASUREMENT

- It is suggested that you take your measurements at the same time each day and use the same arm for consistency.
- Users should wait a minimum of 5 minutes before taking additional measurements. More time may be necessary depending upon your physiology.
- The measurement results that users receive are for reference only. If users have any blood pressure concerns, please consult a physician.
- Once inflation reaches 300 mmHg, the unit will deflate automatically for safety reasons.
- This device may have difficulty determining the proper blood pressure for users with irregular heartbeat, diabetes, poor circulation of blood, kidney problems, or for users who have suffered a stroke.

CARE, MAINTENANCE & CLEANING

- Clean the blood pressure monitor body and cuff carefully with a slightly damp, soft cloth. Do not press. Do not wash cuff or use chemical cleaner on it. Never use thinner, alcohol or petrol (gasoline) as cleaner.
- Make sure the cuff is completely dry before using.
- Leaky batteries can damage the unit. Remove the batteries when the unit will not be used for a long time.
- Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.
- If the unit is stored near freezing, allow it to acclimate to room temperature before use.
- This blood pressure monitor is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems with this device, please contact Homedics Consumer Relations (contact information can be found on the warranty page).
- DO NOT immerse the unit in water as this will result in damage to the unit.
- DO NOT subject the monitor or cuff to extreme temperatures, humidity, moisture, or direct sunlight. Protect from dust.

- DO NOT fold the cuff and tubing tightly.
- DO NOT disassemble the monitor or cuff. If in need of repair, refer to the warranty section of this manual.
- DO NOT subject the monitor to extreme shocks (do not drop on floor).
- DO NOT twist the cuff in any way.
- DO NOT inflate the cuff unless wrapped around arm.
- DO NOT wrap the cuff around body parts other than your arm.
- DO NOT drop or insert any object into any opening or hose.
- To avoid accidental strangulation keep this product away from children. Do not drape tube around neck.
- This monitor may not meet its performance specifications if stored or used outside of these temperature and humidity ranges:

STORAGE/TRANSPORTATION ENVIRONMENT

Temperature: -13°F ~ 158°F (-25°C ~ 70°C)

Humidity: less than 93% RH

OPERATING ENVIRONMENT

Temperature: 41°F ~ 104°F (5°C ~ 40°C)

Humidity: 15% ~ 93% RH

Atmospheric Pressure: 700hPa-1060hPa

- Always store the unit in the storage case after use. It is intended to be transported or stored in a carrying case between uses.
- Do not place the device directly under sunlight, in high temperature, or in humid or dusty places.
- Do not store the device in extremely low (less than -13°F/-25°C) and high (more than 122°F/50°C) temperature, nor in a place its humidity exceeds 93% R.H.

POTENTIAL FOR ELECTROMAGNETIC INTERFERENCE

To avoid inaccurate results caused by electromagnetic interference between electrical and electronic equipment, do not use the device near a cell phone or microwave oven. For most wireless communication devices, it is recommended to maintain a distance of 10.8 feet (3.3m) in order to avoid electromagnetic interference.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

FEDERAL COMMUNICATIONS COMMISSION COMPLIANCE STATEMENT

Changes or modifications to this equipment not expressly approved by the manufacturer could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

1. Reorient or relocate the receiving antenna.
2. Increase the separation between the equipment and receiver.

Connect the equipment into an outlet on a circuit different from that to which the receiver is connected. Consult the dealer or an experienced radio/TV technician for help.

RF exposure:

1. The equipment complies with FCC RF exposure limits set forth for an uncontrolled environment.
2. The equipment must not be co-located or operated in conjunction with any other antenna or transmitter.

FCC Label Compliance Statement:

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) This device must accept any interference received, including interference that may cause undesired operation.

Note: Changes or modifications not expressly approved by the manufacturer could void the user's authority to operate the equipment.

ELECTROMAGNETIC COMPATIBILITY (EMC)

Guidance and manufacturer's declaration – electromagnetic emissions

The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR11	Group 1	RF energy is used only to maintain device's operation. Therefore, its RF emissions are so low that it's not likely to cause any interference in nearby.
RF emissions CISPR11	Class B	
Harmonic Emissions IEC61000-3-2	Class A	The device is suitable for use in all establishments, including domestic establishments, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact discharge ±15 kV air discharge	±8 kV contact ±15 kV air discharge	In the case of air discharge testing, the climatic conditions shall be within the following ranges: Ambient Temperature: 15°C~35°C Relative Humidity: 30%~60%.
Power frequency (50 or 60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 or 60 Hz	30 A/m 50 or 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Electrical fast transient/burst EC 61003.4-4	±2 kV for power supply lines ±2 kV for input/output lines	±2 kV for power supply lines ±2 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.

Emissions Test	Compliance	Electromagnetic Environment - Guidance
Surge IEC 61000-4-5	AC Power port ±1 KV Line to Line	AC Power port ±1 KV Line to Line
Interruptions and voltage variations on power supply input lines IEC 61003 4-11	OR6 UT: 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0% UT: 1 cycles 70 % UT; 25/30 cycles 0% UT; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.


TROUBLESHOOTING

If any abnormality arises during use, please check the following points:

Symptoms	Possible Causes	Correction
Unit does not turn on when the START/STOP button is pushed.	Batteries have run down.	Replace them with six new AA (LR6) alkaline batteries.
	Battery polarities (+/-) have been positioned incorrectly.	Re-insert the batteries in the correct positions.
E E Measuring Error Symbol appears when blood pressure value displayed is excessively low or high.	The Upper Arm Cuff has been positioned on the arm incorrectly.	Re-wrap the cuff properly so that it is positioned correctly. Take new measurement.
	Did you talk or move during measurement?	Keep arm steady during measurement. Measure again. Refer to "Measurement Procedure" instructions.
	Shaking of the arm with the cuff on.	Check cuff connection. See "Using the Upper Arm Cuff" section. Take new measurement.
E I Measuring Error Symbol	Air circuit abnormality. Cuff tube may not be plugged into monitor correctly.	Check cuff connection. See "Using the Upper Arm Cuff" section. Take new measurement.
E 2 Measuring Error Symbol	Inflation pressure exceeding 300 mmHg.	Switch the unit off and then turn back ON. Take new measurement.
E 3 Measuring Error Symbol	Error determining measurement data.	Re-wrap the cuff properly so that it is positioned correctly. Take new measurement.
E 4 Measuring Error Symbol	If BPA-100AHBT-WT cannot be connected to paired Bluetooth-enabled device over 45 seconds, LCD will display error message "E4" and Bluetooth® will be turned off.	Please press STOP button for one time to start Bluetooth® function.
E 5 Try Again Measuring Error Symbol	Finger hasn't be placed on Advanced IHB detection area.	Keep finger gently place and well-covered on Advanced IHB detection area and measure again.
	Finger moved away from the detection area when the measurement has not been completed yet.	Measure again and don't move away your finger before measurement completed.
E 6 Try Again Measuring Error Symbol	Finger press too hard on Advanced IHB detection area.	Measure again. Gently place finger on Advanced IHB detection area.
	Cold Finger and weak pulse signals can't determine Advanced IHB measurement data.	Measure again. Keep finger warm and gently place on Advanced IHB detection area.
E 7 Measuring Error Symbol	Pulse signals could not be detected continuously by the Advanced IHB detection area for a period.	Place the finger on Advanced IHB detection area and keep steady. Measure again.
Excessive Body Motion Detector	Body movement during measurement, especially, the movement on the arm the blood pressure monitor is worn on. e.g. Talking, moving or shaking of the arm with the cuff on while measurement.	Measure again. Keep arm steady during measurement.
Notice:	The measured blood pressure reading may not be accurate if the icon is displayed.	If you have any question about the cuff wearing and/or measurement result, please consult your healthcare professional.
Test specifications for enclosure port immunity to RF wireless communications equipment.		
Test frequency (MHz)	Modulation	Immunity Test Level (V/m)
385	Pulse modulation 18 Hz	27
450	FM± 5 kHz deviation 1kHz sine	28
710	Pulse modulation 217 Hz ^a	9
745		
780		
810	Pulse modulation 18 Hz ^a	28
870		
930		
1720	P	

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*Source: American Heart Association^{MD} (AHA) 2017

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