700 series 5-day trend upper arm blood pressure monitor

IMPORTANT PRODUCT NOTICES AND SAFETY INSTRUCTIONS

When using your blood pressure monitor, basic precautions should always be followed. Please read and follow all instructions and warnings before using this product. Save these instructions for future reference. This device is intended for over-the-counter home use in adults 18 years and older with an arm circumference ranging from 9 inches to 17 inches (approx. 23cm to 43cm). The device can accurately measure blood pressure in pregnant patients including those with known or suspected in preeclampsia condition. This device detects the appearance of irregular heartbeats during measurement; an indication symbol will appear with measuring reading. The BP Category Indicator will show the information with the readings on the screen for the user tracking their blood pressure level. This device features a built-in "Bluetooth" Data Transmission function, which enables the device to automatically transmit measuring results to a paired Bluetooth-enabled device. This function allows users to synchronize to current date and time, and check the battery status on the Homedics Health+ app. When using your blood pressure monitor, basic precautions should always be followed. Please read and follow all instructions and warnings before using this product. Save these instructions for future reference.

- Please note that this is a home healthcare product only and this manual and device are not intended to serve as a substitute for the advice of a physician or medical professional.
- This device uses the oscillometric method to measure systolic and diastolic blood pressure, as well as heart rate.
- DO NOT use this device for diagnosis or treatment of any health problem or disease. Measurement results are for reference only. Consult a healthcare professional for interpretation of pressure measurements. Contact your physician if you have or suspect any medical problem. Do not change your medications without the advice of your physician or healthcare professional.
- Proper cuff size is critical for accurate measurements. Follow the instructions in this manual and printed on the arm cuff to ensure the arm cuff is used properly.
- If you are pregnant, you should pay more attention to your blood pressure changes, because during this time it may change drastically.
- This monitor is clinically validated for use in pregnancy and preeclampsia. When you detect unusual readings in pregnancy, you should measure again after taking some rest. If the reading is still abnormal, consult your doctor or gynecologist.
- This product is not suitable for:
 - People with arrhythmias
 - People undergoing intravenous injection on any limb
 - People currently in dialysis treatment
 - Those who have had a mastectomy surgery (especially when lymph nodes have been removed), it is recommended to take measurements on the unaffected side.
- This device may have difficulty determining proper blood pressure for users with irregular heartbeat, diabetes, liver disease, kidney disease, poor circulation of the blood or for users who have suffered a stroke. Please consult your healthcare provider before using this device.
- Excessive use may result in blood overflow interference, which is likely to cause uncomfortable sensations, such as partial subcutaneous hemorrhage, or temporary numbness to your arm. In general, these symptoms should not last long. However, if you continue to experience these sensations, please seek advice from a medical professional.
- When used along with other medical electronic equipment on the same arm, pressurization of the cuff may cause the other devices to temporarily malfunction.
- The pulse display is not suitable for checking the frequency of heart pacemakers.
- Electromagnetic interference: The device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens). These may lead to temporary impairment of measurement accuracy.
- Consider the electromagnetic compatibility of the device (ex. power disturbance, radio frequency interference etc.). Please use this device indoors only.
- The pulse display is not suitable for checking the frequency of heart pacemakers.
- Electromagnetic interference: This device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. cell phones, microwave ovens). These may lead to temporary impairment of measurement accuracy.
- Consider the electromagnetic compatibility of the device (ex. power disturbance, radio frequency interference etc.). Please use this device indoors only.
- Use blood pressure monitor only for its intended use.
- DO NOT wrap the cuff around body parts other than your arm.
- DO NOT use on persons under 18 years of age.
- DO NOT use this device on infants, children, or those who cannot express their own intention.
- DO NOT plug or unplug the adapter power cord with wet hands.
- Please use only the USB to DC power cable included with this monitor or 1.5V AAA alkaline batteries for power supply.
- Please rest for at least 5-10 minutes before taking a measurement.
- To allow your blood vessels to return to the condition prior to taking the measurement, please wait at least 3-5 minutes between measurements. You may need to adjust the wait time according to your personal physiological situation.
- Wait 30-45 minutes before measurement if you've just consumed caffeinated beverages or smoked cigarettes.
- The applied part is the cuff.
- The patient is the intended operator.

Blood pressure measurements determined with this device are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultatory method, within the accuracy limits prescribed by the American National Standard, Manual, electronic, or automated sphygmomanometers.

In the event that the device needs to be checked for calibration, contact the distributor.

What is blood pressure?
Blood pressure is the measurement of the force of blood pushing against the walls of the arteries. Arterial blood pressure is constantly fluctuating during the course of the cardiac cycle. The pressure measures when the heart contracts and sends blood out of the heart is systolic (highest) blood pressure. The pressure measured when the heart dilates with blood flowing back into the heart is called diastolic (lowest) blood pressure.

Why measure your blood pressure?
Among today's various health problems, those associated with high blood pressure are very common. High blood pressure dangerously correlates with cardiovascular diseases. Therefore, blood pressure monitoring is important for identifying those at risk.

Why do my readings vary?
Blood pressure is a body parameter that is subject to normal variations throughout the day. A single reading that is different from yours or your doctor's readings are not necessarily inaccurate. The average of several readings, taken under similar conditions, using the same arm is preferred for accurate blood pressure readings. Medical physicians generally recommend the "Rule of 3", where you are encouraged to take your blood pressure three times in a row (at 3-5 minute intervals), three times a day for three days. After three days you can average all the results and this will give you an accurate idea of what your blood pressure really is.

Why are my readings different than those at my doctor's office?
Many experience a phenomenon called "White Coat Hypertension" when measured by a doctor. White Coat Hypertension refers to blood pressure that rises above its usual level when measured in a clinical setting, such as a doctor's office.

BLOOD PRESSURE STANDARD

The table below contains defined levels for hypertension that are publicly available from the American Heart Association* (AHA 2017). Users can compare their own blood pressure against these defined levels to determine if they may be potentially at increased risk.

Blood Pressure Category	Systolic mm Hg (upper number)	and	Diastolic mm Hg (lower number)
Normal	<120	and	<80
Elevated	120-129	and	<80
High Blood Pressure (hypertension) Stage 1	130-139	or	80-89
High Blood Pressure (hypertension) Stage 2	≥140	or	≥90
Hypertension Crisis (consult your doctor immediately)	>180	and/or	>120


Source: American Heart Association (AHA) 2017

Blood pressure tends to go up and down, even in people who normally don't have high readings. If your numbers stay above the "normal" range most of the time, you may be at increased risk and should consult your physician. Although one can easily find where their own blood pressure readings fall on this table, this monitor comes equipped with a Risk Category Index that automatically compares each reading to the defined levels and provides a helpful cue if your reading falls into one of the stages that could potentially indicate increased risk. See Risk Category Index section for more information on this feature. Please note that cues provided by this monitor are only intended to assist you in using this table. The table and cues are only provided for convenience to help you understand your non-invasive blood pressure reading as it relates to the American Heart Association* (AHA 2017) information. They are not a substitute for a medical examination or diagnosis by your physician. It is important for you to consult with your physician regularly. Your physician will tell you your normal blood pressure range as well as the point at which you may actually be considered to be at risk.


HOW THIS BLOOD PRESSURE MONITOR WORKS

This monitor uses Smart Measure® inflation technology to detect your blood pressure. Once the cuff is wrapped around your upper arm, with a press of the START/STOP button, the monitor will turn on and the inflation will automatically start, creating pressure around the arteries inside the upper arm. Within the cuff is a gauge which senses the fluctuations (oscillations) in pressure. The fluctuation measured represents the degree of intensity that your arteries are contracting with each heartbeat, which is also a result of the pressure that the cuff has placed on the upper arm. The monitor measures these contractions and converts the information to a digital value. Once the measurement is complete, the cuff will automatically deflate, and the result is displayed on the screen. Please note that any muscle movement during inflation will cause measurement error. When measurement is complete, the monitor will display your systolic pressure, diastolic pressure, and pulse readings. Before measurement, it is suggested that you sit quietly for 15 minutes before measurement as measurements taken in a relaxed state have greater accuracy. The monitor automatically finds where your measurement results fall on the American Heart Association* (AHA 2017) table and provides a cue if your reading falls into one of the stages that could potentially indicate increased risk. Please refer to the Risk Category Index section for more information on this feature. The appearance of the Irregular Heartbeat icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Please refer to the Irregular Heartbeat Detector section for more information on this feature.

NAME/FUNCTION OF EACH PART

DISPLAY SYMBOL EXPLANATIONS

Display Explanation Symbols	User 1: Appears when the monitor is operated by User 1.
no. 1	User 1: Appears when the monitor is operated by User 1.
no. 2	User 2: Appears when the monitor is operated by User 2.
Low Battery Symbol	Appears when batteries should be replaced.
Pulse Symbol	Once pulse is detected, the symbol flashes with each pulse beat. Shows the heart rate per minute.
Risk Category Index	The arrow points out the specific Risk Category that your measurement reading falls in. See Risk Category Index Section for more information.
Irregular Heartbeat Detector	This symbol appears when an irregular heart beat was detected. See Irregular Heartbeat Detector section for more information.
Excessive Body Motion Detector	Displayed if body movement is detected during measurement, especially, the movement on the arm the blood pressure monitor is worn on. If the cuff is worn improperly, or the shape of the upper arm is unusual (for example, the circumference of the upper arm differs largely from the circumference of the forearm), excessive gap might exist between the arm cuff and the arm. NOTE: The measured blood pressure reading may not be accurate if the icon is displayed.
Bluetooth Symbol	Appears when Bluetooth is active, and when data is transmitting to your mobile device.
AVG. 3	The Average of the Last 3 Measurements: This symbol appears when the LCD displays average value of last 3 readings.
5-Day Trend Manage Function Bar Graph	This symbol appears when a measurement is completed or a memory is recalled. It can plot the trend of measurement records in five days. Notice: Please set the correct Date/Time and User-Number before each measurement.

If any of the following letters and numbers appear in the area that systolic pressure should be displayed, an error has occurred with your reading. See Troubleshooting section of this manual for more information.

EE	Measurement Error: Wrap the cuff correctly and keep arm steady during measurement. Take new measurement.
E1	Air Circuit Abnormality: Check cuff connection. Take new measurement.
E2	Pressure Exceeding 300 mmHg: Turn the unit off to clear, then take new measurement.
E3	Error Determining Measurement Data: Ensure the cuff is worn properly, do not move or talk. Take new measurement.
E4	Data Transmission Error: Will occur after 45 seconds, if the monitor cannot connect to the mobile device to transmit data. Make sure Bluetooth* is ON.
EP	System Error: Contact Consumer Relations.

INSTALLING BATTERIES

- Slide the battery cover and lift the cover to open the battery compartment.
- Install or replace 4 AAA (1.5V, LR03) sized alkaline batteries in the battery compartment according to the indications inside the compartment.
- Close the battery cover by sliding the cover back into place.

Replace the batteries if:


- The low battery symbol appears.
- When any button is pressed and nothing is displayed on the screen.

- NOTE:**
- If the batteries are removed or replaced, the date and time will need to be reset either manually or automatically using your Bluetooth* mobile device.
 - Replace all batteries at one time (as simultaneous set). Use only 1.5V AAA alkaline batteries.
 - DO NOT mix alkaline, standard (carbon-zinc) or rechargeable (cadmium) batteries.
 - When the batteries are removed, the measurement values stored in memory are retained. However, the date and time must be reset.
 - Remove batteries when unit is not in use for extended periods of time.
 - Clean contacts on battery and in battery compartment with a soft dry cloth each time you install batteries.
 - Batteries are hazardous waste. Do not dispose of them together with the household garbage.
 - DO NOT dispose of batteries in fire. Batteries may explode or leak.
 - Recycle or dispose of properly in accordance with local, state, province, and country regulations.

USING THE USB TO DC POWER CABLE

This monitor is designed for operating with batteries, a USB to DC cable, or an AC/DC adapter (not included).

- Plug the DC cable connector into the main unit's power supply inlet as shown below:

- Plug the USB cable connector into a USB Wall adapter, then plug the Wall adapter into an outlet, as shown below:

Use only a compatible USB/DC or AC/DC adapter with required voltage and current as indicated in this manual.

Recommend Adapter specification, do not use otherwise:
Model: FranMar International, FRM06-S05-UU
Rating:
Input: 100-240V, 50/60 Hz, 0.2 A
Output: 5V, DC, 1.2A

NOTE:

- When the USB to DC cable, or AC adapter is your main power supply, make sure the adapter plug can be easily removed from the unit.
- Please unload the batteries when operating with a USB to DC cable, or AC adapter for an extended period of time. Leaving the batteries in the compartment for a long time may cause leakage, which may lead to damage to the unit.
- No batteries are needed when operating with a USB to DC cable, or AC adapter.
- Date and time will need to be reset if USB to DC cable, or AC adapter is unplugged even if the unit has the batteries in the battery compartment.

DATE & TIME SET PROCEDURE

Date & time can be set by two methods, either sync automatically using your Bluetooth* mobile device, or set manually using the User/Date & Time set and Memory/Date & Time Select buttons on the monitor

If setting the date & time by syncing with your mobile device, it is important that this is done prior to taking any measurements to ensure the date & time are accurate.

To set using your device with Bluetooth* wireless technology:

- Before taking the first measurement, open the Homedics Health+ app on your device. If you have not already installed the Homedics Health+ app onto your device, it is available to download on the App Store* and on Google Play* Make sure you have the app downloaded and open on your mobile device before trying to sync the date & time.
- Press the Memory/Date & Time Select button for 3 seconds on the monitor to enable the connection.
- Once the Bluetooth* connection is established, the date & time will automatically update on the blood pressure monitor.

NOTE:

- Once the date & time have been successfully synced, future readings will automatically have the correct date & time.
- If you are having trouble automatically syncing the date & time, make sure the Bluetooth* is ON on your device and ON on the blood pressure monitor.

To Set Manually:

- To set the date and time, the monitor needs to be in standby mode. Press START/STOP button and the monitor will automatically turn to standby mode. Then press User/Date & Time set button for 3 seconds and hold.
- The display will show a blinking number showing the YEAR. Change the YEAR by pressing the Memory/Date & Time Select button. Each press will increase the number by one in a cycling manner. Press the User/Date & Time set button again to confirm the entry, and the screen will show a blinking number representing the MONTH.
- Change the DAY, HOUR and MINUTE as described in Step 2 above, using the Memory/Date & Time Select button to change the numbers and the User/Date & Time set button to confirm the entries.


- The date and time will only need to be set manually if the monitor will not be used with the Homedics Health+ app.
- If the Homedics Health+ app is used after manually setting the date and time on the blood pressure monitor, the date and time on the mobile device will override the date and time on the blood pressure monitor.

BLUETOOTH* OPERATION

This monitor has Bluetooth function ON for your convenience. This will allow your readings to automatically transmit to the Homedics Health+ app. By turning this function OFF, the measurements cannot be transmitted.

Turning Bluetooth* Function OFF:

- The default setting for this function is ON.
- Once the date and time has been set (please refer to the "Setting Year, Date and Time" section), press the User/Date & Time set button, then press the Memory/Date & Time Select button to turn the Bluetooth* ON or OFF.

App for Bluetooth

Download and install the Homedics Health+ app on your smart device from Google Play or App store.

- System requirement of the Bluetooth-enabled device
- Bluetooth 5.0 for Android 6.0 or above
 - Bluetooth 5.0 for iOS 9.0 or above

5-DAY TREND MANAGE FUNCTION

The 5-Day Trend Manage Function can help users to acquire a better knowledge of their "5-day" trend of blood pressure variation. It averages all blood pressures of individual day by weighting the Systolic and Diastolic pressures and compares the blood pressure readings of each individual day with each other; it then illustrates the comparative result using a bar graph, depicting blood pressure changing in five successive days.


The bar graph updates when a new measurement completes or when users recall previous records in memory mode. It will only show up when there exists at least 2 days of measuring records among 5 successive days.

- NOTE:**
- The function relies on correct Date/Time and User Number settings of each blood pressure measurement. Please check Date/Time and User Number before taking a measurement.
 - Do not take the measurement result as a replacement for medical diagnosis. Please consult your doctor if you have any concerns regarding the measurement results.


Detailed explanation of the screens:

- There are 5 columns in the bottom of LCD. Each column represents a day.
- "Latest-day" always corresponds with the date of reading that shows on the screen. From "D-1" to "D-4" "D-1" is for one day before the "Latest-day", and so forth.
- Each column has 4 levels. It compares each day's averaged blood pressure with each other and subdivides the largest difference into 4 levels. The upmost level represents the highest blood pressure among 5 days, and the ground level for the lowest. Each day's bar represents the level of that day.


Suppose that the user has taken 7 measurements from 11/1 to 11/6, with two on 11/6 and one on each prior days. When the user had just completed the most recent measurement, the screen would be as follows:

If the user then pressed the Memory/Date & Time Select button and scrolled to the most recent reading in memory mode, the screen would be the same as that one just measured. If user pressed Memory/Date & Time Select button again to view prior reading:

If user pressed Memory/Date & Time Select button again to view prior reading:


USING THE UPPER ARM CUFF


Very Important: Proper cuff size is critical to accurate measurement.

This monitor comes with a One Size Fits Most Upper Arm Cuff that fits arm sizes 9in-17in (23cm - 43cm). The appropriate cuff is suitable for your use if the colored arrow falls within the solid color line as shown below.


If the arrow falls outside of the solid color line, you will need an XL arm cuff model # BPA-CUFFXL A (17in - 22in, 43cm-56cm). If you have any additional questions, or need help on purchasing, please contact Homedics Consumer Relations, please see warranty section.

- If the cuff is not assembled, pass the end of the cuff furthest away from the tubing through the metal D-ring in order to form a loop. The smooth side without the felt material should be on the inside of the cuff loop.

- Plug the cuff tube into the left side of the unit.


- Open the arm cuff so that your arm may be placed through it.
- Wrap the cuff on a bare arm or over thin clothing. Thick clothing or a rolled up sleeve will cause inaccurate blood pressure measurements.
- Position cuff edge 0.8 - 1.2 inches (2-3 cm) above elbow.
- Center tube over middle of arm.

- Pull the end of the cuff so that it tightens evenly around your arm. Press the thick end loop material together to secure. Allow room for 2 fingers to fit between the cuff and your arm. Please make sure the cuff does not slip during measurement, and the arrow falls within the proper fit range.
- Lay your arm on a table (palm upward) so the cuff is at the same height as your heart. Make sure the tube is not kinked.

Note:

- This device should not be used when your wrist is wounded or injured.
- If it is not possible to fit the cuff on your left arm, it can also be placed on your right arm. However, all measurements should be made using the same arm.
- To use on the right arm, you must position the artery symbol over the main artery. Locate the main artery by pressing with two fingers approximately 1in (2 cm) above the bend of your elbow on the inside of your right arm. Identify where the pulse can be felt the strongest. This is your main artery.


MEASUREMENT PROCEDURE

- Important Note:**
- Blood pressure changes with every heartbeat and is in constant fluctuation throughout the day.
 - Blood pressure measurement can be affected by the position of the user, his or her physiologic condition and other factors.
- Before Measurement:**
- To help ensure an accurate reading, follow these instructions before taking a measurement:
- Wait 1 hour after exercising, bathing, eating, drinking beverages with alcohol or caffeine, or smoking to measure blood pressure.
 - Sit quietly and rest for 15 minutes.
 - Stress raises blood pressure. Avoid taking measurements during stressful times.
 - Take your reading in a comfortable environment as measurements can be affected by hot or cold temperatures. Take your blood pressure at normal body temperature.

During Measurement:

- DO NOT talk or move your arm or hand muscles.
- Do not cross your legs. Sit with feet flat on the floor.
- DO NOT touch cuff or monitor during measurement.

If you are using this blood pressure monitor for the first time, please remove the protective film from the screen.

- With the cuff wrapped around your arm, press the START/STOP button to confirm the chosen user and start measurement.
- All display symbols appear on the screen for 1.5 seconds.
- After all symbols appear, the display will show "00". The monitor is "Ready to Measure" and will automatically inflate to the level that is right for you.

- As the cuff inflates, the monitor automatically determines your ideal inflation level. This monitor features Smart Measure® inflation technology, and detects your blood pressure and pulse rate during inflation. The Heartbeat Symbol flashes with every heartbeat. Remain still and do not move until the entire measurement process is completed. The device will detect your pulse and determine the blood pressure.

Note!

If the cuff does not stop inflating, remove the cuff at once. To stop measurement, press START/STOP button. The cuff will deflate immediately after the button is pressed.

- After the monitor has determined your blood pressure and heart rate, the cuff automatically deflates. Your results including systolic pressure, diastolic pressure, heart rate, corresponding Risk Category Indicator, graphic of 5 Day Trend Function, Irregular Heartbeat Detector and Excessive Body Motion Detector (if any) will be displayed on the screen and the cuff will deflate to memory automatically.
- Device automatically shuts off if no operation over 1 minute.

TRANSMIT READINGS TO YOUR MOBILE DEVICE

This monitor syncs your blood pressure readings to the Homedics Health+ app that is free to download on the App Store* and on Google Play*. Make sure you have the Homedics Health+ app downloaded, and open on your mobile device before trying to transmit your blood pressure measurements.

Before attempting to sync the blood pressure monitor with your mobile device, make sure Bluetooth* is turned ON, on both your mobile device and the monitor.

Transmit Reading

There are 2 ways to activate the Bluetooth* function.

- Automatically activate:** When measurement completed, the device activates the Bluetooth* function automatically, and Bluetooth* symbol will flash on the screen.
- Manually activate:** Under standby mode, user can press the button for 3 seconds to activate the Bluetooth* function, and the Bluetooth symbol will flash on the screen.

RISK CATEGORY INDEX

This monitor comes equipped with a Risk Category Index that automatically compares each reading to defined levels established by the American Heart Association* (AHA 2017) as described earlier in this manual, and provides a helpful cue if your reading falls into one of the stages that could potentially indicate increased risk. Please note that cues provided by this monitor are only intended to assist you in using this table. The table and cues are only provided for convenience to help you understand your non-invasive blood pressure reading as it relates to the AHA information. They are not a substitute for a medical examination or diagnosis by your physician. It is important for you to consult with your physician regularly. Your physician will tell you your normal blood pressure range as well as the point at which you may actually be considered to be at risk.

Blood Pressure Category	Systolic mm Hg (upper number)	and	Diastolic mm Hg (lower number)	Indicator Color
Normal	<120	and	<80	Green
Elevated	120-129	and	<80	Yellow
High Blood Pressure (hypertension) Stage 1	130-139	or	80-89	
High Blood Pressure (hypertension) Stage 2	≥140	or	≥90	Red
Hypertension Crisis (consult your doctor immediately)	>180	and/or	>120	

Source: American Heart Association (AHA) 2017

IRREGULAR HEARTBEAT DETECTOR (IHB)

The appearance of the Irregular Heartbeat Detector icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Usually this is not a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. Please note that the device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

Movement, shaking or talking during the measurement can result in pulse irregularities which may cause the appearance of this icon. Therefore, it is of great importance to not move or talk during measurement.

To determine the presence of an irregular heartbeat, the average of the heartbeat intervals is calculated with the first 3 normal effective heartbeats. It is important to note that the average is not a strict mathematical averaging of all recorded intervals. At least 3 beats with 25% or greater difference from the average heartbeat interval will generate the Irregular Heartbeat Detector icon on the screen.

Important Information:


This blood pressure monitor is not designed for use by people with arrhythmias nor for diagnosing or treating an arrhythmia problem. As a safeguard, we recommend that if you have arrhythmias such as atrial or ventricular premature beats and atrial fibrillation or any other special conditions you should check with your physician before using your blood pressure monitor.


USING THE MEMORY FUNCTION

After each measurement, the results including systolic and diastolic pressure, heart rate and corresponding Risk Category indicator, graphic of 5-Day Trend Function, Irregular Heartbeat Detector and Excessive Body Motion Detector (if any) with the time and date will be automatically stored. The monitor features 2 user memory capabilities. The monitor holds the last 120 measurements for each user, and automatically replaces the oldest data with a new one.


RECALLING VALUES FROM MEMORY

- Press the User/Date & Time set  button to select User 1 or User 2.

- Press the Memory/Date & Time Select  button to access memory. If there is no data stored before, nothing (except month, date and time) will appear on the display. If yes, the first reading will be the average of the last 3 measurements.

- Press the Memory/Date & Time Select  button repeatedly to scroll through the previous readings stored in the memory. The 5-Day Trend bar graph will plot the "5-Day" trend of blood pressure variation, spanning 4 days before the date of the reading that shows on the screen.




Note:

For viewing the 5-Day Trend bar graph:

- The rightmost bar, "Latest day", will always correspond with the date of reading that shows on the screen.
- For more detailed explanation about the 5-Day Trend Manage function, please see the "5-Day Trend" section.

- To stop reading the memories, press the START/STOP button to switch to Standby Mode.

CLEARING VALUES FROM MEMORY

- Press the User/Date & Time set  button to select User 1 or User 2.
- Press the Memory/Date & Time Select  button to access the memory recall mode.
- Press and hold User/Date & Time set  and Memory/Date & Time Select  buttons at the same time and the data for the selected user will be erased automatically.

Note: Once deleted, your readings cannot be restored.

IMPORTANT NOTES REGARDING YOUR BLOOD PRESSURE MEASUREMENT

- It is suggested that you take your measurements at the same time each day and use the same arm consistently.
- Users should wait a minimum of 5 minutes before taking additional measurements. More time may be necessary depending upon your physiology.
- The measurement results that users receive are for reference only. If users have any blood pressure concerns, please consult a physician.
- Once inflation reaches 300 mmHg, the unit will deflate automatically for safety reasons.
- This device is not suitable for people with arrhythmias.
- This device may have difficulty determining the proper blood pressure for users with irregular heartbeat, diabetes, poor circulation of blood, kidney problems, or for users who have suffered a stroke.

CARE, MAINTENANCE, & CLEANING

- Clean the blood pressure monitor body and cuff carefully with a slightly damp, soft cloth. DO NOT press. Do not wash cuff or use chemical cleaner on it. Never use thinner, alcohol or petrol (gasoline) as cleaner.
- Make sure the cuff is completely dry before using.
- Leaky batteries can damage the unit. Remove the batteries when the unit will not be used for a long time.
- Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.
- If the unit is stored near freezing, allow it to acclimate to room temperature before use.
- This blood pressure monitor is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems with this device, please contact Homedics Consumer Relations (contact information on warranty page).
- DO NOT immerse the unit in water as this will result in damage to the unit.
- DO NOT subject the monitor or cuff to extreme temperatures, humidity, moisture, or direct sunlight. Protect from dust.
- DO NOT fold the cuff tightly.
- DO NOT disassemble the monitor or cuff. If in need of repair, refer to the warranty section of this manual.
- DO NOT subject the monitor to extreme shocks (DO NOT drop on floor).
- DO NOT twist the cuff in any way.
- DO NOT inflate the cuff unless wrapped around arm.
- DO NOT wrap the cuff around body parts other than your arm.
- DO NOT drop or insert any object into any opening or hose.
- To avoid accidental strangulation keep this product away from children. Do not drape tube around neck.
- This monitor may not meet its performance specifications if stored or used outside of these temperature and humidity ranges:

POTENTIAL FOR ELECTROMAGNETIC INTERFERENCE

To avoid inaccurate results caused by electromagnetic interference between electrical and electronic equipment, Do not use the device near a cell phone or microwave oven. For most wireless communication devices, it is recommended to maintain a distance of 10-8 feet (3.3m) in order to avoid electromagnetic interference.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

FEDERAL COMMUNICATIONS COMMISSION COMPLIANCE STATEMENT

Changes or modifications to this equipment not expressly approved by the manufacturer could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antennas.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

RF exposure warning

- The equipment complies with FCC RF exposure limits set forth for an uncontrolled environment.
 - The equipment must not be co-located or operation in conjunction with any other antenna or transmitter.
- FCC Label Compliance Statement:
This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:
(1) This device may not cause harmful interference, and
(2) This device must accept any interference received, including interference that may cause undesired operation.
Note
Changes or modifications not expressly approved by the manufacturer could void the user's authority to operate the equipment*.

ELECTROMAGNETIC COMPATIBILITY (EMC)

Guidance and manufacturer's declaration – electromagnetic emissions
The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF Emissions CISPR 11	Group 1	RF energy is used only to maintain device's operation. Therefore, its RF emissions are so low that it's not likely to cause any interference in nearby.
RF emissions CISPR 11	Class B	
Harmonic Emissions IEC 61000-3-2	Class A	The device is suitable for use in all establishments, including domestic establishments, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity
The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact discharge ±15 kV air discharge	± 8 kV contact ± 15 kV air discharge	In the case of air discharge testing, the climatic conditions shall be within the following ranges: Ambient Temperature: 15°C – 35°C Relative Humidity: 30% – 60%.
Power frequency (50 or 60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 or 60 Hz	30A/m 50 or 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Electrical fast transient/burst IEC 61003-4	±2 kV for power supply lines ±kV for input/output lines	±2 kV for power supply lines ±kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	AC Power port ±1 kV Line to Line	AC Power port ±1 kV Line to Line	Mains power quality should be that of a typical commercial or hospital environment.
Interruptions and voltage variations on power supply input lines IEC 61003-4-11	OR6 UT: 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°.	OR6 UT: 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°.	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.

Guidance and manufacturer's declaration – electromagnetic immunity
The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Conducted RF IEC 61000-4-6	3V rms At 0.15-80 MHz 6V rms At ISM & Radio Amateur Freq	Not Applicable	Portable and mobile RF communication equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3 Proximity fields from RF wireless communications equipment IEC 61000-4-3	10 V/m at 80-2700 MHz AM Modulation And 9-28V/m at 385-6000MHz, Pulse Mode and other Modulation. The system shall be tested as specified in IEC60601-1-2 Table 9 for proximity fields from RF wireless communications equipment using the test methods specified in IEC 61000-4-3	10 V/m at 80-2700 MHz AM Modulation And 9-28V/m at 385-6000MHz, Pulse Mode and other Modulation. The system shall be tested as specified in IEC60601-1-2 Table 9 for proximity fields from RF wireless communications equipment using the test methods specified in IEC 61000-4-3	Recommended separation distance Considering to reduce the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation: $E=6/d^{0.7}$ where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVELS in V/m. Field strengths from fixed RF transmitters, as determined by an electromagnetic site should be less than the compliance level in each frequency Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a: Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast can-not be predicted theoretically with accuracy. To assess the electro-magnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be ob-served to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.
- b: Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Test specifications for enclosure port immunity to RF wireless communications equipment.

Test frequency (MHz)	Modulation	Immunity Test Level (V/m)
385	Pulse modulation 18 Hz	27
450	FM± 5 kHz deviation 1kHz sine	28
710		
745	Pulse modulation 217 Hz ^{±1}	9
780		
810		
870	Pulse modulation 18 Hz ^{±1}	28
930		
1720		
1845	Pulse modulation 217 Hz ^{±1}	28
1970		
2450	Pulse modulation 217 Hz ^{±1}	28
5240		
5500	Pulse modulation 217 Hz ^{±1}	9
5785		

- a). The carrier shall be modulated using a 50% duty cycle square wave signal.
b). AS an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

TROUBLESHOOTING

If any abnormality arises during use, please check the following points:

Symptoms	Possible Causes	Correction
Unit does not turn on when the START/STOP button is pushed.	Batteries have run down. Battery polarities (+/-) have been positioned incorrectly.	Replace them with four new "AAA" alkaline batteries. Re-insert the batteries in the correct positions.
EE Measurement error symbol shown on display or the blood pressure value is displayed excessively low (or high).	The Upper Arm Cuff has been positioned on the arm incorrectly. Did you talk or move during measurement? Shaking of the arm with the cuff on.	Re-wrap the cuff properly so that it is positioned correctly. Take new measurement. Keep arm steady during measurement. Measure again. Refer to "Measurement Procedure" instructions.
E1 Error symbol shown on display.	Air circuit abnormality. Cuff tube may not be plugged into monitor correctly.	Check cuff connection. See Using the Upper Arm Cuff section. Take new measurement.
E2 Error symbol shown on display.	Inflation pressure exceeding 300 mmHg.	Switch the unit off and then turn back ON. Take new measurement.
E3 Error symbol shown on display.	Error determining measurement data.	Re-wrap the cuff properly so that it is positioned correctly. Take new measurement.
E4 Data transmission error symbol shown on the display.	Bluetooth® function is turned OFF, on your mobile device. Bluetooth® function is turned OFF, on the blood pressure monitor. Mobile device does not pair with the blood pressure monitor.	Turn Bluetooth® ON, on your mobile device, and try again. Turn the Bluetooth® function ON, on your blood pressure monitor, then try again. Make sure Bluetooth® is ON, on the mobile device and the blood pressure monitor, and try again.
	Requires compatible mobile device. This blood pressure monitor is compatible with mobile devices that are running Bluetooth® 5.0 for Android 6.0 or above, Bluetooth® 5.0 for iOS 7.0 or above.	
	The App on the mobile device is not open.	Make sure the app is open on the mobile device, and try again.
	The blood pressure monitor and mobile device are out of transmitting range.	Make sure the mobile device and blood pressure monitor are within the acceptable range of 32 feet (10 meters).
	Unexpected loss of electrical/mechanical integrity.	Make sure your last reading is stored in memory and the App is open and try again. See Transmitting Readings to Your Mobile Device section. Remove batteries, re-insert, and try again. Return the device to your local distributor or importer.



 Excessive Body Motion Detector	Body movement during measurement, especially, the movement on the arm the blood pressure monitor is worn on. e.g. Talking, moving or shaking of the arm with the cuff on while measurement.	Measure again. Keep arm steady during measurement.
	Cuff is worn improperly, or the shape of the upper arm is unusual (for example, the circumference of the upper arm differs largely from the circumference of the forearm), excessive gap might be exist between the arm cuff and the arm.	Wrap the cuff properly and keep steady. Measure again. If you have any question about the cuff wearing and/or measurement result, please consult your healthcare professional.
	The distance between BPM and Bluetooth-enabled device is out of transmitting range.	Please make sure the acceptable distance (±510 meters) with each other.
	Use an incompatible Bluetooth-enabled device, Use non-Bluetooth-enabled device,	Please refer to Bluetooth Operation and RF specification section.
BPM cannot communicate with Bluetooth-enabled device	Paring has not been completed. Bluetooth® function is not turned on.	Please re-pair the BPM and Bluetooth-enabled device with each other. See the "Bluetooth® Transmission" section to turn on Bluetooth® function.
EP error symbol shown on display.	System Error.	Measure again. If error persists, contact Consumer Relations.
	Unexpected loss of electrical/mechanical integrity,	Re-insert the batteries and try again. Return the device to your local distributor or importer.

Note:
If the unit still does not work, contact Homedics Consumer Relations. Under no circumstance should you disassemble or attempt to repair the unit by yourself. Contact information for Homedics Consumer Relations Department can be found on the warranty page.

SPECIFICATIONS

Model Number	HL858A2-BL/BPA-100TRBT
Measurement Method	Oscillometric (Inflation)
Rated Range of Cuff Pressure	0 – 300 mmHg
Rated Range of Determination	40 – 280 mmHg
Measurement Range of Heart Rate	40 ~ 199 beats/minute
Accuracy	Pressure: ± 3 mmHg Pulse: ± 5 % Max.
Inflation	Automatic Inflation (Air Pump)
Deflation	Automatic Air Release Control Valve
Display	Liquid Crystal Display
Memory	240 Memory Total for 2 Users
Unit Dimensions	5.51 x 4.33 x 2.22 inch (LxWxH) 140 x 110 x 56.5 mm (LxWxH)
Unit Weight	276 ± 5 g (9.74 ± 0.18 oz) (Cuff and batteries excluded)
Cuff Size	UC-01: Universal size cuff 9 – 17 inch (23 – 43 cm)
Storage/Transportation Environment	Temperature: -25°C ~ 70°C (-139°F ~ 158°F) Humidity: ≤ 93% R.H.
Operation Environment	Temperature: 5°C ~ 40°C (41°F ~ 104°F) Humidity: 15% ~ 93% R.H. Atmospheric pressure: 700 hPa ~ 1060 hPa
Power Supply	DC 6V AAA "LR03" (1.5V) alkaline Battery x 4 USB to DC Jack cable DC 5V 1A AC/DC adapter (Optional)
Battery Life	Approx. 250 Measurements
Shelf life (battery)	3 years (Temperature: 20 ± 2°C; Relative humidity: 65 ± 20% RH)
Product Life	5 Years (4 times per day)
Sleeping Mode	Without any operation for 1 minute, device automatically shuts off.
Accessories	4 AAA (LR03) Alkaline Batteries, Arm Cuff with Tube, Instruction Manual, Pouch
RF Type	Bluetooth® 5.0 BLE
System requirement of the Bluetooth®-enabled device	Bluetooth® 5.0 for Android 6.0 or above Bluetooth® 5.0 for iOS 9.0 or above

* The contents of this manual and the specifications of the device covered by this manual are subject to change for improvement without notice.

	Follow instruction for use.
	TYPE BF Applied Part.
	To avoid inaccurate results caused by electromagnetic interference. Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30cm (12 inches) to any part of the device. Otherwise, degradation of the performance of this equipment could result.
	Waste of electrical and electronic equipment (WEEE).
	Discard the used product to the recycling collection point according to local regulations.
	Serial number SN YYYMMXXXXXX
IP22	Ingress Protection Rating. First characteristic numeral - Degree of protection against access to hazardous parts and against solid foreign objects N1=2 (Protected against solid foreign objects of 12.5 mm Ø g and greater).
	Second characteristic numeral - Degree of protection against ingress of water. N2=2 (Protected against vertically falling water drops when ENCLOSURE tilted up to 15°).
	Non-ionizing electromagnetic radiation.

* The contents of this manual and the specifications of the device covered by this manual are subject to change for improvement without notice.
Note: These specifications are subject to change without notice.
Device information:
• Internally powered equipment
• Not suitable for use in presence of flammable anesthetic mixture with air or with Oxygen or nitrous oxide
• Continuous operation with short-time loading

Innovation, Science and Economic Development Canada ICES-003 Compliance Label: CAN ICES-3 (B)/NMB-3(B)

Homedics

For service or repair, do not return this unit to the retailer. Contact Homedics Consumer Relations at:

Email: cservice@homedics.com

Phone: 1-800-466-3342
Business Hours: 8:30am-7pm EST
Monday-Friday

Distributed by: Homedics, LLC
3000 N Pontiac Trail
Commerce Township, MI 48390
Printed in China or Vietnam

LIMITED 5-YEAR WARRANTY

Homedics sells its products with the intent that they are free of defects in manufacture and workmanship for a period of 5 years from the date of original purchase, except as noted below. Homedics warrants that its products will be free of defects in material and workmanship under normal use and service. This blood pressure monitor meets the simulated measurement cycles test requirement per EN1060-3, part 8.10. This warranty extends only to consumers and does not extend to Retailers.

To obtain warranty service on your Homedics product, contact a Consumer Relations Representative. Please make sure to have the model number of the product available.

Homedics does not authorize anyone, including, but not limited to, Retailers, the subsequent consumer purchaser of the product from a Retailer or remote purchasers, to obligate Homedics in any way beyond the terms set forth herein. This warranty does not cover damage caused by misuse or abuse; accident; the attachment of any unauthorized accessory; alteration to the product; improper installation; unauthorized repairs or modifications; improper use of electrical/power supply; loss of power; dropped product; malfunction or damage of an operating part from failure to provide manufacturer's recommended maintenance; transportation damage; theft; neglect; vandalism; or environmental conditions; loss of use during the period the product is at a repair facility or otherwise awaiting parts or repair; or any other conditions whatsoever that are beyond the control of Homedics.

This warranty is effective only if the product is purchased and operated in the country in which the product is purchased. A product that requires modifications or adoption to enable it to operate in any other country than the country for which it was designed, manufactured, approved and/or authorized, or repair of product damaged by these modifications is not covered under this warranty.

THE WARRANTY PROVIDED HEREIN SHALL BE THE SOLE AND EXCLUSIVE WARRANTY. THERE SHALL BE NO OTHER WARRANTIES EXPRESS OR IMPLIED INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS OR ANY OTHER GUARANTEE ON THE PART OF THE COMPANY WITH RESPECT TO PRODUCTS COVERED BY THIS WARRANTY. HOMEDICS SHALL HAVE NO LIABILITY FOR ANY INCIDENTAL, CONSEQUENTIAL OR SPECIAL DAMAGES. IN NO EVENT SHALL THIS WARRANTY REQUIRE MORE THAN THE REPAIR OR REPLACEMENT OF ANY PART OR PARTS WHICH ARE FOUND TO BE DEFECTIVE WITHIN THE EFFECTIVE PERIOD OF THE WARRANTY. NO REFUNDS WILL BE GIVEN. IF REPLACEMENT PARTS FOR DEFECTIVE MATERIALS ARE NOT AVAILABLE, HOMEDICS RESERVES THE RIGHT TO MAKE PRODUCT SUBSTITUTIONS IN LIEU OF REPAIR OR REPLACEMENT.

This warranty does not extend to the purchase of opened, used, repaired, repackaged and/or resaled products, including but not limited to sale of such products on internet auction sites and/or sales of such products by surplus or bulk resellers. Any and all warranties or guarantees shall immediately cease and terminate as to any products or parts thereof which are repaired, replaced, altered, or modified, without the prior express and written consent of Homedics. This warranty provides you with specific legal rights. You may have additional rights which may vary from state to state. Because of individual state regulations, some of the above limitations and exclusions may not apply to you.

For more information regarding our product line in the USA, please visit: www.homedics.com

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of these marks by Homedics is under license. Other trademarks and trade names are those of their respective owners. App Store is a service mark of Apple Inc. Android and Google Play are trademarks of Google, LLC. Homedics and the Smart Measure word mark and logo are registered trademarks of Homedics, LLC. American Heart Association is a registered trademark of the American Heart Association, Inc. ©2022 Homedics, LLC. All rights reserved. IB-BPA100TRBT P/N: 323103675

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The 5-Day Trend Management Function can help users to acquire a better knowledge of their "5-day" trend of blood pressure variation. It averages all blood pressures of individual day by weighting the Systolic and Diastolic pressures and compares the blood pressure readings of each individual day with each other; it then illustrates the comparative result using a bar graph, depicting blood pressure changing in five successive days.

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