# Homedics

## 300 series upper arm blood pressure monitor

BPA-800-WT  
Doc # L-0295, Rev.5

**5 year limited warranty**

### WARNINGS AND PRECAUTIONS

Please read the instruction manual carefully before using product. If you have any questions, please contact Homedics Consumer Relations. Save these instructions for future reference.

### WARNINGS

**This product is not suitable for newborns.**

- This product is for home use only. Consult your doctor with any concerns about the measurement results.
- DO NOT use this device for diagnosis or treatment of any health problem or disease. Measurement results are for reference only. Consult a healthcare professional for interpretation of pressure measurements. Contact your physician if you have or suspect any medical problem. DO NOT change your medications without the advice of your physician or healthcare professional.
- When common arrhythmias (such as atria premature beats, premature ventricular and atrial fibrillation) occurs, it will affect the measurement accuracy of blood pressure. If it is only sporadic, please rest for 2 hours before measuring. If it is frequent it is recommended you seek medical advice.
- If the cuff is inflated for a long time, immediately press the "START/STOP" button to stop the measurement. Continued inflation pressure of the arm can cause limb numbness.
- DO NOT allow children to use without supervision. Small parts may pose a choking hazard.
- Those who have had a mastectomy surgery (especially when lymph nodes have been removed), should take measurements on the unaffected side.

### PRECAUTIONS

- Avoid high temperature, humidity, dust, and direct sunlight.
- Avoid damage to the arm band and rubber tube due to folding.
- DO NOT disassemble the monitor or cuff. If in need of repair, refer to the Warranty section of this manual.
- Avoid dropping or shaking the device.
- The normal measurement period is about 1 minute. If the measurement does not stop after 1 minute, please press the "START/STOP" button to stop the measurement. Avoid numbness of the limbs due to prolonged compression of the arm.
- For better reading of the display, please note the following visible conditions:
  - Ambient brightness: 100lm-1500lm
  - Line of sight: less than 30cm
  - Viewing angle: normal display ± 30°
- Extreme temperature, humidity, and altitude conditions can affect the performance of the measurement, and the sphygmomanometer may not meet the stated performance specifications.
- DO NOT put the cuff on a wounded or injured arm, which will cause further damage.
- DO NOT measure blood pressure if there is intravascular access or treatment, or arteriovenous (A-V) shunt.

### Purpose/Intended Use of Device

Upper Arm Blood Pressure Monitor is intended for over-the-counter home use to measure the blood pressure and pulse rate of adults and children at least 12 years of age. The cuff is intended for use on the left upper arm. Use according to instructions.

**NOTE:** NOT suitable for neonate, pregnancy, or preeclampsia.

### INSTRUCTIONS

#### Product Features

This is a home healthcare product only and is not intended to serve as a substitute for the advice of a physician or medical professional. This device uses the oscillometric method to measure systolic and diastolic blood pressure, as well as heart rate.

#### Self-Measurement

- Please keep in mind: This reading is not meant to be a diagnosis or treatment. If you have abnormal blood pressure readings, consult your physician immediately and follow your physician's instructions.
- The pulse displayed by this unit is not suitable to be a fixed frequency detector that identifies heart rate!
- This device may have difficulty determining proper blood pressure for users with irregular heartbeat, diabetes, liver disease, kidney disease, poor circulation of the blood, or for users who have suffered a stroke. Please consult your healthcare professional before using the device.

#### Electromagnetic Interference

Electromagnetic interference: The device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens). These may lead to temporary impairment of measurement accuracy. Consider the electromagnetic compatibility of the device (ex. power disturbance, radio frequency interference, etc.) Please use this device in a home health care environment only.

### ABOUT BLOOD PRESSURE

#### What is blood pressure?

Blood pressure is the measurement of the force of blood pushing against the walls of the arteries. Arterial blood pressure is constantly fluctuating during the course of the cardiac cycle. The pressure measured when the heart contracts and sends blood out of the heart is systolic (highest) blood pressure. The pressure measured when the heart dilates with blood flowing back into the heart is called diastolic (lowest) blood pressure.

#### Why measure your blood pressure?

Among today's various health problems, those associated with high blood pressure are very common. High blood pressure dangerously correlates with cardiovascular diseases. Therefore, blood pressure monitoring is important for identifying those at risk.

#### Why do my readings vary?

Blood pressure is a body parameter that is subject to normal variations throughout the day. A single reading that is different from yours or your doctor's readings are not necessarily inaccurate. The average of several readings, taken under similar conditions, using the same arm is preferred for accurate blood pressure readings.

#### Why Are My Readings Different Than Those Taken at My Doctor's Office?

Many experience a phenomenon called "White Coat Hypertension" when measured by a doctor. White Coat Hypertension refers to blood pressure that rises above its usual level when measured in a clinical setting, such as a doctor's office.

#### What is normal blood pressure?

If blood pressure value is too high during rest: diastolic pressure (DIA) exceeds 80 mmHg, and either systolic pressure (SYS) is higher than 120 mmHg. In this case, please consult your physician immediately. Prolonged hypertension will damage blood vessels and important organs such as kidneys and even the heart. Consult your physician if systolic blood pressure (SYS) exceeds 130mmHg and diastolic pressure (DIA) exceeds 80 mmHg. It is necessary to use this sphygmomanometer for regular blood pressure self-measurement even though your blood pressure is in the normal scope. Through this method, you can find changes in blood pressure earlier and take appropriate measures. If you are using medication to control your blood pressure, please take measurements and record your blood pressure every day. Share these records with your physician. Do not change the prescription or dosage prescribed by your physician according to your measurement. The following is a standard classification established by the American Heart Association® (AHA 2017) for blood pressure, shown in the below table: Unit: mm(Hg).

Blood Pressure Category	Systolic mm Hg (upper number)	and	Diastolic mm Hg (lower number)	Indicator Color
Normal	<120	and	<80	Green
Elevated	120-129	and	<80	Yellow
High Blood Pressure (hypertension) Stage 1	130-139	or	80-89	Red
High Blood Pressure (hypertension) Stage 2	≥140	or	≥90	
Hypertension Crisis (consult your doctor immediately)	>180	and/or	>120	

\*Source: American Heart Association® (AHA) 2017

Blood pressure tends to go up and down, even in people who normally don't have high readings. If your numbers stay above the "normal" range most of the time, you may be at increased risk and should consult your physician. Although one can easily find where their own blood pressure readings fall on this table, this monitor comes equipped with a Risk Category Index that automatically compares each reading to the defined levels and provides a helpful cue if your reading falls into one of the stages that could potentially indicate increased risk.

**PLEASE NOTE** that cues provided by this monitor are only intended to assist you in using this table. The table and cues are only provided for convenience to help you understand your non-invasive blood pressure reading as it relates to the American Heart Association® (AHA 2017) information. They are not a substitute for a medical examination or diagnosis by your physician. It is important for you to consult with your physician regularly. Your physician will tell you your normal blood pressure range as well as the point at which you may actually be considered to be at risk.

**What measure should be taken when blood pressure is too high or low?**  
**Please consult physician.**

Prolonged ascending blood pressure (different types of SYS) will endanger human health. Deposits on the walls of the blood vessels limit the flow of blood (that can result in arteriosclerosis), which is very dangerous. Because it will cause insufficient blood supply to important body parts (heart, brain, muscles, etc.), it may severely destroy the structure of heart. There are many factors that cause high blood pressure. We could divide them into common hypertension and secondary hypertension. Secondary hypertension will lead to organ disorders. If your blood pressure value continues to rise, consult your physician about possible causes. Changing your lifestyle can also prevent or lower hypertension. These healthy lifestyle habits include:

**Dietary Habit**  
Maintain normal weight as guided by doctor. DO NOT eat too much salt, as many "packaged foods" contain more salt. Avoid eating greasy food. (Packaged foods usually contain large amounts of fat).

**Preventing Diseases**  
Adhere to medical guidelines for preventing certain diseases, such as diabetes, fat metabolic disorders and gout.

**Living Habits**  
Don't smoke, avoid drinking too much or high concentrations of drinking; limit the intake of caffeine (coffee, tea, chocolate, etc.).


**Physical Exercise**  
After a medical examination, exercise regularly; choose an exercise program that requires endurance rather than strength.

**PLEASE NOTE:** DO NOT exercise beyond your physical limit. Patients over 40 years old with a medical history, should consult your physician before you start to exercise.

**PRECAUTION:** Consult your physician before using the device for any of the following conditions: common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation, arterial sclerosis, poor perfusion, diabetes, old age, pregnancy, pre-eclampsia, renal diseases.

**NOTE:** that motion, trembling, and shivering may affect the reading.

### DESCRIPTION OF PRODUCT STRUCTURE

Low Battery Alert

Blood Pressure Level Display

Unit kPa

Average Memory

Memory Storage

Date

Time

Date/Time AM PM

SYS

DIA

Unit mmHG

Pulse

Heartbeat


Irregular pulse wave

### OPERATION MANUAL

#### Install the batteries

After opening the package, first install the batteries. The battery case is located on the back of the monitor. Process of battery installation is as follows:

- Remove battery cover.
- Insert batteries, and make sure the positive, negative poles of the batteries are the same as the battery case's positive, negative poles are as shown.
- If LCD displays the low battery signal, please replace with 4 new batteries.

### PRECAUTIONS

- Once the low battery signal is displayed, the unit will not be available unless you replace new battery.
- Please use 4 standard "AA", long-lasting alkaline batteries.
- Remove battery if the unit is to remain unused for an extended period.

**NOTE:** If the batteries are removed or replaced, the date and time will need to be re-set.

- Replace all batteries at one time (as simultaneous set). Use only 1.5V "AA" alkaline batteries. DO NOT mix alkaline, standard (carbon-zinc) or rechargeable (cadmium) batteries.
- When the batteries are removed, the measurement values stored in memory are retained. However, the date and time must be reset.
- Remove batteries when unit is not in use for extended periods of time.
- Clean contacts on battery and in battery compartment with a soft dry cloth each time you install batteries.
- Batteries are hazardous waste. DO NOT dispose of them together with household garbage.**
- DO NOT dispose of batteries in fire. Batteries may explode or leak.**
- Recycle or dispose of properly in accordance with local, state, province, and country regulation.**


#### Connect to cuff

After removing the cuff, connect the tube of the cuff to the tube connection hole on the left side of the unit. The air pipe should be properly connected when connecting to prevent the air path from being bent or blocked.



#### Setting time, date and kPa/mmHg unit selection

The unit automatically records the date and time of each measurement, which is key information. Because the blood pressure of the human body is constantly changing on the same day, it is recommended that you set the correct date and time immediately after installing the battery. Please set the correct date and time as following. (Example: enter time and day :12:08 pm, on Dec. 8th).


- In the shutdown mode, press the [SET] key, the product enters the function setting mode (if not selected, the system defaults to the manufacturing time). The last two digits of the year will flash on the screen. You can enter the year by pressing the [MEM] button (for example, you can increase the year from 22 to 23 by pressing it one time).


- Press the [SET] button, and there will be a jump of the month in the left corner of the screen. At this time, the first digit of the month starts to flash. You can enter the month by pressing the [MEM] button (for example, you can increase the month from January to December by pressing 11 times).

- Press the [SET] button again, the screen will switch to the time setting, at this time, the last two digits of the date will start to flash. You can enter the date by pressing the [MEM] button (for example, you can change the date from 1st to 08th by pressing the button 7 times).


- Press the [SET] button again, the screen will switch to the time setting. At this time, the first digit (indicating hour) starts to flash. Please input the time by pressing the [MEM] button (for example, you can change the time from 0:00 to 12:00 by pressing the [MEM] button 12 times).

- Press the [SET] button again, and the last two digits (in minutes) will start to flash. You can enter the minute by pressing the [MEM] button (for example, pressing the button 8 times will increase the time from 0 minutes to 8 minutes).


- In the last step, press the [SET] button, and the date will be displayed after the setting is completed (the date will be displayed first, and it will automatically jump to the time after 3 seconds). When there is no operation, it will automatically shut down after 30 seconds.

### MEASUREMENT OF BLOOD PRESSURE

#### Preparation before measurement

- Avoid eating, smoking or any form of activity before measuring, all of which will affect the measurement. Try to relax in a quiet environment, rest for 10 minutes, and then take measurements.
- If you have heavy clothing, please remove from the upper arm.
- Select the arm for blood pressure measurement (usually the left arm) and take blood pressure measurements on the same arm and in the same area.
- Blood pressure measurements should be taken at the same time every day.


**Common factors that can lead to erroneous measurements.**  
**NOTE:** Comparative blood pressure measurements should be performed in the same condition (usually referred to quiet conditions). If the arm artery is low (high) relative to the heart, it will result in a high (low) blood pressure value. (Each 15cm height difference will produce 1.3kPa/mmHg error value)!

- This unit is not suitable for newborns.
- NOTE:** Please use the original arm cuff that meets the clinical test requirement!
- Too loose of a cuff will lead to a false measurement.
- The cuff on the arm will compress the blood vessels due to repeated measurements. This condition will also cause a biased blood pressure value. Therefore, when making repeated measurements, be sure to rest for 3-5 minutes or raise your arm for 3 minutes.

#### The correct method of using the cuff:

Make sure the cuff tube connector is fully inserted into the cuff tube connector hole on the left side of the monitor. Please take off your coat, sweater and other thicker clothes, and any accessories, use a bare upper arm or wear a thin shirt for measurement. Wrap the cuff around the left arm. DO NOT wrap the cuff too tightly (it is the best to easily insert a finger) and the lower edge of the cuff should be 0.8-1.2in(2-3)cm above the elbow (as shown). After wrapping the cuff around the upper arm, place the air tube on the inside of the arm and on the same line as the middle finger (as shown). If you cannot measure your left hand blood pressure, measure your right hand blood pressure. When measuring, please sit in the chair, place your feet flat on the floor, place your arms on the table, and put the arms and the heart on the same level to relax and measure in a relaxed and natural posture.

### ERROR WARNING/TROUBLESHOOTING

**False Alarm**  
The LCD will display an error warning if any of the following occurs, as shown in figure marked:

Error code	Description	Reason of error	Solutions
E01	Error of measurement	Too much noise to detect an effective pulse signal	Please rest for 2 minutes, adjust the cuff, keep quiet during measurement, DO NOT move or talk.
		No pulse was detected	
		Test results deviate from normal range	
E02	Abnormal cuff wearing	Static pressure exceeds the set protection point	Measure blood pressure again
		Air pressure/sensor abnormal	1. During the process of returning to zero, the air pressure in the cuff continues to fluctuate; 2. The sensor part of the circuit is abnormal; 3. Sensor damage;
Lo	Power is not enough	The power is lower than the minimum operating voltage.	Replace the battery

**Over range alarm**  
When the test blood pressure value is higher than 280mmHg (excluding 280mmHg), the blood pressure value will continue to flash at a frequency of 2Hz. When the blood pressure value is lower than 280mmHg, the flicker is released.

**Other Information**  
Even a healthy individual, the blood pressure is constantly changing (presenting a jagged line), so when you make comparative measurements, you must be in a fixed state (quiet environment)! If the difference is greater than 2.0 kPa / 15 mmHg if the above conditions are met, or in some cases, in case of irregular heartbeat, please consult your doctor.

**Troubleshooting**  
If any faults (or abnormal conditions) occur during use, you can check and exclude according to the items listed in the following table:

Malfunction	Exclude
When the battery is installed and the switch is turned on, the LCD shows nothing.	1. Check whether the positive and negative poles of the battery are placed correctly. 2. If the fault persists, please reposition or replace the battery.
The air pump has begun to inflate, but there is no rise in arm pressure.	Check the hose connection for air leakage or whether it is fully inserted into the socket.
Sphygmomanometers fail to measure blood pressure frequently, or the blood pressure is abnormally high or low.	1. Re-set the correct cuff. 2. If the left upper arm of the belt is covered with sleeves or other clothing, please take off. Re-measure blood pressure.
Each measurement was different, although the sphygmomanometer was functional or showed normal blood pressure.	Please study the following points: Re-wrap the cuff properly so that it is positioned correctly, and take new measurement. Consult your doctor.
Self-measured value is different from the measured value of the doctor.	Record daily measurements and consult your doctor.
After the sphygmomanometer is pressurized, the air pressure of the cuff is released. And the rate is slow even not released at all.	The air hole connection of the hose in the arm belt has the occurrence of the "plastic ring" falling off. Please put the plastic ring on and measure again.

**NOTE:** If the unit still does not work, contact Homedics Consumer Relations. Under no circumstance should you disassemble or attempt to repair the unit by yourself. Contact information for Homedics Consumer Relations Department can be found on the warranty page.

**Memory—Using the Memory Function**  
This monitor automatically stores up to 60 memories for one user. When the monitor is off, press the [MEM] button to display the average value of the last

- Measurements. Press [MEM] button again, will display (MEM1-60). After that press [SET] button will display the last reading taken with date shown first, press [SET] button again, will display (MEM60-1).

**NOTE:** Each memory will show the date first, then after a few seconds the time will show on the display.

**Stop measurement**  
If there is any discomfort or for some reason, such as when the airbag is in a persistent over-inflated state, there may be a risk. At this time, the blood pressure measurement must be stopped. You can press the [START/STOP] button and the sphygmomanometer will immediately release the air pressure in the cuff.

**Memory-erasing all memory**  
**Please Note:**

- Before you delete the memory, please confirm carefully stored data will not be used in the future. The best way is to keep good records so that you can provide your doctor with the necessary information.
- If you decide to permanently delete all stored values, press and hold the [MEM] and [START/STOP] buttons at the same time until the display position of SYS and DIA and heart rate shows horizontal bars "---" then release the button.
- The machine cannot erase individual numeric readings.

**NOTE:** Once deleted, your readings cannot be restored.

### REFERENCE STANDARD

**This blood pressure monitor utilizes advanced technology that has been clinically proven to meet the following standards:**

Performance safety standard: IEC60601-2-30  
Electromagnetic compatibility: IEC60601-1-2  
Safety standard: IEC60601-1

**Calibration**  
The electronic sphygmomanometer has been calibrated at the time of manufacture. We recommend a static pressure test for this sphygmomanometer every 2 years, with the option of an authorized dealer to calibrate your device. At any time, if you question the accuracy of the measurement, please feel free to contact your dealer or manufacturer for disposal.

### DISPOSAL

If this unit is damaged and needs to be discarded, please dispose of the discarded electronic waste in accordance with the relevant local ordinances, national laws and regulations. Disposal of the battery or product should not be directly placed in the garbage can.

### TECHNICAL SPECS

Model	FDBP-A1/BPA-800-WT
Display	LCD displayer
Measurement Method	Oscillometric method
Measuring scope	Systolic pressure(SYS):60mmHg-250mmHg Diastolic pressure(DIA) 30mmHg-195mmHg Pulse: (40-199) times / minute
Memory	60 sets
Resolution	0.1kPa(1mmHg)
Static pressure scope	0-295mmHg
Accuracy	Static pressure: ±0.4 kPa (±3 mmHg) Pulse: within ±5% of the reading
Power	DC 6.0V(4*1.5V AA Batteries)
Special accessory	Cuff, instruction manual, 4pcs AA batteries 1.5V
Size	(L) 135mm* (W) 110 mm* (H) 75.5 mm
Weight	About 368g (including battery)
Withstand pressure for the cuff	360mmHg
Upper arm circumference	8.7-16.5in
Electric shock protection type	Internally powered ME equipment
Shock protection procedure	BF application part
Expected use lifetime	Body in 5 years, Cuff in 2 years
Application component	Cuff
Operating condition	Temperature: 5°C-40°C Humidity: 15% RH-90% RH. No condensation Atmospheric pressure: 70kPa-106kPa
Transportation and storage condition	Temperature: -25°C-55°C Humidity: 15% RH-95% RH. No condensation Atmospheric pressure: 70kPa-106kPa Please strictly observe the environmental conditions of transportation and storage, otherwise it will affect the accuracy of the equipment.

## INCLUDED IN BOX

Component	Quantity
Body	1Set
Cuff	1pcs (with trachea)
AA Battery	4pcs

## STANDARDIZED SYMBOLIC DESCRIPTION

	<b>NOTE:</b> please refer to the instruction before using
	Attention: see Instructions for use!
	Type BF applied parts
	Disposal in accordance with Directive 2002/96/EC (WEEE)
	Manufacturer: Famidoc Technology Co., Ltd
<b>SYS</b>	The value of SYSTOLIC BLOOD PRESSURE
<b>DIA</b>	The value of DIASTOLIC BLOOD PRESSURE
	Temperature limitation
	Humidity limitation
	Atmospheric pressure limitation
<b>IP21</b>	Level of protection for ingress of water or particulate matter into ME EQUIPMENT

## FEDERAL COMMUNICATIONS COMMISSION COMPLIANCE STATEMENT

### FCC CAUTION

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**NOTE:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

## ELECTROMAGNETIC COMPATIBILITY INFORMATION

### ▲NOTE:

- The Upper Arm Blood Pressure Monitor FDBP-A1/BPA-800-WT meets the electromagnetic compatibility requirements of IEC60601-1-2.
- Users should install and use the electromagnetic compatibility information provided by the random files. Portable and mobile RF communication devices may affect the performance of the Upper Arm Blood Pressure Monitor FDBP-A1/BPA-800-WT. Avoiding strong electromagnetic interference when used, such as close to mobile phones, microwave ovens, etc.
- The instructions for the guide and the manufacturer are detailed in the attachment.

### ▲WARNING:

- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this the FDBP-A1/BPA-800-WT could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the FDBP-A1/BPA-800-WT, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

### Attachment:

Guidance and manufacturer's declaration – electromagnetic emission – for all EQUIPMENT AND SYSTEMS

Guidance and manufacturer's declaration – electromagnetic emission		
The FDBP-A1/BPA-800-WT is intended for use in the electromagnetic environment specified below. The customer or the user of FDBP-A1/BPA-800-WT should assure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment- Guidance
RF emissions CISPR 11	Group 1	The FDBP-A1/BPA-800-WT uses RF energy only for its internal function. There for, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The FDBP-A1/BPA-800-WT suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity – for all EQUIPMENT and SYSTEMS

Guidance and manufacturer's declaration – electromagnetic immunity

The FDBP-A1/BPA-800-WT is intended for use in the electromagnetic environment specified below. The customer or the user of the FDBP-A1/BPA-800-WT should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment- Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output linesww	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	± 1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0,5 cycle UT At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycle	0 % UT; 0,5 cycle UT At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the FDBP-A1/BPA-800-WT requires continued operation during power mains interruptions, it is recommended that the FDBP-A1/BPA-800-WT be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30A/m	30A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

**NOTE:** UT is the a. c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity – for EQUIPMENT and SYSTEM

Guidance and manufacturer's declaration – electromagnetic immunity			
The FDBP-A1/BPA-800-WT is intended for use in the electromagnetic environment specified below. The customer or the user of the FDBP-A1/BPA-800-WT should assure that it is used in such an environment			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment- Guidance
Conducted RF IEC 61000-4-6	3 Vrms	3V	Portable and mobile RF communications equipment should be used no closer to any part of the FDBP-A1/BPA-800-WT, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \left[ \frac{W}{P} \right]^{1/2}$ 80 MHz to 800 MHz: $d = \left[ \frac{3.2}{P} \right]^{1/2}$ 800 MHz to 2.7 GHz where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). <sup>3</sup> Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. <sup>3</sup> Interference may occur in the vicinity of equipment marked with the following symbol: 
	150 kHz to 80 MHz	150 kHz to 80 MHz	
Radiated RF IEC 61000-4-3	6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz	6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz	
	10 V/m	10 V/m	
	80 MHz to 2.7 GHz	80 MHz to 2.7 GHz	
	385MHz-5785MHz Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communication equipment (Refer to table 9 of IEC 60601-1-2:2014)	385MHz-5785MHz Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communication equipment (Refer to table 9 of IEC 60601-1-2:2014)	

**NOTE 1:** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2:** These guidelines may not apply in all situations. Electromagnetic is affected by absorption and reflection from structures, objects and people.

a The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 19,07 MHz to 19,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

b Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the FDBP-A1/BPA-800-WT is used exceeds the applicable RF compliance level above, the FDBP-A1/BPA-800-WT should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the FDBP-A1/BPA-800-WT.

c Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM - for EQUIPMENT and SYSTEM

Recommended separation distances between portable and mobile RF communications equipment and the FDBP-A1/BPA-800-WT.

The FDBP-A1/BPA-800-WT is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the FDBP-A1/BPA-800-WT can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the FDBP-A1/BPA-800-WT as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output of transmitter W	Separation distance according to frequency of transmitter			
	150 kHz to 80 MHz outside ISM and amateur radio bands $d = \left[ \frac{3.2}{P} \right]^{1/2}$	150 kHz to 80 MHz in ISM and amateur radio bands $d = \left[ \frac{17}{P} \right]^{1/2}$	80 MHz to 800 MHz $d = \left[ \frac{3.2}{P} \right]^{1/2}$	800 MHz to 2.7 GHz $d = \left[ \frac{17}{P} \right]^{1/2}$
0.01	0.02	0.20	0.035	0.07
0.1	0.38	0.63	0.11	0.22
1	1.2	2.00	0.35	0.70
10	3.8	6.32	1.10	2.21
100	12	20.00	3.5	7.0

For transmitters rated at a maximum output power not listed above the recommended separation distance in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE 1:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## Homedics

### LIMITED 5-YEAR WARRANTY

Homedics sells its products with the intent that they are free of defects in manufacture and workmanship for a period of 5 years from the date of original purchase, except as noted below. Homedics warrants that its products will be free of defects in material and workmanship under normal use and service. This blood pressure monitor meets the simulated measurement cycles test requirement per EN1060-3, part 8.10. This warranty extends only to consumers and does not extend to Retailers.

To obtain warranty service on your Homedics product, contact a Consumer Relations Representative. Please make sure to have the model number of the product available.

Homedics does not authorize anyone, including, but not limited to, Retailers, the subsequent consumer purchaser of the product from a Retailer or remote purchasers, to obligate Homedics in any way beyond the terms set forth herein. This warranty does not cover damage caused by misuse or abuse; accident; the attachment of any unauthorized accessory; alteration to the product; improper installation; unauthorized repairs or modifications; improper use of electrical/power supply; loss of power; dropped product; malfunction or damage of an operating part from failure to provide manufacturer's recommended maintenance; transportation damage; theft; neglect; vandalism; or environmental conditions; loss of use during the period the product is at a repair facility or otherwise awaiting parts or repair; or any other conditions whatsoever that are beyond the control of Homedics.

This warranty is effective only if the product is purchased and operated in the country in which the product is purchased. A product that requires modifications or adoption to enable it to operate in any other country than the country for which it was designed, manufactured, approved and/or authorized, or repair of products damaged by these modifications is not covered under this warranty.

THE WARRANTY PROVIDED HEREIN SHALL BE THE SOLE AND EXCLUSIVE WARRANTY. THERE SHALL BE NO OTHER WARRANTIES EXPRESS OR IMPLIED INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS OR ANY OTHER OBLIGATION ON THE PART OF THE COMPANY WITH RESPECT TO PRODUCTS COVERED BY THIS WARRANTY. HOMEDICS SHALL HAVE NO LIABILITY FOR ANY INCIDENTAL, CONSEQUENTIAL OR SPECIAL DAMAGES. IN NO EVENT SHALL THIS WARRANTY REQUIRE MORE THAN THE REPAIR OR REPLACEMENT OF ANY PART OR PARTS WHICH ARE FOUND TO BE DEFECTIVE: WITHIN THE EFFECTIVE PERIOD OF THE WARRANTY. NO REFUNDS WILL BE GIVEN. IF REPLACEMENT PARTS FOR DEFECTIVE MATERIALS ARE NOT AVAILABLE, HOMEDICS RESERVES THE RIGHT TO MAKE PRODUCT SUBSTITUTIONS IN LIEU OF REPAIR OR REPLACEMENT.

This warranty does not extend to the purchase of opened, used, repaired, repackaged and/or resaled products, including but not limited to sale of such products on internet auction sites and/or sales of such products by surplus or bulk resellers. Any and all warranties or guarantees shall immediately cease and terminate as to any products or parts thereof which are repaired, replaced, altered, or modified, without the prior express and written consent of Homedics. This warranty provides you with specific legal rights. You may have additional rights which may vary from state to state. Because of individual state regulations, some of the above limitations and exclusions may not apply to you.

For more information regarding our product line in the USA, please visit: [www.homedics.com](http://www.homedics.com)

**For service or repair, do NOT return this unit to the retailer, Contact Homedics Consumer Relations at:**

**Email:**  
[cservice@homedics.com](mailto:cservice@homedics.com)

**Phone:**  
1-800-466-3342  
**Business Hours:**  
8:30am-7pm EST  
Monday-Friday

**In USA Distributed by:**  
**Homedics, LLC**  
3000 N Pontiac Trail  
Commerce Township, MI 48390  
Printed in China

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