Homedics

500 series voice out guide upper arm blood pressure monitor

BPA-900BT-WT
Doc # L-03019, Rev.1
5 year
limited warranty

IMPORTANT PRODUCT NOTICES AND SAFETY INSTRUCTIONS

When using your blood pressure monitor, basic precautions should always be followed. Please read and follow all instructions and warnings before using this product.

SAVE THESE INSTRUCTIONS FOR FUTURE REFERENCE.

This device is intended for over-the-counter home use in adults 18 years and older with an arm circumference ranging from 9 inches to 17 inches (approx. 23cm to 43cm). This device features a built-in Bluetooth® data transmission function, which enables the device to automatically transmit measuring results to a paired Bluetooth®-enabled device. This function allows users to synchronize to current date and time, and check the battery status on the **Homedics Health+** app.

- Please note that this is a home healthcare product only and this manual and device are not intended to serve as a substitute for the advice of a physician or medical professional.
- This device uses the oscillometric method to measure systolic and diastolic blood pressure, as well as heart rate.
- DO NOT use this device for diagnosis or treatment of any health problem or disease. Measurement results are for reference only. Consult a healthcare professional for interpretation of pressure measurements. Contact your physician if you have or suspect any medical problem. Do not change your medications without the advice of your physician or healthcare professional.
- Proper cuff size is critical for accurate measurements. Follow the instructions in this manual and printed on the arm cuff to ensure the arm cuff is used properly.

This product is not suitable for:

- Pregnant women
- People with arrhythmias
- People undergoing intravenous injection on any limb
- People currently in a dialysis treatment
- Pregnant women in pre-eclampsia condition
- Those who have had a mastectomy surgery (especially when lymph nodes have been removed), it is recommended to take measurements on the unaffected side.
- This device may have difficulty determining proper blood pressure for users with irregular heartbeat, diabetes, liver disease, kidney disease, poor circulation of the blood or for users who have suffered a stroke. Please consult your healthcare professional before using this device.
- Excessive use may result in blood flow interference, which is likely to cause uncomfortable sensations, such as partial subcutaneous hemorrhage, or temporary numbness to your arm. In general, these symptoms should not last long. However, if you continue to experience these sensations, please seek advice from a medical professional.
- When used along with other medical electronic equipment on the same arm, pressurization of the cuff may cause the other devices to temporarily malfunction.

- The pulse display is not suitable for checking the frequency of heart pacemakers.
- Electromagnetic interference: The device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens). These may lead to temporary impairment of measurement accuracy.
- Consider the electromagnetic compatibility of the device (ex. power disturbance, radio frequency interference etc.). Please use this device indoors only.
- Use blood pressure monitor only for its intended use.
- DO NOT wrap the cuff around body parts other than your arm.
- Not for use by or on persons under the age of 18.
- DO NOT use this device on infants, children, or those who cannot express their own intention.
- DO NOT plug or unplug the adapter power cord with wet hands.
- Please use only the AC adapter indicated for use with this monitor or 1.5V AAA alkaline batteries for power supply.
- Please rest for at least 5-10 minutes before taking a measurement.
- To allow your blood vessels to return to the condition prior to taking the measurement, please wait at least 3-5 minutes between measurements. You may need to adjust the wait time according to your personal physiological situation.
- Wait 30-45 minutes before measurement if you've just consumed caffeinated beverages or smoked cigarettes.
- The applied part is the cuff.
- The patient is the intended operator.

Blood pressure measurements determined with this device are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultatory method, within the accuracy limits prescribed by the *American National Standard, Manual, electronic, or automatic sphygmometers*.

In the event that the device needs to be checked for calibration, contact the distributor.

ABOUT BLOOD PRESSURE

What is blood pressure?

Blood pressure is the measurement of the force exerted on the artery walls while blood flows through the arteries. The pressure measured when the heart contracts and sends blood out of the heart is systolic (highest) blood pressure. The pressure measured when the heart dilates with blood flowing back into the heart is called diastolic (lowest) blood pressure.

Why measure your blood pressure?

Among today's various health problems, those associated with high blood pressure are very common. High blood pressure dangerously correlates with cardiovascular diseases. Therefore, blood pressure monitoring is important for identifying those at risk.

Why do my readings vary?

Blood pressure is a body parameter that is subject to normal variations throughout the day. A single reading that is different from yours or your doctor's readings are not necessarily inaccurate. The average of several readings, taken under similar conditions, using the same arm is preferred for accurate blood pressure readings. Medical physicians generally recommend the "Rule of 3", where you are encouraged to take your blood pressure three times in a row (at 3 ~ 5 minute intervals), three times a day for three days. After three days you can average all the results and this will give you an accurate idea of what your blood pressure really is.

Why Are My Readings Different Than Those Taken at My Doctor's Office?

Many experience a phenomenon called "White Coat Hypertension" when measured by a doctor. White Coat Hypertension refers to blood pressure that rises above its usual level when measured in a clinical setting, such as a doctor's office.

Bluetooth®

The table below contains defined levels for hypertension that are publicly available from the American Heart Association® (AHA 2017). Users can compare their own blood pressure readings against these defined levels to determine if they may be potentially at increased risk.

Blood Pressure Category	Systolic mm Hg (upper number)	and	Diastolic mm Hg (lower number)	Indicator Color
Normal	<120	and	<80	Green
Elevated	120-129	and	<80	Yellow
High Blood Pressure (hypertension) Stage 1	130-139	or	80-89	Red
High Blood Pressure (hypertension) Stage 2	≥140	or	≥90	Red
Hypertension Crisis (consult your doctor immediately)	>180	and/or	>120	Red

*Source: American Heart Association® (AHA) 2017


Blood pressure tends to go up and down, even in people who normally don't have high readings. If your numbers stay above the "normal" range most of the time, you may be at increased risk and should consult your physician. Although one can easily find where their own blood pressure readings fall on this table, this monitor comes equipped with a Risk Category Index that automatically compares each reading to the defined levels and provides a helpful cue if your reading falls into one of the stages that could potentially indicate increased risk. See *Risk Category Index* section for more information on this feature. Please note that cues provided by this monitor are only intended to assist you in using this table. The table and cues are only provided for convenience to help you understand your non-invasive blood pressure reading as it relates to the American Heart Association® (AHA 2017) information. They are not a substitute for a medical examination or diagnosis by your physician. It is important for you to consult with your pressure range as well as the point at which you may actually be considered to be at risk.

HOW THIS BLOOD PRESSURE MONITOR WORKS


This monitor uses Smart Measure® inflation technology to detect your blood pressure. Once the cuff is wrapped around your upper arm, with a press of the START/STOP button, the monitor will turn on and the inflation will automatically start, creating pressure around the arteries inside the upper arm. Within the cuff is a gauge which senses the fluctuations (oscillations) in pressure. The fluctuation measured represents the degree of intensity that your arteries are contracting with each heartbeat, which is also a result of the pressure that the cuff has placed on the upper arm. The monitor measures these contractions and converts the information to a digital value. Once the measurement is complete, the cuff will automatically deflate, and the result is displayed on the screen. Please note that any muscle movement during inflation will cause measurement error. When measurement is complete, the monitor will display your systolic pressure, diastolic pressure, and pulse readings. Before measurement, it is suggested that you sit quietly for 15 minutes before measurement as measurements taken in a relaxed state have greater accuracy. The monitor automatically finds where your measurement results fall on the American Heart Association® (AHA 2017) table and provides a cue if your reading falls into one of the stages that could potentially indicate increased risk. Please refer to the *Risk Category Index* section for more information on this feature.

The appearance of the icon indicates that a pulse irregularity consistent with an Irregular Heartbeat was detected during measurement. Please refer to *Irregular Heartbeat Detector* section for more information on this feature.

NAME/FUNCTION OF EACH PART

DISPLAY SYMBOL EXPLANATIONS

Display	Explanation Symbols
no. 1	User 1: Appears when the monitor is operated by User 1.
no. 2	User 2: Appears when the monitor is operated by User 2.
	Low Battery Symbol: Appears when batteries should be replaced.
	Pulse Symbol: Once pulse is detected, the symbol flashes with each pulse beat. Shows the heart rate per minute.
	Irregular Heartbeat Detector: This symbol appears when an irregular heart beat was detected. See Irregular Heartbeat Detector section for more information. If symbol appears frequently, please consult your physician.
	Risk Category Index: The arrow points out the specific Risk Category that your measurement reading falls in. See Risk Category Index Section for more information.
	Excessive Body Motion Detector: Displayed if body movement is detected during measurement, especially, the movement on the arm the blood pressure monitor is worn on. If the cuff is worn improperly, or the shape of the upper arm is unusual (for example, the circumference of the upper arm differs largely from the circumference of the forearm), excessive gap might exist between the arm cuff and the arm. NOTE: The measured blood pressure reading may not be accurate if the icon is displayed.
	Bluetooth® Symbol: Appears when Bluetooth is active, and when data is transmitting to your mobile device.
AVG. 3	Average of the Last 3 Measurements: This symbol appears when the LCD displays the average value of last 3 readings.

If any of the following letters and numbers appear in the area that systolic pressure should be displayed, an error has occurred with your reading. See Troubleshooting section of this manual for more information.

EE	Measurement Error: Wrap the cuff correctly and keep arm steady during measurement. Take new measurement.
E1	Air Circuit Abnormality: Check cuff connection. Take new measurement.
E2	Pressure Exceeding 300 mmHg: Turn the unit off to clear, then take new measurement.
E3	Error Determining Measurement Data: Ensure the cuff is worn properly, do not move or talk. Take new measurement.
E4	Data Transmission Error: Will occur after 45 seconds, if the monitor cannot connect to the mobile device to transmit data. Make sure Bluetooth® is ON. Press button for 3 seconds to start Bluetooth® function.
EP	System Error: Contact Consumer Relations

INSTALLING BATTERIES


- Press and lift the battery cover to open the battery compartment.
- Install or replace 4 AAA (1.5V, LR03) alkaline batteries in the battery compartment according to the indications inside the compartment.
- Close the battery cover by pressing the cover back into place.

- Replace the batteries if:**
- The low battery symbol appears.
 - When any button is pressed and nothing is displayed on the screen.

NOTE:

- If the batteries are removed or replaced, the date and time will need to be re-set either manually or automatically using your Bluetooth® mobile device.
- Replace all batteries at one time (as simultaneous set). Use only 1.5V AAA alkaline batteries. DO NOT mix alkaline, standard (carbon-zinc) or rechargeable (cadmium) batteries.
- When the batteries are removed, the measurement values stored in memory are retained. However, the date and time must be reset.
- Remove batteries when unit is not in use for extended periods of time.
- Clean contacts on battery and in battery compartment with a soft dry cloth each time you install batteries.
- Batteries are hazardous waste. **DO NOT** dispose of them together with the household garbage.
- DO NOT** dispose of batteries in fire. Batteries may explode or leak.
- Recycle or dispose of properly in accordance with local, state, province, and country regulations.


- This monitor is designed to operate with a 6V DC AC adapter, sold separately. If you have any additional questions, or need help on purchasing, please contact Homedics Consumer Relations, please see warranty section.
- Connect the AC adapter with the AC adapter jack in the side of the unit as shown below.

- Use only the AC adapter specified.

NOTE:

- When the AC adapter is your main power supply, make sure the adapter plug can be easily removed from the unit.
- Please unload the batteries when operating with an AC adapter for an extended period of time. Leaving the batteries in the compartment for a long time may cause leakage, which may lead to damage to the unit.
- No batteries are needed when operating with an AC adapter.
- Date and time will need to be reset if AC adapter is unplugged even if the unit has the batteries in the battery compartment.

DATE & TIME SET PROCEDURE

Date & Time can be set by two methods, either sync automatically using your Bluetooth® mobile device, or set manually using the and buttons on the monitor.

If setting the Date & Time by syncing with your mobile device, it is important that this is done prior to taking any measurements to ensure the date and time are accurate.

To set using your device with Bluetooth® wireless technology:

- Before taking the first measurement, open the **Homedics Health+** app on your device. If you have not already installed the **Homedics Health+** app onto your device, it is available to download on the App Store and on Google Play™. Make sure you have the app downloaded and open on your mobile device before trying to sync the Date & Time.
- Press and hold the Bluetooth® button for 3 seconds on the monitor to enable the connection.
- Once the Bluetooth® connection is established, the Date & Time will automatically update on the blood pressure monitor.

NOTE:

- Once the Date & Time have been successfully synced, future readings will automatically have the correct date and time.
- If you are having trouble automatically syncing the Date & Time, make sure the Bluetooth® is ON on your device and ON on the blood pressure monitor.

To set Manually:

To set the date and time when monitor is ON, press the button for 3 seconds.

The display will show a blinking number showing the YEAR. Change the YEAR by pressing the button. Each press will increase the number by one in a cycling manner. Press the button again to confirm the entry, and the screen will show a blinking number representing the MONTH.

Change the DAY, HOUR and MINUTE as described in Step 2 above, using the button to change the numbers and the button to confirm the entries.

NOTE:

- The date and time will only need to be set manually if the monitor will not be used with the **Homedics Health+** app.
- If the **Homedics Health+** app is used after manually setting the date and time on the blood pressure monitor, the date and time on the mobile device will override the date and time on the blood pressure monitor.


BLUETOOTH® OPERATION

This monitor has Bluetooth® function ON for your convenience. This will allow your readings to automatically transmit to the **Homedics Health+** app. By turning this function OFF, the measurements cannot be transmitted.

Turning the Bluetooth® Function ON/OFF:

While the screen is on, press and hold the button for 3 seconds to enter date/time set mode.

- Once the date and time has be set press the button again, then press the button to turn the Bluetooth® ON or OFF.

App for Bluetooth®

Download and install the **Homedics Health+** app on your smart device from App Store® and on Google Play™.

System requirement of the Bluetooth-enabled device

- Bluetooth® 5.0 for Android 6.0 or above
- Bluetooth® 5.0 for IOS 9.0 or above

TALKING FUNCTION

The Talking Function will guide you through the measurement procedure and also announce measurement results. This monitor has Talking Function ON for your convenience. The Talking Function is not a substitute for reading and understanding the instruction manual.

The Talking Function has Two Modes:

Instructions ON Mode - Announces instructions and results

Instructions OFF Mode - No voice or announcements

This function default setting is on. When monitor is in standby mode (Monitor is ON) user can switch the function on or off. By pressing button for 3 seconds, the device will enter to setting year, date and time first. When these settings are done, press the button again, use button to turn Talking Function ON or OFF.


USING THE UPPER ARM CUFF

Very Important: Proper cuff size is critical to an accurate measurement.


This monitor comes with a one size fits most upper arm cuff that fits arm sizes 9"-17" (23cm-43cm).

The appropriate cuff is suitable for your use if the colored arrow falls within the solid color line as shown below.


If the arrow falls outside of the solid color line you will need to purchase an XL arm cuff, model# BPA-CUFFFLA. If you have any additional questions, or need help on purchasing, please contact Homedics Consumer Relations, please see warranty section.

- If the cuff is not assembled, pass the end of the cuff furthest away from the tubing through the metal D-ring in order to form a loop. The smooth side without the felt material should be on the inside of the cuff loop.


- Plug the cuff tube into the left side of the unit.

- Open the arm cuff so that your arm may be placed through it.

- Wrap the cuff on a bare arm or over thin clothing. Thick clothing or a rolled up sleeve will cause inaccurate blood pressure measurements.

- Position cuff edges 0.8-1.2in (2-3cm) above elbow.


- Center tube over middle of arm.

- Pull the end of the cuff so that it tightens evenly around your arm. Press the hook and loop material together to secure. Allow room for 2 fingers to fit between the cuff and your arm. Please make sure the cuff does not slip during measurement, and the arrow falls within the proper fit range.


- Lay your arm on a table (palm upward) so the cuff is at the same height as your heart. Make sure the tube is not kinked.

NOTE:

- This device should not be used when your arm is injured or wounded.
- If it is not possible to fit the cuff on your left arm, it can also be placed on your right arm. However, all measurements should be made using the same arm.
- To use on the right arm, you must position the artery symbol over the main artery. Locate the main artery by pressing with two fingers approximately 1" (2cm) above the bend of your elbow on the inside of your right arm. Identify where the pulse can be felt the strongest. This is your main artery.

MEASUREMENT PROCEDURE

Important Notes:

Blood pressure changes with every heartbeat and is in constant fluctuation throughout the day.

Blood pressure measurement can be affected by the position of the user, his or her physiologic condition and other factors.

Before Measurement:

To help ensure an accurate reading, follow these instructions before taking a measurement:


- Wait 1 hour after exercising, bathing, eating, drinking beverages with alcohol or caffeine, or smoking to measure blood pressure.
- Sit quietly and rest for 15 minutes.
- Stress raises blood pressure. Avoid taking measurements during stressful times.
- Take your reading in a comfortable environment as measurements can be affected by hot or cold temperatures. Take your blood pressure at normal body temperature.

During Measurement:

- DO NOT talk or move your arm or hand muscles.
- DO NOT cross your legs. Sit with feet flat on the floor.
- DO NOT touch cuff or monitor during measurement.

Taking a Measurement:

- Under Standby Mode, press button to select a USER 1 or 2.

- With the cuff wrapped around your arm, press start/stop button to confirm the chosen user and start measurement.

- All display symbols appear on the screen for 1.5 seconds. After all symbols disappear, the display will show "00". The monitor is "Ready to Measure" and will automatically inflate to the level that is right for you.

- As the cuff inflates, the monitor automatically determines your ideal inflation level. This monitor features Smart Measure® inflation technology, and detects your blood pressure and pulse rate during inflation. The Heartbeat Symbol flashes at every heartbeat. Remain still and do not move until the entire measurement process is completed. The device will detect your pulse and determine the blood pressure.

NOTE

- If the cuff does not stop inflating, remove the cuff at once.
- To stop measurement, press START/STOP button. The cuff will deflate immediately after the button is pressed.
- After the monitor has determined your blood pressure and heart rate, the cuff automatically deflates. Your systolic rate, diastolic rate, heart rate and corresponding Risk Category Indicator, Irregular Heartbeat Detector and Excessive Body Motion Detector (if any) are displayed with date and time for 1 minute and save results to memory automatically.
- Device automatically shuts off if no operation over 1 minute.

TRANSMIT READINGS FROM YOUR MOBILE DEVICE

This monitor syncs your blood pressure readings to the **Homedics Health+** app that is free to download on the App Store® and on Google Play™. Make sure you have the **Homedics Health+** app downloaded, and open on your mobile device before trying to transmit your blood pressure measurements.

Before attempting to sync the blood pressure monitor with your mobile device, make sure Bluetooth® is turned ON, on both your mobile device and the monitor.

Transmit Readings

There are 2 ways to activate Bluetooth® function.

- Automatically activate:** When measurement is completed, the device activates Bluetooth® function automatically, and the Bluetooth® symbol will flash on the screen.
- Manually activate:** Under standby mode, user can press the button for 3 seconds to activate Bluetooth function, and Bluetooth® symbol will flash on the screen.
- If BPA-900BT-WT is connected successfully to your smart device, Bluetooth® symbol will appear on the screen.
- If the monitor cannot be connected to paired Bluetooth®-enabled device over 45 seconds, LCD will display Error message "E4" and Bluetooth® will be turned off.


RISK CATEGORY INDEX


This monitor comes equipped with a Risk Category Index that automatically compares each reading to defined levels established by the American Heart Association® (AHA 2017) as described earlier in this manual, and provides a helpful cue if your reading falls into one of the stages that could potentially indicate increased risk. Please note that cues provided by this monitor are only intended to assist you in using this table. The table and cues are only provided for convenience to help you understand your non-invasive blood pressure reading as it relates to the AHA® information. They are not a substitute for a medical examination or diagnosis by your physician. It is important for you to consult with your physician regularly. Your physician will tell you your normal blood pressure range as well as the point at which you may actually be considered to be at risk.

Blood Pressure Category	Systolic mm Hg (upper number)	and	Diastolic mm Hg (lower number)	Indicator Color
Normal	<120	and	<80	Green
Elevated	120-129	and	<80	Yellow
High Blood Pressure (hypertension) Stage 1	130-139	or	80-89	Red
High Blood Pressure (hypertension) Stage 2	≥140	or	≥90	Red
Hypertension Crisis (consult your doctor immediately)	>180	and/or	>120	Red

*Source: American Heart Association® (AHA) 2017

IRREGULAR HEARTBEAT DETECTOR (IHB)

The appearance of the  icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Usually this is not a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. Please note that the device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage. Movement, shaking or talking during the measurement can result in pulse irregularities that may cause the appearance of this icon. Therefore, it is of great importance to not move or talk during measurement. To determine the presence of an irregular heartbeat, the average of the heartbeat intervals is calculated with the first 3 normal effective heartbeat values.

It is important to note that the average is not a strict mathematical averaging of all recorded intervals. At least 3 beats with 25% or greater difference from the average heartbeat interval will generate the  icon on the screen.

Important Information:

This blood pressure monitor is not designed for use by people with arrhythmias not for diagnosing or treating an arrhythmia problem. As a safeguard, we recommend that if you have arrhythmias, such as atrial or ventricular premature beats and atrial fibrillation or any other special conditions, you should check with your physician before using your blood pressure monitor.


USING THE MEMORY FUNCTION

This monitor can be used by 2 users. Each user can store up to 120 measurements, and will automatically replace the oldest data with the new one.

After each measurement, the Systolic and Diastolic Pressure, Heart Rate and Risk Category Indicator will be automatically stored. If an Irregular Heartbeat, or Excessive Body Motion was detected the icon will also be automatically stored with that measurement result in memory.

RECALLING VALUES FROM MEMORY

1. Press the **USER** button to select User 1 or User 2.


2. Press **M/** button to access memory. If there is no stored memory, nothing except month, date, and time will appear on the display.

3. The monitor will first display the calculated average applied to the last 3 memories. 4. Every press of the **M/** button will recall a previous reading. The last reading will be recalled first.

CLEARING VALUES FROM MEMORY

1. Press the **USER** button to select USER 1 or USER 2. 2. Press the **M/** button to enter into memory recall mode. 3. Press and hold the **USER** and **M/** buttons at the same time and the data for the selected user will be erased automatically.

NOTE: Once deleted, your readings cannot be restored.

IMPORTANT NOTES REGARDING YOUR BLOOD PRESSURE MEASUREMENT

- It is suggested that you take your measurements at the same time each day and use the same arm for consistency.
- Users should wait a minimum of 5 minutes before taking additional measurements. More time may be necessary depending upon your physiology.
- The measurement results that users receive are for reference only. If users have any blood pressure concerns, please consult a physician.
- Once inflation reaches 300 mmHg, the unit will deflate automatically for safety reasons.
- This product is not suitable for people with arrhythmias.
- This device may have difficulty determining the proper blood pressure for users with irregular heartbeat, diabetes, poor circulation of blood, kidney problems, or for users who have suffered a stroke.

CARE, MAINTENANCE & CLEANING

- Clean the blood pressure monitor body and cuff carefully with a slightly damp, soft cloth. DO NOT press, DO NOT wash cuff or use chemical cleaner on it. Never use thinner, alcohol or petrol (gasoline) as cleaner.
- Make sure the cuff is completely dry before using.
- Leaky batteries can damage the unit. Remove the batteries when the unit will not be used for a long time.
- Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.
- If the unit is stored near freezing, allow it to acclimate to room temperature before use.
- This blood pressure monitor is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems with this device, please contact Homedics Consumer Relations (contact information can be found on the *Warranty* section).
- DO NOT immerse the unit in water as this will result in damage to the unit.
- DO NOT subject the monitor or cuff to extreme temperatures, humidity, moisture, or direct sunlight. Protect from dust.
- DO NOT fold the cuff and tubing tightly.
- DO NOT disassemble the monitor or cuff. If in need of repair, refer to the warranty section of this manual.
- DO NOT subject the monitor to extreme shocks (do not drop on floor).
- DO NOT inflate the cuff unless wrapped around arm.
- DO NOT wrap the cuff around body parts other than your arm
- DO NOT drop or insert any object into any opening or hose.
- To avoid accidental strangulation keep this product away from children.
- DO NOT drape tube around neck.
- This monitor may not meet its performance specifications if stored or used outside of these temperature and humidity ranges:

STORAGE/TRANSPORTATION ENVIRONMENT

Temperature: -13°F ~ 158°F (-25°C ~ 70°C)
Humidity: less than 93% RH

OPERATING ENVIRONMENT

Temperature: 41°F ~ 104°F (5°C ~ 40°C)
Humidity: 15% ~ 93% RH
Atmospheric Pressure: 700hPa-1060hPa

POTENTIAL FOR ELECTROMAGNETIC INTERFERENCE

To avoid inaccurate results caused by electromagnetic interference between electrical and electronic equipment, do not use the device near a cell phone or microwave oven. For most wireless communication devices, it is recommended to maintain a distance of 10.8 feet (3.3m) in order to avoid electromagnetic interference.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

FEDERAL COMMUNICATIONS COMMISSION COMPLIANCE STATEMENT

Changes or modifications to this equipment not expressly approved by the manufacturer could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, it can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

ELECTROMAGNETIC COMPATIBILITY (EMC)

Guidance and manufacturer's declaration – electromagnetic emissions
The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF emissions CISPR 11	Group 1	RF energy is used only to maintain device's operation. Therefore, its RF emissions are so low that it's not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments, and those directly connected to the public low-voltage power supplies network that supplies buildings used for domestic purposes.
Harmonic Emissions IEX 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic emissions

The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	In the case of air discharge testing, the climatic conditions shall be within the following ranges: Ambient Temperature: 15°C ~ 35°C, Relative Humidity: 30% ~ 60%
Power frequency (50 or 60 Hz) Magnetic field IEC 61000-4-8	30 A/m 50 or 60 Hz	30 A/m 50 or 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	AC Power port ± 1 kV line to line	AC Power port ± 1 kV line to line	Mains power quality should be that of a typical commercial or hospital environment.
Interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycles 70% UT, 25/30 cycles 0% UT; 250/300 cycles	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycles 70% UT; 25 cycles 0% UT; 250 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.

Guidance and manufacturer's declaration – electromagnetic immunity
The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Conducted RF IEC 61000-4-6	3V rms At 0.15-800 MHz 6V rms At ISM & Radio Amateur Freq.	3V rms At 0.15-80 MHz 6V rms At ISM & Radio Amateur Freq.	Portable and mobile RF communications equipment should be separated to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000 4-3	10 V/m at 80-2700 MHz AM Modulation And 9-28V /m at 385-6000 MHz; Pulse Mode and other Modulation. The system shall be tested as specified in IEC60601-1-2 Table 9 for proximity fields from RF wireless communications equipment using the test methods specified in IEC 61000 4-3	10 V/m at 80-2700 MHz AM Modulation and 9-28V /m at 385-6000MHz; Pulse Mode and other Modulation. The system shall be tested as specified in IEC60601-1-2 Table 9 for proximity fields from RF wireless communications equipment using the test methods specified in IEC 61000 4-3	Recommended separation distance considering to reduce the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the distance for higher IMMUNITY TEST LEVELS shall be calculated using the following equation: $E=6/d^2P$ where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVELS in V/m. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.
- Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Test specifications for enclosure port immunity to RF wireless communications equipment.

Test frequency (MHz)	Modulation	Immunity Test Level (V/m)
385	Pulse modulation 18 Hz ^{a)}	27
450	FM± 5 kHz deviation 1kHz sine ^{b)}	28
710	Pulse modulation 217 Hz ^{a)}	9
745		
780	Pulse modulation 18 Hz ^{a)}	28
810		
870		
930	Pulse modulation 217 Hz ^{a)}	28
1720		
1845		
1970	Pulse modulation 217 Hz ^{a)}	28
2450		
5240		
5500	Pulse modulation 217 Hz ^{a)}	9
5785		

NOTE:

If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

- The carrier shall be modulated using a 50% duty cycle square wave signal.
- As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

If any abnormality arises during use, please check the following points:



Symptoms	Possible Causes	Correction
Unit does not turn on when the START/STOP button is pushed.	Batteries have run down.	Replace them with four new AAA (LR03) alkaline batteries.
	Battery polarities (+/-) have been positioned incorrectly.	Re-insert the batteries in the correct positions.
EE measurement error symbol shown on display or the blood pressure value is displayed excessively low (or high).	The Upper Arm Cuff has been positioned on the arm incorrectly.	Re-wrap the cuff properly so that it is positioned correctly. Take new measurement.
E1 error symbol shown on display.	Air circuit abnormality. Cuff tube may not be plugged into monitor correctly.	Check cuff connection. See Using the Upper Arm Cuff section. Take new measurement.
E2 error symbol shown on display.	Inflation pressure exceeding 300 mmHg.	Switch the unit off and then turn On. Take new measurement.
E3 error symbol shown on display.	Error determining measurement data.	Re-wrap the cuff properly so that it is positioned correctly. Take new measurement. If you have any question about the cuff wearing and/or measurement result, please consult your healthcare professional.
	Cuff is worn improperly, or the shape of the upper arm is unusual (for example, the circumference of the upper arm differs largely from the circumference of the forearm), excessive gap might exist between the arm cuff and the arm.	
E4- Data transmission error symbol shown on the display.	Bluetooth® function is turned OFF, on your mobile device.	Turn Bluetooth® ON, on your mobile device, and try again.
	Bluetooth® function is turned OFF, on the blood pressure monitor.	Turn the Bluetooth® ON, on the mobile device and the blood pressure monitor, then try again.
		Make sure Bluetooth® is ON, on the mobile device and the blood pressure monitor, and try again.
	Mobile device does not pair with the blood pressure monitor.	Requires compatible mobile device. This blood pressure monitor is compatible with mobile devices that are running Bluetooth® 5.0 for Android 6.0 or above, Bluetooth® 5.0 for iOS 9.0 or above.
	The App on the mobile device is not open.	Make sure the app is open on the mobile device, and try again.
	The blood pressure monitor and mobile device are out of transmitting range.	Make sure the mobile device and blood pressure monitor are within the acceptable range of 32 feet (10 meters).
	Unexpected loss of electrical/mechanical integrity.	Remove batteries, re-insert, and try again. Return the device to your local distributor or importer.
	Body movement during measurement, especially, the movement on the arm the blood pressure monitor is worn on, e.g. talking, moving or shaking of the arm with the cuff on while measurement.	Measure again. Keep arm steady during measurement.
	Cuff is worn improperly, or the shape of the upper arm is unusual (for example, the circumference of the upper arm differs largely from the circumference of the forearm), excessive gap might exist between the arm cuff and the arm.	Wrap the cuff properly and keep steady. Measure again. If you have any question about the cuff wearing and/or measurement result, please consult your healthcare professional.
BPM cannot communicate with Bluetooth-enabled device	Paring has not been completed.	Please re-pair the BPM and Bluetooth-enabled device with each other.
	Bluetooth® function is not turned on.	See the "Bluetooth® Transmission" section to turn on Bluetooth® function.
	The distance between BPM and Bluetooth-enabled device is out of transmitting range.	Please make sure the acceptable distance (±10 meters) with each other.
	Use an incompatible Bluetooth-enabled device.	Please refer to <i>RF Specification</i> section
	Use non-Bluetooth-enabled device.	
	Unexpected loss of electrical/mechanical integrity.	Re-insert the batteries and try again. Return the device to your local distributor or importer.
EP error symbol shown on display.	System Error.	Measure again. If error persists, contact Consumer Relations.

NOTE: If the unit still does not work, contact Homedics Consumer Relations. Under no circumstance should you disassemble or attempt to repair the unit by yourself. Contact information for Homedics Consumer Relations Department can be found in the *Warranty* section.

SPECIFICATIONS

Model Number	HL85A3-BT/ BPA-900BT-WT
Measurement Method	Oscillometric (inflation)
Rated Range of Cuff Pressure	0 ~ 300 mmHg
Rated Range of Determination	40 ~ 280 mmHg
Measurement Range of Heart Rate	40-199 beats/minute
Accuracy	Pressure: ± 3 mmHg Pulse: ± 5 % Max.
Inflation	Automatic Inflation (Air Pump)
Deflation	Automatic Air Release Control Valve
Display	Liquid Crystal Display
Memory	240 Memory Total for 2 Users
Unit Dimensions	5.51 × 4.33 × 2.22 inch (L x W x H) 140 × 110 × 56.5 mm (L x W x H)
Unit Weight	0.51lbs (230 grams) (Cuff and batteries excluded)
Cuff Size	UC-01: Universal size cuff (9"-17" / 23-43cm)
Storage/Transportation Environment	Temperature: -25°C ~ 70°C (-13°F ~ 158°F) Humidity: ± 93% R.H.
Operation Environment	Temperature: 5°C ~ 40°C (41°F ~ 104°F) Humidity: 15% ~ 93% R.H. Atmospheric pressure: 700hPa ~ 1060hPa
Power Supply	1. DC 6V, AAA "LR03" (1.5V) alkaline battery x 4 2. 5V 1A AC/DC adapter (Optional)
Battery Life	Approx. 250 measurements
Shelf life (battery)	3 years (Temperature: 20 ± 2°C; Relative humidity: 65 ± 20%RH)
Product Life	5 Years (4 times per day)
Sleeping Mode	Without any operation for 1 minute, device automatically shuts off.
Accessories	4 AAA (LR03) Alkaline Batteries, Arm Cuff with Tube, Instruction Manual, Pouch
RF Type	Bluetooth® 5.0 BLE
System Requirement of the Bluetooth enabled device	Bluetooth® 5.0 for Android 6.0 or above Bluetooth® 5.0 for iOS 7.0 or above

*The contents of this manual and the specifications of the device covered by this manual are subject to change for improvement without notice

	Follow instruction for use.
	TYPE BF Applied Part.
	To avoid inaccurate results caused by electromagnetic interference. WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30cm (12 inches) to any part of the device. Otherwise, degradation of the performance of this equipment could result.
	Wase of electrical and electronic equipment (WEEE). Discard the used product to the recycling collection point according to local regulations.
	Serial number  YMMXXXXX
	Ingress Protection Rating. First characteristic numeral - Degree of protection against access to hazardous parts and against solid foreign objects N=2 (Protected against solid foreign objects of 12.5 mm (0.9 and greater). Second characteristic numeral - Degree of protection against ingress of water. N=2 (Protected against vertically falling water drops when ENCLOSURE tilted up to 15°).
	Non-ionizing electromagnetic radiation.

NOTE:

These specifications are subject to change without notice.

DEVICE INFORMATION:

- Internally powered equipment
- Not suitable for use in presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide
- Continuous operation with short-time loading

Homedics

For service or repair, do not return this unit to the retailer, Contact Homedics Consumer Relations at:

Email:
cservice@homedics.com

Phone:
1-800-466-3342
Business Hours:
9:30am-7pm EST
Monday-Friday

Distributed by:
Homedics, LLC
3000 N Pontiac Trail
Commerce Township, MI 48390

Printed in China or Vietnam

To obtain warranty service on your Homedics product, contact a Consumer Relations Representative. Please make sure to have the model number of the product available.

Homedics does not authorize anyone, including, but not limited to, Retailers, the subsequent consumer purchaser of the product from a Retailer or remote purchaser, to obligate Homedics in any way beyond the terms set forth herein. This warranty does not cover damage caused by misuse or abuse; accident; the attachment of any unauthorized accessory; alteration to the product; improper installation; unauthorized repairs or modifications; improper use of electrical/power supply; loss of power; dropped product; malfunction or damage of an operating part from failure to provide manufacturer's recommended maintenance; transportation damage; theft; neglect; vandalism; or environmental conditions; loss of use during the period the product is at a repair facility or otherwise awaiting parts or repair; or any other conditions whatsoever that are beyond the control of Homedics.

This warranty is effective only if the product is purchased and operated in the country in which the product is purchased. A product that requires modifications or adaption to enable it to operate in any other country than the country for which it was designed, manufactured, approved and/or authorized, or repair of product damaged by these modifications is not covered under this warranty.

THE WARRANTY PROVIDED HEREIN SHALL BE THE SOLE AND EXCLUSIVE WARRANTY. THERE SHALL BE NO OTHER WARRANTIES EXPRESS OR IMPLIED INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS OR ANY OTHER OBLIGATION ON THE PART OF THE COMPANY WITH RESPECT TO PRODUCTS COVERED BY THIS WARRANTY. HOMEDICS SHALL HAVE NO LIABILITY FOR ANY INCIDENTAL, CONSEQUENTIAL OR SPECIAL DAMAGES. IN NO EVENT SHALL THIS WARRANTY REQUIRE MORE THAN THE REPAIR OR REPLACEMENT OF ANY PART OR PARTS WHICH ARE FOUND TO BE DEFECTIVE WITHIN THE EFFECTIVE PERIOD OF THE WARRANTY. NO REFUNDS WILL BE GIVEN. IF REPLACEMENT PARTS FOR DEFECTIVE MATERIALS ARE NOT AVAILABLE, HOMEDICS RESERVES THE RIGHT TO MAKE PRODUCT SUBSTITUTIONS IN LIEU OF REPAIR OR REPLACEMENT.

This warranty does not extend to the purchase of opened, used, repaired, repackaged and/or resold products, including but not limited to sale of such products on internet auction sites and/or sales of such products by surplus or bulk resellers. Any and all warranties or guarantees shall immediately cease and terminate as to any products or parts thereof which are repaired, replaced, altered, or modified, without the prior express and written consent of Homedics. This warranty provides you with specific legal rights. You may have additional rights which may vary from state to state. Because of individual state regulations, some of the above limitations and exclusions may not apply to you.

For more information regarding our product line in the USA, please visit: www.homedics.com

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Homedics is under license. App Store is a service trademark of Apple Inc. Android and Google Play are trademarks of Google, LLC. Homedics, Smart Measure and the Smart Measure logo are registered trademarks of Homedics, LLC. American Heart Association is a registered trademark of the American Heart Association, Inc. ©2022 Homedics, LLC. All rights reserved. IE-BPA900BTWT P/N: 323103673

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