## 400 series automatic upper arm blood pressure monitor

INSTRUCTION MANUAL

BPA-H100  
Doc# L-03603, Rev.1  
**5 year limited warranty**

\*Source: American Heart Association® (AHA) 2017


Blood pressure tends to go up and down, even in people who normally don't have high readings. Individual blood pressures vary greatly both on a daily and a seasonal basis. These variations are even more pronounced in hypertensive patients. Normally the blood pressure rises while at work and is at its lowest during sleeping period. If your numbers stay above the "normal" range most of the time, you may be at increased risk and should consult your physician. Although one can easily find where their own blood pressure readings fall on this table, this monitor comes equipped with a Risk Category Index that automatically compares each reading to the defined levels and provides a helpful cue if your reading falls into one of the stages that could potentially indicate increased risk. See Risk Category Index section for more information on this feature.

### HOW THIS BLOOD PRESSURE MONITOR WORKS


This monitor uses Smart Measure® inflation technology to detect your blood pressure. Once the cuff wrapped around your upper arm, with a press of the start/stop button, the monitor will turn on and the inflation will automatically start, creating pressure around the arteries inside the upper arm. Within the cuff is a gauge which senses the fluctuations (oscillations) in pressure. The fluctuation measured represents the degree of intensity that your arteries are contracting with each heartbeat, which is also a result of the pressure that the cuff has placed on the upper arm. The monitor measures these contractions and converts the information to a digital value. Once the measurement is complete, the cuff will automatically deflate, and the result is displayed on the screen. Please note that any muscle movement during inflation will cause measurement error. When measurement is complete, the monitor will display your systolic pressure, diastolic pressure, and pulse readings. Before measurement, it is suggested that you sit quietly for 15 minutes before measurement as measurements taken in a relaxed state have greater accuracy. The monitor automatically finds where your measurement results fall on the American Heart Association® (AHA 2017) table and provides a cue if your reading falls into one of the stages that could potentially indicate increased risk. Please refer to the Risk Category Index section for more information on this feature.

The appearance of the icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Please refer to Irregular Heartbeat Detector section for more information on this feature.

### NAME/FUNCTION OF EACH PART

### DISPLAY SYMBOL EXPLANATIONS

### BLOOD PRESSURE STANDARD

The table below contains defined levels for hypertension that are publicly available from the American Heart Association® (AHA 2017). Users can compare their own blood pressure readings against these defined levels to determine if they may be potentially at increased risk.

BLOOD PRESSURE CATEGORY	SYSTOLIC mm Hg (upper number)	DIASTOLIC mm Hg (lower number)
NORMAL	LESS THAN 120	and LESS THAN 80
ELEVATED	120–129	and LESS THAN 80
HIGH BLOOD PRESSURE (HYPERTENSION) STAGE 1	130–139	or 80–89
HIGH BLOOD PRESSURE (HYPERTENSION) STAGE 2	140 OR HIGHER	or 90 OR HIGHER
HYPERTENSION CRISIS (consult your doctor immediately)	HIGHER THAN 180	and/or HIGHER THAN 120

\*Source: American Heart Association® (AHA) 2017

### Display Explanation Symbols

	<b>Low Battery Symbol:</b> Appears when battery power is excessively low, or the polarities are reversed.
	<b>Pulse Symbol:</b> Once pulse is detected, the symbol flashes with each pulse beat. Shows the heart rate per minute.
	<b>Irregular Heartbeat Detector:</b> This symbol appears when an irregular heart beat was detected. See Irregular Heartbeat Detector section for more information. If symbol appears frequently, please consult your physician.
	<b>Risk Category Index:</b> The arrow points out the specific Risk Category that your measurement reading falls in. See Risk Category Index Section for more information.
	<b>Excessive Body Motion Detector:</b> This symbol appears if body movement is detected during measurement, especially, the movement on the arm the blood pressure monitor is worn on. If the cuff is worn improperly, or the shape of the upper arm is unusual (for example, the circumference of the upper arm differs largely from the circumference of the forearm), excessive gap might exist between the arm cuff and the arm.
	<b>NOTE:</b> The measured blood pressure reading may not be accurate if the icon is displayed.
	This symbol appears when the LCD displays the average value of the last 3 readings.

- Clean contacts on battery and in battery compartment with a soft, dry cloth each time you install batteries.
- Batteries are hazardous waste. DO NOT dispose of them together with the household garbage.**
- DO NOT dispose of batteries in fire. Batteries may explode or leak.**
- Recycle or dispose of properly in accordance with local, state, province, and country regulations.

### DATE & TIME SET PROCEDURE

To set the date and time when monitor is ON, press the **START | STOP** button for 3 seconds. The display will show a blinking number showing the MONTH. Change the MONTH by pressing the **M/+** button. Each press will increase the number by one in a cycling manner. Press the **START | STOP** button again to confirm the entry, and the screen will show a blinking number representing the DAY. Change the YEAR, HOUR and MINUTE as described in Step 2 above, using the **M/+** button to change the numbers and the **START | STOP** button to confirm the entries. When settings are done, press **START | STOP** button to confirm the entries. The device is ready to measure.

### MEASUREMENT PROCEDURE

#### Important Notes:

Blood pressure changes with every heartbeat and is in constant fluctuation throughout the day. Blood pressure measurement can be affected by the position of the user, his or her physiologic condition and other factors.

#### Before Measurement:


- To help ensure an accurate reading, follow these instructions before taking a measurement:
  - Wait 30 minutes after exercising, eating, drinking beverages with alcohol or caffeine, or smoking to measure blood pressure.
  - Sit quietly and rest for 5 minutes.
  - Empty your bladder, as a full one can temporarily raise blood pressure.
  - Stress raises blood pressure. Avoid taking measurements during stressful times.
  - Take your reading in a comfortable environment as measurements can be affected by hot or cold temperatures. Take your blood pressure at normal body temperature.

#### During Measurement:

- DO NOT talk or move your arm or hand muscles.
- DO NOT cross your legs. Sit with feet flat on the floor.
- DO NOT touch cuff or monitor during measurement.

#### Taking a Measurement:

- With the cuff wrapped around your upper arm, press **START | STOP** button to start measurement. All display symbols appear on the screen for 1.5 seconds.


#### USING THE UPPER ARM CUFF

##### Very Important: Proper cuff size is critical to an accurate measurement.

This monitor comes with a one size fits most upper arm cuff that fits arm sizes 9" ~17" (23cm ~ 43 cm).

The appropriate cuff is suitable for your use if the colored arrow falls within the solid color line as shown below.

1. If the cuff is not assembled, pass the end of the cuff furthest away from the tubing through the metal D-ring in order to form a loop. The smooth side without the felt material should be on the inside of the cuff loop.

2. Plug the cuff tube into the left side of the unit. You may hear a soft click once fully inserted.

3. Open the arm cuff so that your arm may be placed through it.

4. Wrap the cuff on a bare arm or over thin clothing. Thick clothing or a rolled up sleeve will cause inaccurate blood pressure measurements.

5. The bottom edge of the cuff should be approximately 1 inch (2-3 cm) above the inner elbow.


6. Center tube over middle of arm.

7. Pull the end of the cuff so that it tightens evenly around your arm. Press the hook and loop material together to secure. Allow room for 2 fingers to fit between the cuff and your arm. Please make sure the cuff does not slip during measurement, and the arrow falls within the proper fit range.

8. Lay your arm on a table (palm upward) so the cuff is at the same height as your heart. Make sure the tube is not kinked.

- After all symbols disappear, the display will show "00". The monitor is "Ready to Measure" and will automatically inflate to the level that is right for you.
- As the cuff inflates, the monitor automatically determines your ideal inflation level.

This monitor features Smart Measure® inflation technology, and detects your blood pressure and pulse rate during inflation. The Pulse Symbol () flashes at every heartbeat. Remain still and do not move until the entire measurement process is completed. The device will detect your pulse and determine the blood pressure.

4. After the monitor has determined your blood pressure and heart rate, the cuff automatically deflates. Your systolic rate, diastolic rate, heart rate and corresponding Risk Category Indicator, Irregular Heartbeat Detector and Excessive Body Motion Detector (if any) are displayed with date and time for 1 minute and save results to memory automatically.

5. Device automatically shuts off if no operation over 1 minute.

- If the cuff does not stop inflating, remove the cuff at once.
- To stop measurement, press **START | STOP** button. The cuff will deflate immediately after the button is pressed.


### USING THE MEMORY FUNCTION

The monitor can store up to 60 memory sets, and automatically replace the oldest data with new one.

After each measurement, the systolic and diastolic pressure, heart rate and corresponding Risk Category Indicator, Irregular Heartbeat Detector and Excessive Body Motion Detector (if any) with date and time will be automatically stored.

### RECALLING VALUES FROM MEMORY

- Press **M/+** button to access memory. If there is no stored memory, nothing except month, date, and time will appear on the display.
- If yes, the first reading will be the average of the last 3 measurements.


- The monitor will first display the calculated average applied to the last 3 memories.
- Every press of the button will recall a previous reading. The last reading will be recalled first.

4. To exit Memory Mode, press the **START | STOP** button


### CLEARING VALUES FROM MEMORY

- Press the **M/+** button to enter into memory recall mode.
- Press and hold the **START | STOP** and **M/+** buttons at the same time and the data for the selected user will be erased automatically.
- To confirm deletion, press the **M/+** button and no data should appear.

#### NOTE:

- This device should not be used when your arm is injured or wounded.
- If it is not possible to fit the cuff on your left arm, it can also be placed on your right arm. However, all measurements should be made using the same arm.
- To use on the right arm, you must position the artery symbol " over the main artery. Locate the main artery by pressing with two fingers approximately 1" (2cm) above the bend of your elbow on the inside of your right arm. Identify where the pulse can be felt the strongest. This is your main artery.

NOTE: Once deleted, your readings cannot be restored.

## IMPORTANT NOTES REGARDING YOUR BLOOD PRESSURE MEASUREMENT

- It is suggested that you take your measurements at the same time each day and use the same arm for consistency.
- Users should wait a minimum of 5 minutes before taking additional measurements. More time may be necessary depending upon your physiology.
- The measurement results that users receive are for reference only. If users have any blood pressure concerns, please consult a physician.
- Once inflation reaches 300 mmHg, the unit will deflate automatically for safety reasons.
- This product is not suitable for people with arrhythmias.
- This device may have difficulty determining the proper blood pressure for users with irregular heartbeat, diabetes, poor circulation of blood, kidney problems, or for users who have suffered a stroke.

## CARE, MAINTENANCE & CLEANING

- Clean the blood pressure monitor body and cuff carefully with a slightly damp, soft cloth. DO NOT press. DO NOT wash cuff or use chemical cleaner on it. Never use thinner, alcohol or petrol (gasoline) as cleaner.
- Make sure the cuff is completely dry before using.
- Leaky batteries can damage the unit. Remove the batteries when the unit will not be used for a long time.
- Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.
- If the unit is stored near freezing, allow it to acclimate to room temperature before use.
- This blood pressure monitor is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems with this device, please contact Homedics Consumer Relations (contact information can be found on the *Warranty* section).
- DO NOT immerse the unit in water as this will result in damage to the unit.
- DO NOT subject the monitor or cuff to extreme temperatures, humidity, moisture, or direct sunlight. Protect from dust.
- DO NOT fold the cuff and tubing tightly.
- DO NOT disassemble the monitor or cuff. If in need of repair, refer to the warranty section of this manual.
- DO NOT subject the monitor to extreme shocks (do not drop on floor).
- DO NOT inflate the cuff unless wrapped around arm.
- DO NOT wrap the cuff around body parts other than your arm
- DO NOT drop or insert any object into any opening or hose.
- To avoid accidental strangulation keep this product away from children.
- DO NOT drap tube around neck.
- This monitor may not meet its performance specifications if stored or used outside of these temperature and humidity ranges:

## STORAGE/TRANSPORTATION ENVIRONMENT

Temperature: -13°F ~ 158°F (-25°C ~ 70°C)

Humidity: less than 93% RH

## OPERATING ENVIRONMENT

Temperature: 41°F ~ 104°F (5°C ~ 40°C)

Humidity: 15% ~ 93% RH

Atmospheric Pressure: 700hPa~1060hPa

## POTENTIAL FOR ELECTROMAGNETIC INTERFERENCE

To avoid inaccurate results caused by electromagnetic interference between electrical and electronic equipment, do not use the device near a cell phone or microwave oven. For most wireless communication devices, it is recommended to maintain a distance of 10.8 feet (3.3m) in order to avoid electromagnetic interference.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

## FEDERAL COMMUNICATIONS COMMISSION COMPLIANCE STATEMENT

Changes or modifications to this equipment not expressly approved by the manufacturer could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

## ELECTROMAGNETIC COMPATIBILITY (EMC)

### • Guidance and manufacturer's declaration - electromagnetic emissions

The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	Group 1	RF energy is used only to maintain device's operation. Therefore, its RF emissions are so low that it's not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments, and those directly connected to the public low-voltage power supplies network that supplies buildings used for domestic purposes.
Harmonic Emissions IEX 61000-3-2	Not Applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not Applicable	

### • Guidance and manufacturer's declaration - electromagnetic emissions

The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment-Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	In the case of air discharge testing, the climatic conditions shall be within the following ranges: Ambient Temperature: 15°C ~ 35°C, Relative Humidity: 30% ~ 60%
Power frequency (50 or 60 Hz) Magnetic field IEC 61000-4-8	30 A/m 50 or 60 Hz	30 A/m 50 or 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

### • Guidance and manufacturer's declaration - electromagnetic immunity

The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment-Guidance
Conducted RF IEC 61000-4-6	3V rms At 0.15-80 MHz 6V rms At ISM & Radio Amateur Freq.	Not Applicable	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment-Guidance
Radiated RF IEC 61000-4-3	10 V/m at 80-2700 MHz AM Modulation And 9-28V / m at 385-6000 MHz, Pulse Mode and other Modulation. The system shall be tested as specified in IEC60601-1-2 Table 9 for proximity fields from RF wireless communications equipment using the test methods specified in IEC 61000-4-3	10 V/m at 80-2700 MHz AM Modulation And 9-28V / m at 385-6000 MHz, Pulse Mode and other Modulation. The system shall be tested as specified in IEC60601-1-2 Table 9 for proximity fields from RF wireless communications equipment using the test methods specified in IEC 61000-4-3	Recommended separation distance considering to reduce the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation: $E=6/d^{1/p}$ where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVELS in V/m.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment-Guidance
			Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,* should be less than the compliance level in each frequency range <sup>b</sup> . Interference may occur in the vicinity of equipment marked with the following symbol:

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment-Guidance
			NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment-Guidance
			NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment-Guidance
			a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment-Guidance
			b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

## Test specifications for enclosure port immunity to RF wireless communications equipment.

Test frequency (MHz)	Modulation	Immunity Test Level (V/m)
385	Pulse modulation 18 Hz <sup>a)</sup>	27
450	FM± 5 kHz deviation 1kHz sine <sup>b)</sup>	28
710		
745	Pulse modulation 217 Hz <sup>a)</sup>	9
780		
810	Pulse modulation 18 Hz <sup>a)</sup>	28
870		
930	Pulse modulation 217 Hz <sup>a)</sup>	28
1720		
1845	Pulse modulation 217 Hz <sup>a)</sup>	28
1970		
2450	Pulse modulation 217 Hz <sup>a)</sup>	28
5240		
5500	Pulse modulation 217 Hz <sup>a)</sup>	9
5785		

**NOTE:** If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m, The 1 m test distance is permitted by IEC 61000-4-3.

a) The carrier shall be modulated using a 50% duty cycle square wave signal.  
b) As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

## TROUBLESHOOTING

If any abnormality arises during use, please check the following points:

Symptoms	Possible Causes	Correction
Unit does not turn on when the START/STOP button is pushed.	Batteries have run down.	Replace them with four new AAA (LRO3) alkaline batteries.
Battery polarities (+/-) have been positioned incorrectly.	Re-insert the batteries in the correct positions.	

EE measurement error symbol shown on display or the blood pressure value is displayed excessively low (or high).	The Upper Arm Cuff has been positioned on the arm incorrectly.	Re-wrap the cuff properly so that it is positioned correctly. Take new measurement. Keep arm steady during measurement.
	Did you talk or move during measurement?	

E1 error symbol shown on display.	Air circuit abnormality. Cuff tube may not be plugged into monitor correctly.	Check cuff connection. See Using the Upper Arm Cuff section. Take new measurement.
	Inflation pressure exceeding 300 mmHg.	Switch the unit off and then turn back ON. Take new measurement.

E2 error symbol shown on display.	Error determining measurement data.	
	Cuff is worn improperly, or the shape of the upper arm is unusual (for example, the circumference of the upper arm differs largely from the circumference of the forearm), excessive gap might exist between the arm cuff and the arm.	Re-wrap the cuff properly so that it is positioned correctly. Take new measurement. If you have any question about the cuff wearing and/or measurement result, please consult your healthcare professional.

Excessive Body Motion Detector	Body movement during measurement, especially, the movement on the arm the blood pressure
--------------------------------	--

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