600 series automatic upper arm blood pressure monitor

BPA-H150BT
Doc # L-03633, Rev.1
5 year limited warranty

IMPORTANT PRODUCT NOTICES AND SAFETY INSTRUCTIONS

When using your blood pressure monitor, basic precautions should always be followed. Please read and follow all instructions and warnings before using this product.

SAVE THESE INSTRUCTIONS FOR FUTURE REFERENCE.

This device is intended for over-the-counter home use in adults 18 years and older with an arm circumference ranging from 9 inches to 17 inches (approx. 23cm to 43cm). The device can accurately measure blood pressure in pregnant patients including those with known or suspected in pre-eclampsia condition. This device features a built-in "Bluetooth® data transmission" function, which enables the device to automatically transmit measuring results to a paired Bluetooth®-enabled device. This function allows users to synchronize to current date and time, and check the battery status on the Homedics Health+ app.

• Please note that this is a home healthcare product only and this manual and device are not intended to serve as a substitute for the advice of a physician or medical professional.

• This device uses the oscillometric method to automatically measure systolic and diastolic blood pressure, as well as heart rate. All values can be read out on one LCD Display.

• DO NOT use this device for diagnosis or treatment of any health problem or disease. Measurement results are for reference only. Consult a healthcare professional for interpretation of pressure measurements. Contact your physician if you have or suspect any medical problem. Do not change your medications without the advice of your physician or healthcare professional.

• We recommend you using the same arm (preferably the left arm) and measuring around the same time each day.

• Perform measurements in a quiet and relaxed environment at room temperature.

• If you are pregnant, you should pay more attention to your blood pressure changes, because during this time, it may change drastically.

• This monitor is clinically validated for use in pregnancy and pre-eclampsia.

• When you detect unusual readings in pregnancy, you should measure again after taking some rest. If the reading is still abnormal, consult your doctor or gynecologist.

• Proper cuff size is critical for accurate measurements. Follow the instructions in this manual and printed on the arm cuff to ensure the arm cuff is used properly.

• This product is not suitable for:

- People with arrhythmias
- People undergoing intravenous injection on any limb
- People currently in a dialysis treatment

• Those who have had a mastectomy surgery (especially when lymph nodes have been removed), it is recommended to take measurements on the unaffected side.

• This device may have difficulty determining proper blood pressure for users with irregular heartbeat, atherosclerosis, diabetes, liver disease, kidney disease, poor circulation of the blood or for users who have suffered a stroke. Please consult your healthcare professional before using this device.

• Excessive use may result in blood flow interference, which is likely to cause uncomfortable sensations, such as partial subcutaneous hemorrhage, or temporary numbness to your arm. In general, these symptoms should not last long. However, if you continue to experience these sensations, please seek advice from a medical professional.

• When used along with other medical electronic equipment on the same arm, pressurization of the cuff may cause the other devices to temporarily malfunction.

• The pulse display is not suitable for checking the frequency of heart pacemakers.

• Electromagnetic interference: The device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g., mobile telephones, microwave ovens). These may lead to temporary impairment of measurement accuracy.

• Consider the electromagnetic compatibility of the device (ex. power disturbance, radio frequency interference etc.). Please use this device indoors only.

• Use blood pressure monitor only for its intended use.

• DO NOT wrap the cuff around body parts other than your arm.

• Not for use by or on persons under the age of 18.

• DO NOT use this device on infants, children, or those who cannot express their own intention.

• DO NOT plug or unplug the adapter power cord with wet hands.

• Please use only the AC adapter indicated for use with this monitor or 1.5V AAA alkaline batteries for power supply.

• Do not take a measurement in a low (less than 41°F/5°C) and high (more than 104°F/40°C) temperature, nor in a place outside humidity ranges (15% ~ 93% R.H.), and atmospheric pressure ranges (700 ~ 1060 hPa), or you may get inaccurate readings.

• Please rest for at least 5-10 minutes before taking a measurement.

• To allow your blood vessels to return to the condition prior to taking the measurement, please wait at least 3-5 minutes between measurements. You may need to adjust the wait time according to your personal physiological situation.

• Wait 30-45 minutes before measurement if you've just consumed caffeinated beverages or smoked cigarettes.

• The applied part is the cuff.

• The patient is the intended operator.

Blood pressure measurements determined with this device are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultatory method, within the accuracy limits prescribed by the American National Standard, Manual, electronic, or automated sphygmomanometers.

In the event that the device needs to be checked for calibration, contact consumer relations.

ABOUT BLOOD PRESSURE

What is blood pressure?

Blood pressure is the measurement of the force exerted on the artery walls while blood flows through the arteries. The pressure measured when the heart contracts and sends blood out of the heart is systolic (highest) blood pressure. The pressure measured when the heart dilates with blood flowing back into the heart is diastolic (lowest) blood pressure.

Why measure your blood pressure?

Among today's various health problems, those associated with high blood pressure are very common. High blood pressure dangerously correlates with cardiovascular diseases. Therefore, blood pressure monitoring is important for identifying those at risk.

Why do my readings vary?

Blood pressure is a body parameter that is subject to normal variations throughout the day. A single reading that is different from yours or your doctor's readings are not necessarily inaccurate. The average of several readings, taken under similar conditions, using the same arm is preferred for accurate blood pressure readings.

Why Are My Readings Different Than Those Taken at My Doctor's Office?
Many experience a phenomenon called "White Coat Hypertension" when measured by a doctor. White Coat Hypertension refers to blood pressure that rises above its usual level when measured in a clinical setting, such as a doctor's office.

BLOOD PRESSURE STANDARD

The table below contains defined levels for hypertension that are publicly available from the American Heart Association® (AHA 2017). Users can compare their own blood pressure readings against these defined levels to determine if they may be potentially at increased risk.

BLOOD PRESSURE CATEGORY	SYSTOLIC mm Hg (upper number)	DIASTOLIC mm Hg (lower number)
NORMAL	LESS THAN 120	and LESS THAN 80
ELEVATED	120-129	and LESS THAN 80
HIGH BLOOD PRESSURE (HYPERTENSION) STAGE 1	130-139	or 80-89
HIGH BLOOD PRESSURE (HYPERTENSION) STAGE 2	140 OR HIGHER	or 90 OR HIGHER
HYPERTENSION CRISIS (consult your doctor immediately)	HIGHER THAN 180	and/or HIGHER THAN 120

*Source: American Heart Association® (AHA) 2017

DISPLAY SYMBOL EXPLANATION

	Low Battery Symbol: Appears when battery power is excessively low, or the polarities are reversed.
	Pulse Symbol: Once pulse is detected, the symbol flashes with each pulse beat. Shows the heart rate per minute.
	Irregular Heartbeat Detector: This symbol appears when an irregular heart beat was detected. See Irregular Heartbeat Detector section for more information. If symbol appears frequently, please consult your physician.
	Risk Category Index: The arrow points out the specific Risk Category that your measurement falls in. See Risk Category Index Section for more information.

	Excessive Body Motion Detector: Displayed if body movement is detected during measurement, especially, the movement on the arm the blood pressure monitor is worn on. If the cuff is worn improperly, or the shape of the upper arm is unusual (for example, the circumference of the upper arm differs largely from the circumference of the forearm), excessive gap might exist between the arm cuff and the arm.
	Cuff Fitness Detection: Displayed if the cuff was wrapped incorrectly, which is too tight or too loose. This is the function aid in detecting if the cuff is wrapped properly.

	Bluetooth® Symbol: Appears when Bluetooth is active, and when data is transmitting to your mobile device.
	Average of the Last 3 Measurements: This symbol appears when the LCD displays the average value of last 3 readings.
	AM/PM Averaging: This symbol appears when LCD displays an average from the last 3 morning or last 3 evening measurements.

	Using the AC/DC Adapter (NOT INCLUDED): This monitor is designed to operate with an AC/DC adapter, sold separately. If you have any additional questions, or need help on purchasing, please contact Homedics Consumer Relations, see warranty section.
	• Connect the AC/DC adapter with the DC jack on the top of the unit. • Use only the AC/DC adapter specified.

	Measurement Error: Wrap the cuff correctly and keep arm steady during measurement. Take new measurement.
	Air Circuit Abnormality: Check cuff connection. Take new measurement.
	Pressure Exceeding 300 mmHg: Turn the unit off to clear, then take new measurement.
	Error Determining Measurement Data: Ensure the cuff is worn properly, do not move or talk. Take new measurement.
	Data Transmission Error: Will occur after 45 seconds, if the monitor cannot connect to the mobile device to transmit data. Make sure Bluetooth® is ON. Press M+/<#x2f;> button for 3 seconds to start Bluetooth® function.
	System Error: Contact Consumer Relations

	AC/DC Adapter: Input: 100-240V, 50/60 Hz, 0.2 A Output: 5V, DC, 1.2A
--	--

DATE & TIME SET PROCEDURE

Date & Time can be set by two methods, either sync automatically using your Bluetooth® mobile device, or set manually using the and buttons on the monitor.

If setting the Date & Time by syncing with your mobile device, it is important that this is done prior to taking any measurements to ensure the date and time are accurate.

To set using your device with Bluetooth® wireless technology:


1. Before taking the first measurement, open the Homedics Health+ app on your device. If you have not already installed the Homedics Health+ app onto your device, it is available to download on the App Store and on Google Play™. Make sure you have the app downloaded and open on your mobile device before trying to sync the Date & Time. Make sure the Bluetooth® is on your mobile device and your App is open.
2. Follow the prompts in the app to select your device to sync with the app.
3. Once the Bluetooth® connection is established with the device follow the prompt on the app to sync the date & time on the device.

NOTE:

- Once the Date & Time have been successfully synced, future readings will automatically have the correct date and time.
- If you are having trouble automatically syncing the Date & Time, make sure the Bluetooth® is ON on your device and ON on the blood pressure monitor.

To Set Manually:

1. To set the date and time when monitor is ON, press and hold the for 3 seconds.
2. The display will show a blinking number showing the MONTH. Change the MONTH by pressing the button. Each press will increase the number by one in a cycling manner. Press the again to confirm the entry, and the screen will show a blinking number representing the DAY.
3. Change the YEAR, HOUR and MINUTE, and select OFF, or leave ON Bluetooth function, as described above, using the to change the numbers and the to confirm the entries.

Repeat process to set DAY, YEAR, HOUR and MINUTE, and select OFF, or leave ON Bluetooth function.

Turning the Bluetooth® Function ON/OFF:

The LCD screen shows 'ON' and 'OFF' with arrows indicating a transition. Below the screen, it says 'Press M+/<#x2f;> button'.

Bluetooth® function ON

Bluetooth® function OFF

Bluetooth® function ON

IMPORTANT NOTES REGARDING YOUR BLOOD PRESSURE MEASUREMENT

- It is suggested that you take your measurements at the same time each day and use the same arm for consistency.
- Users should wait a minimum of 5 minutes before taking additional measurements. More time may be necessary depending upon your physiology.
- The measurement results that users receive are for reference only. If users have any blood pressure concerns, please consult a physician.
- Once inflation reaches 300 mmHg, the unit will deflate automatically for safety reasons.
- This product is not suitable for people with arrhythmias.
- This device may have difficulty determining the proper blood pressure for users with irregular heartbeat, diabetes, poor circulation of blood, kidney problems, or for users who have suffered a stroke.

CARE, MAINTENANCE & CLEANING

- Clean the blood pressure monitor body and cuff carefully with a slightly damp, soft cloth. Do not press. Do not wash cuff or use chemical cleaner on it. Never use thinner, alcohol or petrol (gasoline) as cleaner.
- Make sure the cuff is completely dry before using.
- Leaky batteries can damage the unit. Remove the batteries when the unit will not be used for a long time.
- Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.
- If the unit is stored near freezing, allow it to acclimate to room temperature before use.
- This blood pressure monitor is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems with this device, please contact Homedics Consumer Relations (contact information can be found on the Warranty section).
- DO NOT immerse the unit in water as this will result in damage to the unit.
- DO NOT subject the monitor or cuff to extreme temperatures, humidity, moisture, or direct sunlight. Protect from dust.
- DO NOT fold the cuff and tubing tightly.
- DO NOT disassemble the monitor or cuff. If in need of repair, refer to the warranty section of this manual.
- DO NOT subject the monitor to extreme shocks (do not drop on floor).
- DO NOT inflate the cuff unless wrapped around arm.
- DO NOT wrap the cuff around body parts other than your arm.
- DO NOT drop or insert any object into any opening or hose.
- To avoid accidental strangulation keep this product away from children.
- DO NOT drape tube around neck.
- This monitor may not meet its performance specifications if stored or used outside of these temperature and humidity ranges:

STORAGE/TRANSPORTATION ENVIRONMENT

Temperature: -13°F ~ 158°F (-25°C ~ 70°C)

Humidity: less than 93% RH

OPERATING ENVIRONMENT

Temperature: 41°F ~ 104°F (5°C ~ 40°C)

Humidity: 15% ~ 93% RH

Atmospheric Pressure: 700hPa-1060hPa

POTENTIAL FOR ELECTROMAGNETIC INTERFERENCE

To avoid inaccurate results caused by electromagnetic interference between electrical and electronic equipment, do not use the device near a cell phone or microwave oven. For most wireless communication devices, it is recommended to maintain a distance of 10.8 feet (3.3m) in order to avoid electromagnetic interference.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

FEDERAL COMMUNICATIONS COMMISSION COMPLIANCE STATEMENT
Changes or modifications to this equipment not expressly approved by the manufacturer could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

ELECTROMAGNETIC COMPATIBILITY (EMC)

• Guidance and manufacturer's declaration - electromagnetic emissions

The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	Group 1	RF energy is used only to maintain device's operation. Therefore, its RF emissions are so low that it's not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments, and those directly connected to the public low-voltage power supplies network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

• Guidance and manufacturer's declaration - electromagnetic emissions

The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment-Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	In the case of air discharge testing, the climatic conditions shall be within the following ranges: Ambient Temperature: 15°C ~ 35°C, Relative Humidity: 30% ~ 60%
Power frequency (50 or 60 Hz)	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Magnetic field IEC 61000-4-8	50 or 60 Hz	50 or 60 Hz	
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	AC Power port ± 1 kV line to line	AC Power port ± 1 kV line to line	Mains power quality should be that of a typical commercial or hospital environment.
Interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT: 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 1 cycles 70% UT; 25/30 cycles 0% UT; 250/300 cycles	0% UT: 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 1 cycles 70% UT; 25 cycles 0% UT; 250 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.

• Guidance and manufacturer's declaration - electromagnetic immunity

The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment-Guidance
Conducted RF IEC 61000-4-6	3V rms At 0.15-80 MHz 6V rms At ISM & Radio Amateur Freq.	3V rms At 0.15-80 MHz 6V rms At ISM & Radio Amateur Freq.	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

NOTE: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1m, The 1m test distance is permitted by IEC 61000-4-3.

a) The carrier shall be modulated using a 50% duty cycle square wave signal.

b) As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Test specifications for enclosure port immunity to RF wireless communications equipment.

Test frequency (MHz)	Modulation	Immunity Test Level (V/m)
385	Pulse modulation 18 Hz ^{a)}	27
450	FM± 5 kHz deviation 1kHz sine ^{b)}	28
710		
745	Pulse modulation 217 Hz ^{a)}	9
780		
810		
870	Pulse modulation 18 Hz ^{a)}	28
930		
1720		
1845	Pulse modulation 217 Hz ^{a)}	28
1970		
2450	Pulse modulation 217 Hz ^{a)}	28
5240		
5500	Pulse modulation 217 Hz ^{a)}	9
5785		

NOTE: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1m, The 1m test distance is permitted by IEC 61000-4-3.

a) The carrier shall be modulated using a 50% duty cycle square wave signal.

b) As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

TROUBLESHOOTING

If any abnormality arises during use, please check the following points:

Symptoms	Possible Causes	Correction
Unit does not turn on when the START/STOP button is pushed.	Batteries have run down. Battery polarities (+/-) have been positioned incorrectly.	Replace them with four new AAA (LR03) alkaline batteries. Re-insert the batteries in the correct positions.
EE measurement error symbol shown on display or the blood pressure value is displayed excessively low (or high).	The Upper Arm Cuff has been positioned on the arm incorrectly. Did you talk or move during measurement? Shaking of the arm with the cuff on.	Re-wrap the cuff properly so that it is positioned correctly. Take new measurement. Keep arm steady during measurement.
E1 error symbol shown on display.	Air circuit abnormality. Cuff tube may not be plugged into monitor correctly.	Check cuff connection. See Using the Upper Arm Cuff section. Take new measurement.
E2 error symbol shown on display.	Inflation pressure exceeding 300 mmHg.	Switch the unit off and then turn back ON. Take new measurement.
E3 error symbol shown on display.	Error determining measurement data. Cuff is worn improperly, or the shape of the upper arm is unusual (for example, the circumference of the upper arm differs largely from the circumference of the forearm), excessive gap might exist between the arm cuff and the arm.	Re-wrap the cuff properly so that it is positioned correctly. Take new measurement. If you have any question about the cuff wearing and/or measurement result, please consult your healthcare professional.

Bluetooth® function is turned OFF, on your mobile device, and try again.

Turn the Bluetooth® function ON, on your blood pressure monitor, then try again.

Make sure Bluetooth® is ON, on the mobile device and the blood pressure monitor, and try again.

Requires compatible mobile device. This blood pressure monitor is compatible with mobile devices that are running Bluetooth® 5.0 for Android 6.0 or above, Bluetooth® 5.0 for iOS 9.0 or above.

Mobile device does not pair with the blood pressure monitor.

The App on the mobile device is not open.

The blood pressure monitor and mobile device are out of the acceptable range of 32 feet (10 meters).

Make sure your last reading is stored in memory and the App is open and try again. See Transmitting Readings to Your Mobile Device section.

Unexpected loss of electrical/mechanical integrity.

Remove batteries, re-insert, and try again.

Return the device to your local distributor or importer.

Body movement during measurement, especially, the movement on the arm the blood pressure monitor is worn on. E.g. talking, moving or shaking of the arm with the cuff on while measurement.

Measure again. Keep arm steady during measurement.

Cuff is worn improperly, or the shape of the upper arm is unusual (for example, the circumference of the upper arm differs largely from the circumference of the forearm), excessive gap might exist between the arm cuff and the arm.

Wrap the cuff properly and keep steady. Measure again.

If you have any question about the cuff wearing and/or measurement result, please consult your healthcare professional.

Please reference "applying the Cuff" section to wrap the cuff correctly.

Please re-pair the BPM and Bluetooth-enabled device with each other.

See the "Bluetooth® Transmission" section to turn on Bluetooth® function.

Please make sure the acceptable distance (<10 meters) with each other.

Please refer to RF Specification section.

Use non-Bluetooth-enabled device.

Unexpected loss of electrical/mechanical integrity.

Return the device to your local distributor or importer.

Measure again. If error persists, contact Consumer Relations.

System Error.

Turn Bluetooth® ON, on your mobile device, and try again.

Turn the Bluetooth® function ON, on your blood pressure monitor, then try again.

Make sure Bluetooth® is ON, on the mobile device and the blood pressure monitor, and try again.

Requires compatible mobile device. This blood pressure monitor is compatible with mobile devices that are running Bluetooth® 5.0 for Android 6.0 or above, Bluetooth® 5.0 for iOS 9.0 or above.

Mobile device does not pair with the blood pressure monitor.

The App on the mobile device is not open.

The blood pressure monitor and mobile device are out of the acceptable range of 32 feet (10 meters).

Make sure your last reading is stored in memory and the App is open and try again. See Transmitting Readings to Your Mobile Device section.

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BPA-HIS60T
Doc # L-03633, Rev.1
5 year
limited warranty

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