800 series automatic upper arm blood pressure monitor

INSTRUCTION MANUAL

BPA-H190BT
Doc # L-03645, Rev.1
5 year limited warranty

IMPORTANT PRODUCT NOTICES AND SAFETY INSTRUCTIONS

When using your blood pressure monitor, basic precautions should always be followed. Please read and follow all instructions and warnings before using this product.

SAVE THESE INSTRUCTIONS FOR FUTURE REFERENCE.

This device is intended for over-the-counter home use in adults 18 years and older with an arm circumference ranging from 9 inches to 17 inches (approx. 23cm to 43cm). The device can accurately measure blood pressure in pregnant patients including those with known or suspected in pre-eclampsia condition. This device features a built-in "Bluetooth® data transmission" function, which enables the device to automatically transmit measuring results to a paired Bluetooth®-enabled device. This function allows users to synchronize to current date and time, and check the battery status on the Homedics Health+ app.

• Please note that this is a home healthcare product only and this manual and device are not intended to serve as a substitute for the advice of a physician or medical professional.

• This device uses the oscillometric method to automatically measure systolic and diastolic blood pressure, as well as heart rate. All values can be read out on one LCD Display.

• DO NOT use this device for diagnosis or treatment of any health problem or disease. Measurement results are for reference only. Consult a healthcare professional for interpretation of pressure measurements. Contact your physician if you have or suspect any medical problem. Do not change your medications without the advice of your physician or healthcare professional.

• We recommend you using the same arm (preferably the left arm) and measuring around the same time each day.

• Perform measurements in a quiet and relaxed environment at room temperature.

• When you detect unusual readings in pregnancy, you should measure again after changes, because during this time, it may change drastically.

• This monitor is clinically validated for use in pregnancy and pre-eclampsia.

• When you detect unusual readings in pregnancy, you should measure again after some rest. If the reading is still abnormal, consult your doctor or gynecologist.

• Proper cuff size is critical for accurate measurements. Follow the instructions in this manual and printed on the arm cuff to ensure the arm cuff is used properly.

• This product is not suitable for:

- People with arrhythmias
- People undergoing intravenous injection on any limb
- People currently in dialysis treatment

• Those who have had a mastectomy surgery (especially when lymph nodes have been removed), it is recommended to take measurements on the unaffected side.

• This device may have difficulty determining proper blood pressure for users with irregular heartbeat, atherosclerosis, diabetes, liver disease, kidney disease, poor circulation of the blood or for users who have suffered a stroke. Please consult your healthcare professional before using this device.

• Excessive use may result in blood flow interference, which is likely to cause uncomfortable sensations, such as partial subcutaneous hemorrhage, or temporary numbness to your arm. In general, these symptoms should not last long. However, if you continue to experience these sensations, please seek advice from a medical professional.

• When used along with other medical electronic equipment on the same arm, pressurization of the cuff may cause the other devices to temporarily malfunction.

• The pulse display is not suitable for checking the frequency of heart pacemakers.

• Electromagnetic interference: The device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens). These may lead to temporary impairment of measurement accuracy.

• Consider the electromagnetic compatibility of the device (ex. power disturbance, radio frequency interference etc.). Please use this device indoors only.

• Use blood pressure monitor only for its intended use.

• DO NOT wrap the cuff around body parts other than your arm.

• Not for use by or on persons under the age of 18.

• DO NOT use this device on infants, children, or those who cannot express their own intention.

• DO NOT plug or unplug the adapter power cord with wet hands.

• Please use only the AC adapter indicated for use with this monitor or 1.5V AAA alkaline batteries for power supply.

• Please rest for at least 5-10 minutes before taking a measurement.

• To allow your blood vessels to return to the condition prior to taking the measurement, please wait at least 3-5 minutes between measurements. You may need to adjust the wait time according to your personal physiological situation.

• Wait 30-45 minutes before measurement if you've just consumed caffeinated beverages or smoked cigarettes.

• The applied part is the cuff.

• The patient is the intended operator, who can operate the device by him or herself, not necessarily by a physician or operator.

• Do not take a measurement in a too low (less than 41°F / 5°C) or too high (more than 104°F / 40°C) temperature, nor in a place outside humidity ranges (15% ~ 93% R.H.) and atmospheric pressure ranges (700 ~ 1060hPa), or you may get inaccurate readings.

Blood pressure measurements determined with this device are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultatory method, within the accuracy limits prescribed by the American National Standard, Manual, electronic, or automated sphygmomanometers.

In the event that the device needs to be checked for calibration, contact consumer relations.

ABOUT BLOOD PRESSURE

What is blood pressure?
Blood pressure is the measurement of the force exerted on the artery walls while blood flows through the arteries. The pressure measured when the heart contracts and sends blood out of the heart is systolic (highest) blood pressure. The pressure measured when the heart dilates with blood flowing back into the heart is called diastolic (lowest) blood pressure.

Why measure your blood pressure?

Among today's various health problems, those associated with high blood pressure are very common. High blood pressure dangerously correlates with cardiovascular diseases. Therefore, blood pressure monitoring is important for identifying those at risk.

Why do my readings vary?

Blood pressure is a body parameter that is subject to normal variations throughout the day. A single reading that is different from yours or your doctor's readings are not necessarily inaccurate. The average of several readings, taken under similar conditions, using the same arm is preferred for accurate blood pressure readings.

Why Are My Readings Different Than Those Taken at My Doctor's Office?
Having one's blood pressure measured by a doctor in a hospital or a clinic is often associated with an event called "White Coat Hypertension" where the patient becomes nervous or anxious, thus raising his/her blood pressure. There are also numerous other factors that might cause your blood pressure to be raised at a specific time of day. This is why medical physician recommend home monitoring as it is important to get readings of blood pressure during different times of the day to really get an idea of your real blood pressure.

BLOOD PRESSURE STANDARD

The table below contains defined levels for hypertension that are publicly available from the American Heart Association® (AHA 2017). Users can compare their own blood pressure readings against these defined levels to determine if they may be potentially at increased risk.

BLOOD PRESSURE CATEGORY	SYSTOLIC mm Hg (upper number)	DIASTOLIC mm Hg (lower number)	
NORMAL	LESS THAN 120	and	LESS THAN 80
ELEVATED	120–129	and	LESS THAN 80
HIGH BLOOD PRESSURE (HYPERTENSION) STAGE 1	130–139	or	80–89
HIGH BLOOD PRESSURE (HYPERTENSION) STAGE 2	140 OR HIGHER	or	90 OR HIGHER
HYPERTENSION CRISIS (consult your doctor immediately)	HIGHER THAN 180	and/or	HIGHER THAN 120

*Source: American Heart Association® (AHA) 2017


Blood pressure tends to go up and down, even in people who normally don't have high readings. Individual blood pressures vary greatly both on a daily and seasonal basis. These variations are even more pronounced in hypertensive patients. Normally the blood pressure rises while at work and is at its lowest during sleeping period. If your numbers stay above the "normal" range most of the time, you may be at increased risk and should consult your physician. Although one can easily find where their own blood pressure readings fall on this table, this monitor comes equipped with a Risk Category Index that automatically compares each reading to the defined levels and provides a helpful cue if your reading falls into one of the stages that could potentially indicate increased risk. See Risk Category Index section for more information on this feature.

HOW THIS BLOOD PRESSURE MONITOR WORKS


This monitor uses Smart Measure® inflation technology to detect your blood pressure. Once the cuff is wrapped around your upper arm, with a press of the START/STOP button, the monitor will turn on and the inflation will automatically start, creating pressure around the arteries inside the upper arm. Within the cuff is a gauge which senses the fluctuations (oscillations) in pressure. The fluctuation measured represents the degree of intensity that your arteries are contracting with each heartbeat, which is also a result of the pressure that the cuff has placed on the upper arm. The monitor measures these contractions and converts the information to a digital value. Once the measurement is complete, the cuff will automatically deflate, and the result is displayed on the screen. Please note that any muscle movement during inflation will cause measurement error. When measurement is complete, the monitor will display your systolic pressure, diastolic pressure, and pulse readings. Before measurement, it is suggested that you sit quietly for 15 minutes before measurement as measurements taken in a relaxed state have greater accuracy. The monitor automatically finds where your measurement results fall on the American Heart Association® (AHA 2017) table and provides a cue if your reading falls into one of the stages that could potentially indicate increased risk. Please refer to the Risk Category Index section for more information on this feature.

The appearance of the icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Please refer to the Irregular Heartbeat Detector section for more information on this feature.

NAME/FUNCTION OF EACH PART

DISPLAY SYMBOL EXPLANATIONS

*Source: American Heart Association® (AHA) 2017

TRIPLE-CHECK (ADVANCED) AVERAGE FUNCTION

The Advanced Average Function automatically takes and averages 3 readings in a row, with 1 minute rest intervals between each measurement. To deactivate this feature and take only one single reading, slide the switch on the side to the OFF position.

Display Explanation Symbols	
	Low Battery Symbol: Appears when battery power is excessively low, or the polarities are reversed.
	Pulse Symbol: Once pulse is detected, the symbol flashes with each pulse beat. Shows the heart rate per minute.
	Irregular Heartbeat Detector: This symbol appears when an irregular heart beat was detected. See Irregular Heartbeat Detector section for more information. If symbol appears frequently, please consult your physician.
	Risk Category Index: The arrow points out the specific Risk Category that your measurement reading falls in. See Risk Category Index Section for more information.


pulse irregularities that may cause the appearance of this icon. Therefore, it is of great importance to not move or talk during measurement. To determine the presence of an irregular heartbeat, the average of the heartbeat intervals is calculated with the first 3 normal effective heartbeat values. It is important to note that the average is not a strict mathematical averaging of all recorded intervals. At least 3 beats with 25% or greater difference from the average heartbeat interval will generate the icon on the screen.

Important Information:


This blood pressure monitor is not designed for use by people with arrhythmias nor for diagnosing or treating an arrhythmia problem. As a safeguard, we recommend that if you have arrhythmias, such as atrial or ventricular premature beats and atrial fibrillation or any other special conditions, you should check with your physician before using your blood pressure monitor.

INSTALLING BATTERIES


1. Press and lift the battery cover to open the battery compartment.

2. Install or replace 4 AAA (1.5V, LR03) alkaline batteries in the battery compartment according to the indications inside the compartment.

3. Close the battery cover by pressing the cover back into place.

Replace the batteries if:

- The low battery symbol appears.
- When any button is pressed and nothing is displayed on the screen.

NOTE:

The date and time will only need to be set manually if the monitor will not be used with the Homedics Health+ app.

- If the Homedics Health+ app is used after manually setting the date and time on the blood pressure monitor, the date and time on the mobile device will override the date and time on the blood pressure monitor.

App for Bluetooth®

Download and install the Homedics Health+ app on your smart device from App Store® and on Google Play™.

System requirement of the Bluetooth-enabled device


- Bluetooth® 4.2 for Android 6.0 or above
- Bluetooth® 4.2 for iOS 7.0 or above

BLUETOOTH® OPERATION

This monitor has Bluetooth® function ON for your convenience. This will allow your readings to automatically transmit to the Homedics Health+ app. By turning this function OFF, the measurements cannot be transmitted.

Turning the Bluetooth® Function ON/OFF:

- Press **USER/O** button to enter the Bluetooth setting. The default setting of this function is ON.
- Press **M+/S** button to turn Bluetooth function ON or OFF.

NOTES:

The date and time will only need to be set manually if the monitor will not be used with the Homedics Health+ app.


- If the Homedics Health+ app is used after manually setting the date and time on the blood pressure monitor, the date and time on the mobile device will override the date and time on the blood pressure monitor.

Using the Upper Arm Cuff

Very Important: Proper cuff size is critical to an accurate measurement.

This monitor comes with one size fits most upper arm cuff that fits arm sizes 9"-17" (23cm-43cm).


The appropriate cuff is suitable for your use if the colored arrow falls within the solid color line as shown below.


Using the USB to DC Power Cable

This monitor is designed for operating with batteries, a USB to DC cable (included), or an AC/DC adapter (not included).

1. Plug the DC cable connector into the main unit's power supply inlet as shown below:

2. Plug the USB cable connector into a USB Wall adapter, then plug the Wall adapter into an outlet, as shown below:


Use only a compatible USB/DC or AC/DC adapter with required voltage and current as indicated in this manual.

- Recommend Adapter specification, do not use otherwise: Model: FranMar International, FRM06-S05-UU


• Rating: Input: 100-240V, 50/60Hz, 0.2A - Output: 5V, DC, 1.

RECALLING VALUES FROM MEMORY

1. Press **M+/\$** button to access memory to see previous measuring results, including average of last 3 measurements, and Triple-Check (Advanced) Average measurement, AM/PM average and individual measurement. If there is no stored memory, nothing except month, date, and time will appear on the display. If yes, the first reading will be the average of the last 3 measurements.

2. Press the **M+/\$** button again, LCD displays the average of latest, last 3 AM measurements, and when pressed again will display the average of latest, last 3 PM measurements. If fewer than three measurements in AM or PM, then nothing will be displayed.

3. Press the **M+/\$** button again, LCD displays the single measurement reading or the average of the Triple-Check measurement (if any). Press the **USER/C** button to read the 3 results individually of advanced average measurement.

4. Every new press of the **M+/\$** button will recall a previous reading. The last reading will be recalled first. Press **M+/\$** button repeatedly to scroll through the previous readings stored in the memory.

5. To stop reading the memories, press **START / STOP** button to switch to sleep mode.


✓ AVERAGE		
AVERAGE MEASUREMENT 1	AVERAGE MEASUREMENT 2	AVERAGE MEASUREMENT 3
LAST 3 AVG.		
AM		An average of Last 3 Measurements
PM		An average from the last 3 morning or last 3 evening measurements.
AM		An average of Last 3 Morning Measurements (4:00AM-12:00PM)
PM		An average of Last 3 Nighttime Measurements (6:00PM-2:00AM)

NOTE: The Triple-Check (Advanced) average and its 3 individual measurements are counted as 4 readings when stored in memory.

CLEARING VALUES FROM MEMORY

1. Press the **M+/\$** button to enter into memory recall mode.

2. Press and hold the **USER/C** and **M+/\$** buttons at the same time for 3 seconds and the data for the selected user will be erased automatically.

NOTE: Once deleted, your readings cannot be restored.

IMPORTANT NOTES REGARDING YOUR BLOOD PRESSURE MEASUREMENT

- It is suggested that you take your measurements at the same time each day and use the same arm for consistency.
- Users should wait a minimum of 5 minutes before taking additional measurements. More time may be necessary depending upon your physiology.
- The measurement results that users receive are for reference only. If users have any blood pressure concerns, please consult a physician.
- Once inflation reaches 300 mmHg, the unit will deflate automatically for safety reasons.
- This product is not suitable for people with arrhythmias.
- This device may have difficulty determining the proper blood pressure for users with irregular heartbeat, diabetes, poor circulation of blood, kidney problems, or for users who have suffered a stroke.

CARE, MAINTENANCE & CLEANING

- Clean the blood pressure monitor body and cuff carefully with a slightly damp, soft cloth. DO NOT press. DO NOT wash cuff or use chemical cleaner on it. Never use detergent or any strong chemicals to clean the device.
- Make sure the cuff is completely dry before using.
- Leaky batteries can damage the unit. Remove the batteries when the unit will not be used for a long time.
- Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.
- If the unit is stored near freezing, allow it to acclimate to room temperature before use.
- This blood pressure monitor is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems with this device, please contact Homedics Consumer Relations (contact information can be found on the Warranty section).
- DO NOT immerse the unit in water as this will result in damage to the unit.
- DO NOT subject the monitor or cuff to extreme temperatures, humidity, moisture, or direct sunlight. Protect from dust.
- DO NOT fold the cuff and tubing tightly.
- DO NOT disassemble the monitor or cuff. If in need of repair, refer to the warranty section of this manual.
- DO NOT subject the monitor to extreme shocks (do not drop on floor).
- DO NOT inflate the cuff unless wrapped around arm.
- DO NOT wrap the cuff around body parts other than your arm.
- DO NOT drop or insert any object into any opening or hose.
- To avoid accidental strangulation keep this product away from children.
- DO NOT drap tube around neck.
- This monitor may not meet its performance specifications if stored or used outside of these temperature and humidity ranges:

STORAGE/TRANSPORTATION ENVIRONMENT

Temperature: -13°F ~ 158°F (-25°C ~ 70°C)

Humidity: less than 93% RH

OPERATING ENVIRONMENT

Temperature: 41°F ~ 104°F (5°C ~ 40°C)

Humidity: 15% ~ 93% RH

Atmospheric Pressure: 700hPa-1060hPa

POTENTIAL FOR ELECTROMAGNETIC INTERFERENCE

To avoid inaccurate results caused by electromagnetic interference between electrical and electronic equipment, do not use the device near a cell phone or microwave oven. For most wireless communication devices, it is recommended to maintain a distance of 10.8 feet (3.3m) in order to avoid electromagnetic interference.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

FEDERAL COMMUNICATIONS COMMISSION COMPLIANCE STATEMENT

Changes or modifications to this equipment not expressly approved by the manufacturer could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

ELECTROMAGNETIC COMPATIBILITY (EMC)

Guidance and manufacturer's declaration - electromagnetic emissions

The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	Group 1	RF energy is used only to maintain device's operation. Therefore, its RF emissions are so low that it's not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments, and those directly connected to the public low-voltage power supplies network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration - electromagnetic emissions

The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	In the case of air discharge testing, the climatic conditions shall be within the following ranges: Ambient Temperature: 15°C ~ 35°C, Relative Humidity: 30% ~ 60%
Power frequency (50 or 60 Hz) Magnetic field IEC 61000-4-8	30 A/m 50 or 60 Hz	30 A/m 50 or 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	AC Power port ± 1 kV line to line	AC Power port ± 1 kV line to line	Mains power quality should be that of a typical commercial or hospital environment.
Interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 0.5 cycles 45°, 90°, 135°, 180°, 225°, 270° and 315° 70% UT; 25/30 cycles 0% UT; 250/300 cycles	0% UT; 1 cycles 70% UT; 25 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterrupted power supply or a battery.

Guidance and manufacturer's declaration - electromagnetic immunity

The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Test frequency (MHz)	Modulation	Immunity Test Level (V/m)		
5240	Pulse modulation 217 Hz ^{a)}	9		
5500				
5785				
NOTE: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m, The 1 m test distance is permitted by IEC 61000-4-3.				
a) The carrier shall be modulated using a 50% duty cycle square wave signal.				
b) As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.				

Test frequency (MHz)	Modulation	Immunity Test Level (V/m)
5240		
5500		
5785		

	Excessive Body Motion Detector	Body movement during measurement, especially, the movement on the arm the blood pressure monitor is worn on. e.g. talking, moving or shaking of the arm with the cuff on while measurement.
		Measure again. Keep arm steady during measurement.

	Cuff	Cuff is worn improperly, or the shape of the upper arm is unusual (for example, the circumference of the upper arm differs largely from the circumference of the forearm), excessive gap might exist between the arm cuff and the arm.
		Wrap the cuff properly and keep steady. Measure again.

	Notice: The measured blood pressure reading may not be accurate if the icon is displayed.	If you have any question about the cuff wearing and/or measurement result, please consult your healthcare professional.
		If you have any question about the cuff wearing and/or measurement result, please consult your healthcare professional.

	SN	Serial number SN YYMMXXXXXX
		Ingress Protection Rating.

	IP22	First characteristic numeral - Degree of protection against access to hazardous parts and against solid foreign objects of 12.5 mm Ø g and greater.
		Second characteristic numeral - Degree of protection against ingress of water. N2=2 (Protected against vertically falling water drops when ENCLOSURE tilted up to 15°).

	Non-ionizing electromagnetic radiation.	Non-ionizing electromagnetic radiation.
		NOTE: These specifications are subject to change without notice.

	Homedics	For service or repair, do not return this unit to the retailer, Contact Homedics Consumer Relations at:
		Email: cservice@homedics.com

	Phone:	1-800-466-3342

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