

## HOW THIS BLOOD PRESSURE MONITOR WORKS

This monitor uses SmartMeasure® inflation technology to detect your blood pressure. With one touch of the start/stop button, the monitor will turn on and the inflation will automatically start, creating pressure around the arteries inside the wrist. Within the cuff is a gauge which senses the fluctuations (oscillations) in pressure. The fluctuation measured represents the degree of intensity that your arteries are contracting with each heartbeat, which is also a result of the pressure that the cuff has placed on the wrist. The monitor measures these contractions and converts the information to a digital value. Once the measurement is complete, the cuff will automatically deflate, and the result is displayed on the screen. To help guide your wrist into the ideal position, this blood pressure monitor comes equipped with an optional wrist positioning guide to help you take an accurate blood pressure reading. Please see wrist position guide section for more information on this feature. Please note that any muscle movement during inflation will cause measurement error. When measurement is complete, the monitor will display your systolic pressure, diastolic pressure, and pulse readings. The monitor automatically finds where your measurement results fall on the American Heart Association® (AHA 2017) table and provides a cue if your reading falls into one of the stages that could potentially indicate increased risk. Please see the Risk Category Index section for more information on this feature. The appearance of the icon indicates that a pulse irregularly consistent with an irregular heartbeat was detected during measurement. Refer to Irregular Heartbeat Detector section for more information on this feature.

## NAME/FUNCTION OF EACH PART

## DISPLAY SYMBOL EXPLANATIONS

Display Explanation Symbols:	Low Battery Symbol:
	Appears when battery power is excessively low, or the polarities are reversed.
	Pulse Symbol: Shows the heart rate per minute. Once pulse is detected, the symbol flashes with each pulse beat.
	Risk Category Index: The arrow points out the specific Risk Category that your measurement reading falls in. See Risk Category Index section for more information.
	Irregular Heartbeat Detector: Appears when an irregular heart beat was detected. See Irregular Heartbeat Detector section for more information.
	Memory Average: This symbol appears when LCD displays average value of last 3 readings.
	AM Averaging: Indicates the reading being displayed is an average from the last 3 mornings measurements. AM is defined as 4:00 AM – 12:00 PM.
	PM Averaging: Indicates the reading being displayed is an average from the last 3 evening measurements. PM is defined as 6:00 PM – 2:00 AM.
	Excessive Body Motion Detector: Displayed if body movement is detected during measurement, especially, the movement on the arm the blood pressure monitor is worn on. If the cuff is worn improperly, or the shape of the upper arm is unusual (for example, the circumference of the upper arm differs largely from the circumference of the forearm), excessive gap might exist between the arm cuff and the arm. <b>NOTE:</b> The measured blood pressure reading may not be accurate if the icon is displayed.
	Cuff Fitness Detection: Displayed if the cuff was wrapped incorrectly, which is too tight or too loose. This is the function aid in detecting if the cuff is wrapped properly.
	Wrist Positioning Guide: The Wrist Position Guide is used as an aid in determining if the device is at the appropriate angle and height. See Wrist Positioning Guide section for more information.
	Bluetooth® Symbol: LCD displays this symbol when Bluetooth® is active and data is transmitting to your mobile device.

If any of the following letters and numbers appear in the area that systolic pressure should be displayed, an error has occurred with your reading. See Troubleshooting section of this manual for more information.

	Measurement Error: Wrap the cuff correctly and keep wrist steady during measurement. Measure again.
	Air Circuit Abnormality: Measure again.
	Pressure Exceeding 300 mmHg: Switch the unit off to clear, then measure again.
	Error Determining Measurement Data: Rewrap cuff properly, keep steady, and measure again.
	Data Transmission Error: Monitor cannot connect to the mobile device to transmit data. Make sure Bluetooth® is ON.
	System Error: Contact Consumer Relations.

## RISK CATEGORY INDEX

This monitor comes equipped with a Risk Category Index that automatically compares each reading to defined levels established by the American Heart Association® (AHA 2017) as described earlier in this manual, and provides a helpful cue if your reading falls into one of the stages that could potentially indicate increased risk. Please note that the cues provided by this monitor are only intended to assist you in using this table. The table and cues are only provided for convenience to help you understand your non-invasive blood pressure reading as it relates to the AHA 2017 information. They are not a substitute for a medical examination or diagnosis by your physician. It is important for you to consult your physician regularly. Your physician will tell you your normal blood pressure range as well as the point at which you may actually be considered to be at risk.

BLOOD PRESSURE CATEGORY	SYSTOLIC mm Hg (upper number)	DIASTOLIC mm Hg (lower number)	Indicator Color
NORMAL	LESS THAN 120	and LESS THAN 80	Green
ELEVATED	120–129	and LESS THAN 80	Yellow
HIGH BLOOD PRESSURE (HYPERTENSION) STAGE 1	130–139	or 80–89	Red
HIGH BLOOD PRESSURE (HYPERTENSION) STAGE 2	140 OR HIGHER	or 90 OR HIGHER	Red
HYPERTENSION CRISIS (consult your doctor immediately)	HIGHER THAN 180	and/or HIGHER THAN 120	Red

\*Source: American Heart Association® (AHA) 2017

## IRREGULAR HEARTBEAT DETECTOR

The appearance of the icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Usually this is not a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. Please note that the device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage. Movement, shaking or talking during the measurement can result in pulse irregularities that may cause the appearance

of this icon. Therefore, it is of great importance to not move or talk during measurement. To determine the presence of an irregular heartbeat, the average of the heartbeat intervals is calculated with the first 3 normal effective heartbeat values. It is important to note that the average is not a strict mathematical averaging of all recorded intervals. At least 3 beats with 25% or greater difference from the average heartbeat interval will generate the icon on the screen.  
**IMPORTANT INFORMATION:** This blood pressure monitor is not designed for use by people with arrhythmias nor for diagnosing or treating an arrhythmia problem. As a safeguard, we recommend that if you have arrhythmias such as atrial or ventricular premature beats and atrial fibrillation or any other special conditions you should check with your physician before using your blood pressure monitor.

## INSTALLING BATTERIES

- The battery cover is located on the back of the monitor. Remove battery cover by pressing down and pulling away from the monitor (Fig. 1).
- Install or replace batteries according to the polarity indications inside the compartment (Fig. 2).  
**NOTE:** For easy battery removal, place batteries on top of the ribbon located inside of the battery compartment.
- Replace the cover by inserting it back into place (Fig. 3). Replace the batteries if:  
The low battery symbol appears. When any button is pressed and nothing is displayed on the screen.

- NOTE:**
  - If batteries are removed or replaced, the date and time will need to be re-set either manually or automatically using your Bluetooth® mobile device.
  - Replace all batteries at one time (as simultaneous set). Use only 1.5V "AAA" alkaline batteries. DO NOT mix alkaline, standard (carbon-zinc) or rechargeable (cadmium) batteries.
  - When the batteries are removed, the measurement values stored in memory are retained. However, the date and time must be reset.
  - Remove batteries when unit is not in use for extended periods of time.
  - Clean contacts on battery and in battery compartment with a soft, dry cloth each time you install batteries.
- Batteries are hazardous waste. DO NOT dispose of them together with the household garbage.**
- DO NOT dispose of batteries in fire. Batteries may explode or leak.**
- Recycle or dispose of properly in accordance with local, state, province, and country regulations.**

## DATE & TIME SET PROCEDURE

Date and time can be set by two methods; either sync automatically using your Bluetooth® mobile device, or set manually using the **START/STOP** button on the monitor. If setting the date and time by syncing with your mobile device, it is important that this is done prior to taking any measurements to ensure the date and time are accurate.

- To set using your Bluetooth® mobile device:**
- Before taking the first measurement, open the **Homedics Health+** app on your device. If you have not already installed the **Homedics Health+** app onto your device, it is available to download on the App Store® and on Google Play™. Make sure you have the app downloaded and open on your mobile device before trying to sync the date and time.
  - Follow the prompts in the app to select your device to sync with the app.
  - Once the Bluetooth® connection is established with the device follow the prompt on the app to sync the date & time on the device.

- NOTE:**
- Once the date and time have been successfully synced, future readings will automatically have the correct date and time.
  - If you are having trouble automatically syncing the date and time, make sure Bluetooth® is ON on your device and ON on the blood pressure monitor.


- To set manually:**
- To set the date and time, in standby mode, press and hold the **START/STOP** button for 3 seconds.
  - The display will show a blinking number showing the YEAR. Change this by pressing the **M+/%** button. Each press will increase the number by one in a cycling manner. Press the **START/STOP** button again to confirm the entry, and the screen will show a blinking number representing the MONTH.
  - Change the MONTH, DAY, HOUR, and MINUTE as described in Step 2 above, using the **M+/%** button to change the numbers and **START/STOP** button to confirm the entries.

- NOTE:**
- The date and time will only need to be set manually if the monitor will not be used with the **Homedics Health+** app.
  - If the **Homedics Health+** app is used after manually setting the date and time on the blood pressure monitor, the date and time on the mobile device will override the date and time on the blood pressure monitor.

## TALKING FUNCTION

The Talking Function will guide you through the measurement procedure and also announce measurement results. This monitor has Talking Function ON for your convenience. The Talking Function is not a substitute for reading and understanding the instruction manual.


- The Talking Function has Two Modes:
- Talking Function ON Mode - Announces instructions and results
  - Talking Function OFF Mode - No voice or announcements
- Turning the Talking Function ON/OFF:
- After setting year, date and time (please refer to the setting year, date and time section), press the **START/STOP** button to enter the Talking function.
  - Press the **M+/%** button to turn the Talking function ON or OFF.
  - Press the **START/STOP** button to switch to the Bluetooth function setting.

## BLUETOOTH® OPERATION

This monitor has Bluetooth® function ON for your convenience. This will allow your readings to automatically transmit to the **Homedics Health+** app. By turning this function OFF, the measurements cannot be transmitted. Before attempting to sync the device with your smart device, make sure Bluetooth® function is turned ON in both your smart device and the monitor, and make sure your Bluetooth-enabled device has the **Homedics Health+** app downloaded and open.

- Turning Bluetooth® Function ON/OFF:**
- After setting year, date and time (please refer to the setting year, date and time section) and talking function (please refer to the talking function section), users can press the **START/STOP** button to enter the Bluetooth® function.
  - Press the **M+/%** button to turn the Bluetooth® function ON or OFF.
  - Press the **START/STOP** button to switch to Wrist Position Guide setting.
- NOTE:**
- When Bluetooth® is OFF, measurements cannot be transmitted to the **Homedics Health+** app.
  - The **Homedics Health+** app must be open on your mobile device in order to establish a connection with the monitor.

- Bluetooth® compatibility with blood pressure monitor for Bluetooth-enabled device is:
  - Bluetooth® 4.2 for Android 6.0 or above
  - Bluetooth® 4.2 for iOS 7.0 or above

## WRIST POSITIONING GUIDE SETTINGS

- NOTE:** By default, the Wrist Positioning Guide is ON. In standby mode, press the **START/STOP** button. The device will enter to setting YEAR, DATE, MONTH, DAY, HOUR, MINUTE, Talking, and Bluetooth® function first. When these setting are done, press the **START/STOP** button again, use the **M+/%** button to turn the Wrist Positioning Guide ON and OFF. When the settings are complete, press **START/STOP** button to confirm the entries. The device will switch to sleep mode and is ready to measure.

## WRIST POSITIONING GUIDE

- This model has a built-in Wrist Positioning Guide function used as an aid to help the user determine if the device is at the appropriate height.
- When using wrist blood pressure monitors, it is important to take the measurement while relaxing your wrist at heart level. This will help ensure an accurate reading will be taken. Because it is difficult to find the ideal location for a wrist measurement, we have included an optional Wrist Positioning Sensor that will assist in directing your wrist to the ideal location for a blood pressure measurement.
- After you press **START/STOP** button, the display will illuminate with different icons that are designed to help you move your wrist. Once the appropriate height is found, the Pulse Symbol will flash and measurement will begin. Keep your arm still until measurement is completed.

Symbol	Action
	Move your wrist up.
	Move your wrist down.

**NOTE:** Even if the monitor is not in the ideal position after 10 seconds, the measurement will still begin. If this occurs, the wrist symbol will appear and be stored in the memory with the measurement.  
**NOTE:** Due to the differences in individual sizes and physique, this feature may not be helpful in all cases and you may wish to turn this feature OFF. If you feel the suggested wrist position does not match your heart level, please turn this feature OFF and consult your health care provider.

## APPLYING THE CUFF

- Remove all watches, wrist jewelry, etc. prior to attaching the wrist monitor. Clothing sleeves should be rolled up and the cuff should be wrapped around bare skin for correct measurements.
  - Apply cuff to left wrist with palm facing up as shown in Fig. A.
  - Make sure the edge of the cuff is about 0.4"(10mm) from the palm as shown in Fig. B.
  - In order to ensure accurate measurements, fasten the cuff securely around your wrist so there is no extra space between the cuff and the wrist as shown in Fig. C. If the cuff is not wrapped tight enough, the measurement values will not be accurate.
  - If your physician has diagnosed you with poor circulation in your left arm, place the cuff around your right wrist as shown in Fig. D.
- NOTE:**
- This device should not be used when your wrist is wounded or injured.
  - If it is not possible to use the cuff on your left wrist, it can also be used on your right wrist. However, all measurements should be made using the same wrist. If the cuff does not stop inflating, unwrap the cuff at once.

## CORRECT MEASUREMENT POSTURE

It is extremely important that the cuff be at the same height as the heart. Having the cuff higher or lower may cause inaccurate results.

- Place your elbow on a table so that the cuff is at the same level as your heart as shown in Fig. 1. Relax your entire body, especially between your elbow and fingers.  
**NOTE:** Your heart is located slightly below your armpit.
- If the cuff is not at the same level as your heart or if you cannot keep your arm completely still throughout the reading, use a soft object such as a folded towel to support your wrist as shown in Fig. 2.
- Turn your palm upwards.
- Sit upright in a chair, and take 5–6 deep breaths.  
**NOTE:** Avoid leaning back while the measurement is being taken as shown in Fig. 3.

## MEASUREMENT PROCEDURE

- NOTE:**
- Position the monitor at the same level as your heart during measurement to ensure accurate readings.
  - Blood pressure changes with every heartbeat and is in constant fluctuation throughout the day.
  - Blood pressure measurement can be affected by the position of the user, his or her physiologic condition and other factors.

- Before Measurement:**
- To help ensure an accurate reading, follow these instructions before taking a measurement:
    - Wait 30 minutes after exercising, eating, drinking beverages with alcohol or caffeine, or smoking to measure blood pressure.
    - Sit quietly and rest for 5 minutes.
    - Empty your bladder, as a full one can temporarily raise blood pressure
    - Stress raises blood pressure. Avoid taking measurements during stressful times.
  - Take your reading in a comfortable environment as measurements can be affected by hot or cold temperatures. Take your blood pressure at normal body temperature.

- During Measurement:**
- DO NOT talk or move your arm or hand muscles.
  - DO NOT cross your legs. Sit with feet flat on the floor.
  - DO NOT touch cuff or monitor during measurement.
- With the cuff wrapped around your wrist, press the **START/STOP** button. DO NOT start the measurement process unless the cuff is wrapped around your wrist. All digits will light up, checking the display functions. The checking procedure will be completed after about 1.5 seconds.

If you are using this blood pressure monitor for the first time, please remove the protective film from the screen.

- Adjust the position of your wrist according to the Wrist Positioning Guide symbol that appears on the display. When the device senses that your wrist is in the appropriate height, the Pulse Symbol will flash with beep sound 3 times and measurement will begin.  
**NOTE:** When the Wrist Positioning Guide Function is turned OFF, this instruction will not appear.
- Even if the monitor is not in the ideal position after 10 seconds, the measurement will still begin. If this occurs, the wrist symbol will appear and be stored in the memory with the measurement.
- As the cuff inflates, the monitor automatically determines your ideal level. This monitor detects your blood pressure and pulse rate during inflation. The Pulse Symbol flashes at every pulse beat. Remain still and DO NOT move until the entire measurement process is completed. The device will detect your pulse and determine the blood pressure.
- After the monitor has determined your blood pressure and heart rate, the cuff automatically deflates. Your systolic rate, diastolic rate, heart rate, corresponding Risk Category Indicator, Irregular Heartbeat Detector, Excessive Body Motion Detector and Wrist Position Guide (if any) are displayed with date and time for 1 minute and save results to memory automatically.
- If Bluetooth is connected successfully to the the device and the App is open, the device will activate the Bluetooth function automatically and the Bluetooth symbol will flash on the screen." If the monitor cannot be connected to paired Bluetooth-enabled device over 4.5 seconds, the LCD will display the Error message "E4" and Bluetooth will be turned off.  
**NOTE:**
  - This monitor automatically turns off approximately 1 minute after last operation. You may also press the **START/STOP** button to turn the unit off.
  - To interrupt the measurement, you may press the **START/STOP** button. The cuff will deflate immediately after this button is pressed.
  - If the cuff does not stop inflating, remove the cuff at once.

## TRANSMIT READING TO YOUR MOBILE DEVICE

This monitor syncs your blood pressure readings to the **Homedics Health+** app that is free to download on the App Store® and on Google Play™. Make sure you have the **Homedics Health+** app downloaded, and open on your mobile device before trying to transmit your blood pressure measurements. Before attempting to sync the blood pressure monitor with your mobile device, make sure Bluetooth® is turned ON on both your mobile device and the monitor.  
**Automatically transmit readings:** After a measurement is taken, the Bluetooth® icon will appear on the screen as the monitor automatically transmits your blood pressure readings to the App.  
**Manually transmit readings:** In sleep mode, press and hold the **M+/%** button for 3 seconds to activate the Bluetooth® function, and the Bluetooth® icon will flash on the screen. The readings will automatically be transmitted to the App. The Bluetooth® icon will appear on the screen as the monitor transmits your blood pressure readings to the app. If the transmission(s) were unsuccessful, the E4 symbol will be displayed and Bluetooth® will be turned off.  
**NOTE:**

- The measurement is stored in the monitor's memory even if the reading is not transmitted to your mobile device.
- The **Homedics Health+** app must be open on your mobile device in order to transmit your measurements.
- Only new readings will be accepted by the App.
- This monitor can only pair up with one Bluetooth®-enabled device at a time.

To ensure readings transmitted to the App have the correct date and time, it is important that the correct date and time is set on the blood pressure monitor before taking measurements.

- Measurements transmitted to the App cannot be edited.

### RECALLING VALUES FROM MEMORY

The monitor features 1 user memory capability, storing the last 120 measurements, and automatically replacing the oldest measurement with a new one. This monitor features a memory mode to provide you with a variety of options for reviewing your measurement history, including AM/PM averaging. Reviewing your morning (AM) and nighttime (PM) measurements can provide important information about how your blood pressure changes throughout the day.

- Press the **M+/3** button to enter Memory Mode. The monitor will first display the calculated average applied to the last 3 measurements.

**NOTE:** If there is no data stored, nothing (except month, date, and time) will appear on the display. If the records are less than 3, it will display the latest record.

- Press the **M+/3** button again and the monitor will display the average of the last 3 AM readings. Press the **M+/3** button again to display the average of the last 3 PM readings.

- Every new press of the **M+/3** button will recall a previous reading. The latest reading will be recalled first.

- To exit memory mode, press **START | STOP** button, and switch to sleep mode.

SYMBOL	STATUS
<b>AVG. 3</b>	An average of last 3 measurements.
<b>AM</b>	An average of last 3 morning measurements. (4:00AM - 12:00PM)
<b>PM</b>	An average of last 3 nighttime measurements. (6:00PM - 2:00AM)

**NOTE:** If the wrist positioning guide is ON before the measurement and the monitor is not in the ideal position, the  icon will be displayed with the saved readings.

### CLEARING VALUES FROM MEMORY

- Press the **M+/3** button to access the Memory Recall Mode.
- Press and hold the **M+/3** and **START | STOP** buttons at the same time and the data for the selected user will be erased automatically.
- To confirm deletion, press the **M+/3** button and no data should appear.

**NOTE:** Once deleted, your readings cannot be restored.

### IMPORTANT NOTES REGARDING YOUR BLOOD PRESSURE MEASUREMENT

- It is suggested that you take your measurements at the same time each day and use the same wrist for consistency.
- Users should wait a minimum of 5 minutes before taking additional measurements. More time may be necessary depending upon your physiology.
- The measurement results that users receive are for reference only. If users have any blood pressure concerns, please consult a physician.
- Once inflation reaches 300 mmHg, the unit will deflate automatically for safety reasons.
- This device is not suitable for people with arrhythmias.
- This device may have difficulty determining the proper blood pressure for pregnant women and for users with irregular heartbeat, diabetes, poor circulation of blood, kidney problems, or for users who have suffered a stroke.

### CARE, MAINTENANCE, & CLEANING

- Clean the device and cuff with a soft dry cloth, or Use a dry cloth with water to clean the device (do not flush it directly, do not soak it in water, and hold the device dry), do not use detergent or any strong chemicals to clean the device
- Make sure the cuff is completely dry before using.
- Leaky batteries can damage the unit. Remove the batteries when the unit will not be used for a long time.
- Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.
- If the unit is stored near freezing, allow it to acclimate to room temperature before use.
- This blood pressure monitor is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems with this device, please contact Homedics Consumer Relations (contact information on warranty page).
- DO NOT immerse the unit in water as this will result in damage to the unit.
- DO NOT subject the monitor or cuff to extreme temperatures, humidity, moisture, or direct sunlight. Protect from dust.
- DO NOT fold the cuff tightly.
- DO NOT disassemble the monitor or cuff. If in need of repair, refer to the warranty section of this manual.
- DO NOT subject the monitor to extreme shocks (DO NOT drop on floor).
- DO NOT inflate the cuff unless wrapped around wrist.
- DO NOT wrap the cuff around body parts other than your wrist.
- DO NOT drop or insert any object into any opening.
- This monitor may not meet its performance specifications if stored or used outside of these temperature and humidity ranges:

**Storage/Transportation Environment Temperature:** -13°F ~ 158°F(-25°C ~ 70°C)

**Humidity:** Less than 93% RH

**Operation Environment Temperature:** 41°F ~ 104°F(5°C ~ 40°C)

**Humidity:** 15 ~ 93% RH

**Atmospheric Pressure:** 700hPa-1060hPa

### POTENTIAL FOR ELECTROMAGNETIC INTERFERENCE

To avoid inaccurate results caused by electromagnetic interference between electrical and electronic equipment, Do not use the device near a cell phone or microwave oven. For most wireless communication devices, it is recommended to maintain a distance of 10.8 feet (3.3m) in order to avoid electromagnetic interference. This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

**FEDERAL COMMUNICATIONS COMMISSION COMPLIANCE STATEMENT**  
Changes or modifications to this equipment not expressly approved by the manufacturer could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

FCC RF Radiation Exposure Statement:

- The equipment complies with FCC RF exposure limits set forth for an uncontrolled environment.
- The equipment must not be co-located or operating in conjunction with any other antenna or transmitter.

- HL168AZ/BPW-H190BT essential performance per IC 80601-2-30 additional essential performance requirements:
- 201.12.1.102 Limits of the error of the manometer from environmental conditions Over the temperature range of 5°C to 40°C (41°F ~ 104°F) and the relative humidity range of 15% to 93% (non-condensing), the maximum error for the measurement of the CUFF pressure at any point of the NOMINAL measurement range shall be less than or equal to ± 3 mmHg (± 0.4 kPa) or 2 % of the reading, whichever is greater.
  - 201.12.1.107 Limits of the change in error of the blood pressure determination The laboratory limits of the change in error of the BLOOD PRESSURE DETERMINATION of the AUTOMATED SPHYGMOMANOMETER shall be less than 3 mmHg (0.4 kPa).

### ELECTROMAGNETIC COMPATIBILITY (EMC)

**Guidance and manufacturer's declaration – electromagnetic emissions**  
The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	Group 1	RF energy is used only to maintain device's operation. Therefore, its RF emissions are so low that it's not likely to cause any interference in nearby.
RF emissions CISPR 11	Group B	
Harmonic Emissions IEX 61000-3-2	Not Applicable	The device is suitable for use in all establishments, including domestic establishments, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not Applicable	

**Guidance and manufacturer's declaration – electromagnetic immunity**  
The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment- Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	In the case of air discharge testing, the climatic conditions shall be within the following ranges: Ambient Temperature: 15°C ~ 35°C Relative Humidity: 30% ~ 60%.
Power frequency magnetic field IEC 61000-4-8	30 A/m 50 or 60 Hz	30A/m 50 or 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

**Guidance and manufacturer's declaration – electromagnetic immunity**  
The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment- Guidance
Conducted RF IEC 61000-4-6	3V rms At 0.15-80 MHz 6V rms At ISM & RadioAmateur Freq	Not Applicable	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3 Proximity fields from RF wireless communications equipment IEC 61000-4-3	10 V/m at 80-2700 MHz AM Modulation And 9-28V/m at 385-6000MHz,- Pulse Mode and other Modulation. The system shall be tested as specified in IEC60601-1-2 Table 9 for proximity fields from RF wireless communications equipment using the test methods specified in IEC 61000-4-3	10 V/m at 80-2700 MHz AM Modulation And 9-28V/m at 385-6000MHz, Pulse Mode and other Modulation. The system shall be tested as specified in IEC60601-1-2 Table 9 for proximity fields from RF wireless communications equipment using the test methods specified in IEC 61000-4-3	Recommended separation distance Considering to reduce the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation: $E = 6/dVP$ where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVELS in V/m. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. <sup>3</sup> Interference may occur in the vicinity of equipment marked with the following symbol: 

**NOTE 1:** At 80 MHz and 800 MHz, the higher frequency range applies.  
**NOTE 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a: Field strengths from fixed trans-mitters, such as base stations for radio (cellular/cordless) tele-phones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast can-not be predicted theoretically with accuracy. To assess the electro-magnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be ob-served to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.
- b: Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

### TROUBLESHOOTING

If any abnormality arises or "EP" appears, during use, please check the following points:

Symbols/ Symptoms	Conditions/ Causes	Indication/Correction	
Unit does not turn on when START/STOP button is pushed.	Worn-out batteries	Replace them with 2 new AAA (LR03) alkaline batteries.	
	Battery polarities have been positioned incorrectly.	Re-insert the batteries in the correct positions.	
	Cuff has been placed incorrectly	Wrap the cuff properly so that it is positioned correctly. Keep wrist and arm steady during measurement.	
<b>EE</b> Measuring Error Symbol appears when blood pressure value displayed is excessively low or high.	Did you talk or move during measurement?	Wrap the cuff properly so that it is positioned correctly. Keep wrist and arm steady during measurement.	
	Shaking of the wrist with the cuff on.		
	Air circuit abnormality. Cuff tube may not be plugged into monitor correctly.		Check cuff connection. Measure again.
<b>E2</b> Measuring Error Symbol	Inflation pressure exceeding 300 mmHg.	Switch the unit off, then measure again.	
<b>E3</b> Measuring Error Symbol	Can't determine blood pressure measurement data.	Wrap the cuff snugly so that it is positioned correctly. If you have any question about the cuff wearing and/or measurement result, please consult your healthcare professional.	
	Cuff is worn improperly		
<b>E4</b> Measuring Error Symbol	If the device cannot be connected to paired Bluetooth-enabled device over 45 seconds, LCD will display [Error message "E4" and Bluetooth® will be turned off.	Please press <b>[ ]</b> button for 3 seconds to start Bluetooth® function	
 Excessive Body Motion Detector	Body movement during measurement, especially, the movement on the wrist the blood pressure monitor is worn on. e.g. Talking, moving or shaking of the arm with the cuff on while measurement.	Measure again. Keep steady during measurement.	
	Notice: The measured blood pressure reading may not be accurate if the icon is displayed.		
 Cuff Fitness detection Symbol	Cuff is worn improperly, or the shape of the wrist is unusual, excessive gap might be exist between the cuff and the wrist.	Wrap the cuff properly and keep steady. Measure again. If you have any question about the cuff wearing and/or measurement result, please consult your healthcare professional.	
	The cuff was wrapped incorrectly (for example too loosely or too tightly).	Please reference "Applying the Cuff" section to wrap the cuff correctly.	
	BPM cannot communicate with Bluetooth®-enabled device	Pairing has not been completed.	Please reference "Applying the Cuff" section to wrap the cuff correctly.
	Bluetooth® function is not turned on.	See the "Bluetooth® Operation" section to turn on Bluetooth® function.	
	The distance between BPM and Bluetooth®-enabled device is out of transmitting range.	Please make sure the acceptable distance (<10 meters) with each other.	
 Use an incompatible Bluetooth®-enabled device.	Use an incompatible Bluetooth®-enabled device.	Please refer to "Bluetooth® Operation" and "RF Specification" sections.	
	Use non-Bluetooth®-enabled device.		
	Unexpected loss of electrical/mechanical integrity.		Re-insert the batteries and try again.
EP error symbol shown on display.	System Error.	Measure again. If error persists, contact Consumer Relations.	

**NOTE:** If the unit still does not work, contact Homedics Consumer Relations. Under no circumstance should you disassemble or attempt to repair the unit by yourself. Contact information for Homedics Consumer Relations Department can be found in the Warranty page.

### SPECIFICATIONS

Model Number	HL168AZ / BPW-H190BT
Measurement Method	Oscillometric
Rated Range of Cuff Pressure	0 ~ 300 mmHg
Rated Range of Determination	40 ~ 280 mmHg
Measurement Range of Heart Rate	Pulse: 40 ~ 199 Beats / Minute
Accuracy	Pressure: ± 3 mmHg Pulse: ± 5 % Max.
Inflation	Automatic Inflation (Air Pump)
Deflation	Automatic (Passive Exhaust Valve)
Display	Liquid Crystal Display
Memory	120 Memories Total for 1 User
Unit Dimensions	73.79 mm x 60.15 mm x 22.94 mm (LxWxH) 2.9 x 2.36 x 0.9 inch (LxWxH)
Unit Weight	53.1 g ±3 g (0.11lbs±0.006 lbs) (Without Cuff & Batteries)
Cuff Size	5.3 ~ 7.7 inch (approx.135 ~ 195 mm)
Storage/Transportation Environment	Temperature: -25°C ~ 70°C (-13°F ~ 158°F) Humidity: ≤ 93% R.H.
Operation Environment	Temperature: 5°C ~ 40°C (41°F ~ 104°F) Humidity: 15% ~ 93% R.H. Atmospheric pressure: 700 hPa ~ 1060 hPa
Power Supply	DC 3V AAA "LR03" (1.5V) Alkaline Battery x 2
Battery Life	Approx. 250 Measurements

Shelf life (battery)	3 years (Temperature: 20 ± 2°C; Relative humidity: 65 ± 20% RH)
Product Life	5 Years (4 times per day)
Sleeping Mode	Without any operation for 1 minute, device automatically shuts off.
Accessories	Instruction Manual, 2 AAA (LR03) Alkaline Batteries, Storage Pouch
RF Type	Bluetooth® 5.0 BLE
System requirement of the Bluetooth-enabled device	Bluetooth® 4.2 for Android 6.0 or above Bluetooth® 4.2 for IOS 7.0 or above

\* The contents of this manual and the specifications of the device covered by this manual are subject to change for improvement without notice.

### EXPLANATION OF SYMBOLS

	Follow instruction for use.
	TYPE BF Applied Part.
	To avoid inaccurate results caused by electromagnetic interference. <b>WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30cm (12 inches) to any part of the device. Otherwise, degradation of the performance of this equipment could result.</b>
	Waste of electrical and electronic equipment (WEEE). Discard the used product to the recycling collection point according to local regulations.
	Serial number [SN] YYYMMXXXXXX
<b>IP22</b>	Ingress Protection Rating.First characteristic numeral - Degree of protection against access to hazardous parts and against solid foreign objects N1=2 (Protected against solid foreign objects of 12.5 mm ø and greater). Second characteristic numeral Degree of protection against ingress of water. N2=2 (Protected against vertically falling water drops when ENCLOSURE tilted up to 15°).
	Non-ionizing electromagnetic radiation.

Innovation, Science and Economic Development Canada ICES-003 Compliance  
Label: CAN ICES-3 (B)/NMB-3(B)

### Homedics

For service or repair, DO NOT return this unit to the retailer. Contact Homedics Consumer Relations at:

Email: [cservice@homedics.com](mailto:cservice@homedics.com)  
Phone: 1-800-466-3342  
Business Hours: 8:30am-7pm EST Monday-Friday

Distributed by: Homedics, LLC  
3000 N Pontiac Trail  
Commerce Township, MI 48390  
Printed in Vietnam

**LIMITED 5-YEAR WARRANTY**  
Homedics sells its products with the intent that they are free of defects in manufacture and workmanship for a period of 5 years from the date of original purchase, except as noted below. Homedics warrants that its products will be free of defects in material and workmanship under normal use and service. This blood pressure monitor meets the simulated measurement cycles test requirement per EN1060-3, part 8.10. This warranty extends only to consumers and does not extend to Retailers. To obtain warranty service on your Homedics product, contact a Consumer Relations Representative. Please make sure to have the model number of the product available. Homedics does not authorize anyone, including, but not limited to, Retailers, the subsequent consumer purchaser of the product from a Retailer or remote purchasers, to obligate Homedics in any way beyond the terms set forth herein. This warranty does not cover damage caused by misuse or abuse; accident; the attachment of any unauthorized accessory; alteration to the product; improper installation; unauthorized repairs or modifications; improper use of electrical/power supply; loss of power; dropped product; malfunction or damage of an operating part from failure to provide manufacturer's recommended maintenance; transportation damage; theft; neglect; vandalism; or environmental conditions; loss of use during the period the product is at a repair facility or otherwise awaiting parts or repair; or any other conditions whatsoever that are beyond the control of Homedics. This warranty is effective only if the product is purchased and operated in the country in which the product is purchased. A product that requires modifications or adoption to enable it to operate in any other country than the country for which it was designed, manufactured, approved and/or authorized, or repair of products damaged by these modifications is not covered under this warranty.

THE WARRANTY PROVIDED HEREIN SHALL BE THE SOLE AND EXCLUSIVE WARRANTY. THERE SHALL BE NO OTHER WARRANTIES EXPRESS OR IMPLIED INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS OR ANY OTHER OBLIGATION ON THE PART OF THE COMPANY WITH RESPECT TO PRODUCTS COVERED BY THIS WARRANTY. HOMEDICS SHALL HAVE NO LIABILITY FOR ANY INCIDENTAL, CONSEQUENTIAL OR SPECIAL DAMAGES. IN NO EVENT SHALL THIS WARRANTY REQUIRE MORE THAN THE REPAIR OR REPLACEMENT OF ANY PART OR PARTS WHICH ARE FOUND TO BE DEFECTIVE: WITHIN THE EFFECTIVE PERIOD OF THE WARRANTY, NO REFUNDS WILL BE GIVEN. IF REPLACEMENT PARTS FOR DEFECTIVE MATERIALS ARE NOT AVAILABLE, HOMEDICS RESERVES THE RIGHT TO MAKE PRODUCT SUBSTITUTIONS IN LIEU OF REPAIR OR REPLACEMENT.

This warranty does not extend to the purchase of opened, used, repaired, repackaged and/or resold products, including but not limited to sale of such products on internet auction sites and/or sales of such products by surplus or bulk resellers. Any and all warranties or guarantees shall immediately cease and terminate as to any products or parts thereof which are repaired, replaced, altered, or modified, without the prior express and written consent of Homedics. This warranty provides you with specific legal rights. You may have additional rights which may vary from state to state. Because of individual state regulations, some of the above limitations and exclusions may not apply to you.

For more information regarding our product line in the USA, please visit: [www.homedics.com](http://www.homedics.com)

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©2022 Homedics, LLC. Reservados todos los derechos. Homedics y la marca denominativa y el logotipo de Smart Measure son marcas comerciales registradas de Homedics, LLC. American Heart Association es una marca registrada de American Heart Association, Inc. La marca denominativa y los logotipos de Bluetooth® son marcas comerciales registradas propiedad de Bluetooth SIG, Inc y cualquier uso de dichas marcas por parte de Homedics, LLC se realiza bajo licencia. App Store es una marca de servicio de Apple Inc. Google Play y el logotipo de Google Play son marcas comerciales de Google LLC. B-BPW-H190BT P/N: 323103768
DISTRIBUIDO por: Homedics, LLC Sendero Pontiac 3000 N Municipio de Comercio, MI 48390